# sta-seal f

DE

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# sta-seal f

**GB**  
Instructions for use

**Precision impression material, silicone based, condensation curing, heavy bodied consistency**

Fig. 1

## 1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive** spray or **sili Adhesive** fluid onto the tray and let dry for 60 sec (fig. 1).

## 2. Mixing

Extrude **sta-seal f** from the dosing syringe onto mixing pad (fig. 2). Dosage: 1 cm **sta-seal f** to 1 drop of liquid catalyst. Mixing 30 sec. with spatula and then hand-kneading with fingers until a homogeneous mix is achieved. Load tray. For easy application of individual border use easy application syringe (fig. 3).

Fig. 2

## 3. Disinfection

After removal from mouth rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

## 4. Model casting

Casting not before 2 hours at the latest after 24 hours. Recommended materials: Dental plasters of class III and IV.

Fig. 3

## 5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

## 6. Cleaning of trays

Remove set material with a blunt instrument. **sili adhesive** film can be easily removed with **sili spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.


## Important working hints

- Overdosing of catalyst or increased temperatures accelerate, underdosing or decreased temperatures retard the setting time.
- Catalyst bottles should be closed tightly immediately after use.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes. Wear suitable eye protection and gloves! It may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.
- Don't leave any residual material in the patient's mouth.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

## Indications for use:

- Special impression material for functional impression with extended setting
- Functional border moulding

## Technical Data:

DIN EN ISO 4823 – Type 1

- Mixing ratio:** **sta-seal f / sta-seal catf**  
1 cm sta-seal f = 1.70 g / 1 drop catf = 0.02 g
- Colour code:** **sta-seal f:** green **sta-seal catf:** light blue
- Mixing time:** approx. 30 sec.
- Working time:** approx. 2 min.\*
- Setting time:** approx. 6 min.\*
- Time in mouth:** approx. 4 min. incl. 2 min. 30 sec. plastic phase
- Strain in compression:** approx. 12.7 %
- Recovery from deformation:** approx. 99 %
- Linear dimensional change:** approx. 0.4 %
- Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

## Ordering information:

**sta-seal f**  
Standard packing **02012**  
3 x 80 ml dosing syringes  
10 ml sta-seal catf bottle

Large packing **02013**  
9 x 80 ml dosing syringes  
3 x 10 ml sta-seal catf bottles

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