# Detaseal® hydroflow putty

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# Detaseal<sup>®</sup> hydroflow putty

**GB**  
Instructions for use

VPS precision impression silicone, kneadable consistency  
– fast set

Fig. 1

## 1. Tray preparation

Select a tray type suited for desired application (standard-size or custom, in metal or hard plastics). To avoid teeth contact with the tray, place stops from wax or light-curing resin into the tray at non-relevant areas of the impression. For optimal adhesion of impression material, we recommend applying a thin layer of **trayloc<sup>®</sup> A** to tray (Fig. 1). Let dry for approx. 3 minutes.


Fig. 2

## 2. Mixing and dispensing

Dispense in a mixing ratio of 1:1, using the enclosed measuring scoops. Mix base and catalyst thoroughly by manual kneading until a uniform colour is achieved (Fig. 2).

## 3. Using

First apply the mixed **Detaseal<sup>®</sup> hydroflow putty** onto the impression tray, giving the teeth area a distinctly trough-like shape (Fig. 3). Then apply a layer of **Detaseal<sup>®</sup> hydroflow lite** or **-Xlite** into the shaped depression (Fig. 4). Insert the filled tray into patient's mouth and position it within 3 – 4 seconds, under light pressure. Then seat the impression without pressure, until **Detaseal<sup>®</sup> hydroflow** has set.

Fig. 3

## 4. Disinfecting

After removing tray from the mouth, rinse impression under running, lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.


Fig. 4

## 5. Pouring the model

We recommend waiting at least 30 minutes before pouring the model. The impression will retain its dimensional stability for up to 14 days. All class III and IV dental plasters and standard modelling resins can be used.

Fig. 5

## 6. Cleaning the tray

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc<sup>®</sup> A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.


## Important working hints

- Substances such as retraction liquids can adversely affect the setting reaction and should therefore be thoroughly rinsed off before taking the impression.
- Do not combine with condensation curing silicones.
- Contact with latex examination gloves or other latex material may impair the setting reaction of **Detaseal<sup>®</sup> hydroflow putty**.
- Be sure to remove all traces of impression material from patient's mouth.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- We recommend use of **Detaseal<sup>®</sup> hydroflow lite** or **-Xlite** as correction material in the double mix technique for the one time impression.
- In order to enhance easy mouth removal in case of one-step impressions and to avoid blocking in strong undercuts, fill undercuts beforehand with **Detaseal<sup>®</sup> antilock**.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

**DETAX**

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 07243/510-0 · Fax: 07243/510-100  
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

Made in Germany  
**PG** **CE**  
IM05 08/2017

## Indications for use:

- Preliminary impressions for the double impression technique
- Impressions with sandwich and double-mix technique
- Situation impressions

## Technical data:

DIN EN ISO 4823 - Type 0

- **Mixed volume:** 2 x 250 ml jars
- **Mixing ratio:** 1:1
- **Product colour:** base: green catalyst: white
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 1 min. 30 sec.\*
- **Setting time:** approx. 3 min. 30 sec.\*
- **Time in mouth:** approx. 2 min.
- **Strain in compression:** approx. 3 %
- **Recovery from deformation:** approx. 99,8 %
- **Linear dimensional change:** approx. 0,2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.
- **Storage:** 15°C / 59°F to 25°C / 77°F

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperature accelerate, decreased temperatures retard a. m. times.

## Ordering information:

**Detaseal<sup>®</sup> hydroflow putty Standard pack 02727**

2 jars of 250 ml  
2 measuring scoops

**Multi packing 02728**

4 jars of 250 ml  
2 measuring scoops

**Detaseal<sup>®</sup> hydroflow lite Standard packing**

2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set 02730**  
**fast set 02733**

**Detaseal<sup>®</sup> hydroflow Xlite Standard packing**

2 cartridges of 50 ml  
12 mixing cannulas, yellow  
**regular set 02741**  
**fast set 02744**

**Detaseal<sup>®</sup> antilock Standard pack**

1 x 50 ml cartridge  
6 Mixing cannulas, yellow  
6 Intra-oral tips, yellow  
**03535**

**trayloc<sup>®</sup> A 03098**

17 ml bottle with brush  
15 ml spray bottle  
**03392**

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