DE	Gebrauchsanweisung .....	3
GB	Instructions for use .....	6
FR	Mode d'emploi .....	9
ES	Instrucciones de uso .....	12
IT	Istruzioni per l'uso .....	15
PT	Instruções de uso .....	18
NL	Gebruiksaanwijzing .....	21
SE	Bruksanvisning .....	24
DK	Brugervejledning .....	27
NO	Bruksanvisning .....	30
FI	Käyttöohjeet .....	33
GR	Οδηγίες χρήσης .....	36
CZ	Návod k použití .....	39

HU	Használati utasítás .....	42
LV	Lietošanas instrukcija .....	45
LT	Naudojimo instrukcija .....	48
PL	Instrukcja obsługi .....	51
EE	Kasutusjuhend .....	54
HR	Upute za uporabu .....	57
RO	Instrucțiuni de utilizare .....	60
TR	Kullanma talimatı .....	63
KR	<b>사용 설명서</b> .....	66
BG	Инструкции за употреба ..	69
ID	Petunjuk penggunaan .....	72
SI	Navodila za uporabo .....	75
RS	Uputstva za upotrebu .....	78

### Signum metal bond

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Stand: 2013-07

## Signum metal bond

Signum metal bond I+II is a metal resin bonding system for crown, bridge and CrCo framework.

### Composition:

Inorganic pigments, SiO<sub>2</sub> (silanised), methyl methacrylate, acetone, bifunctional methacrylate-based components

### Indication:

- Bonding methacrylate based resins / composites to metal framework surfaces in conjunction with macromechanical retention
- Bonding methacrylate based resins / composites to CrCo framework surfaces
- Bonding agent for repairing composite facings (intraoral and extraoral)
- Bonding agent for repairing porcelain bonded restorations
- Bonding between root posts, Maryland bridges and composite cement

### Technique

Sandblast the surface to be faced with min. 110 µm Al<sub>2</sub>O<sub>3</sub> at 3 bar pressure (non-recyclable abrasive). Then clean with uncontaminated (oil free) compressed air or clean tap water (do not steam clean) and dry thoroughly.

Dispense an ample amount of Signum metal bond I into a mixing dish. Apply to the prepared surface with a fine brush and allow to air dry. Apply a second coat if necessary.

Shake Signum metal bond II well to attain an optimal distribution of the colour pigments on the surface during application. Apply 1–2 uniform layers of Signum metal bond II with another brush and polymerize for 90 sec. with the Heraflash, HiLite power or UniXS. For intraoral use polymerize Signum metal bond II for 40 seconds with a dental polymerization lamp (e.g. Translux Power Blue).

For an optimum bond apply the relevant opaquer without undue delay and without contaminating the surface by touching it. Then apply the composite or acrylic.

#### **Instructions and contraindications:**

Only for use by dentists and dental technicians for the above indications.

Close the Signum metal bond I and II bottles immediately after use.

Signum metal bond II should not cover the surface completely. The shade is masked by the application of Signum opaque F.

#### **Hazard/Safety information:**

Store Signum metal bond I + II at room temperature (max. 25°C (77°F)).

The product is highly flammable. Keep away from naked flames, sources of ignition and direct sunlight.

The product can irritate the eyes, respiratory organs and skin. On contact with the eyes or skin, immediately rinse thoroughly with water and consult a doctor.

Avoid direct contact with the mucosa when carrying out repairs intraorally. Use a rubber dam if necessary.

On contact with the skin, clean immediately with soap and water.

The acetone and MMA content may cause an allergic reaction with hypersensitive patients.

This product or one of its components may in particular case cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

**Shelf life:**

The use by date is printed on the label of each bottle. Do not use after the expiry date.

**Disposal instructions:**

Observe the safety data sheet or national regulations for disposal.

**Supplied as:**

Signum metal bond Set I + II, Mixing dish

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Heraeus Kulzer GmbH