CHARISMA® flow

CHARISMA® flow

- Gebrauchsanweisung
- Instructions for use
- Mode d'emploi
- Instrucciones de uso
- Istruzioni per l'uso
- Instruções de uso
- Gebruiksaanwijzing
- Bruksanvisning
- Brugervejledning

- Bruksanvisning
- Käyttöohjeet
- Οδηγίες Χρήσης
- Návod k použití
- Használati utasítás
- Lietošanas instrukcija
- Naudojimo instrukcija
- Instrukcja obsługi
- Инструкция по применению

66034065/02 - 07.09

CHARISMA® flow

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CHARISMA® flow

Charisma flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Charisma flow provides excellent coverage of the tooth structure.

Composition: Charisma Flow is based on multifunctional methacrylate monomers (EBADMA/TEGDMA) and contains approximately 62 % by weight or 38 % by volume inorganic fillers such as Ba-Al-F silicate glass and SiO₂. The filler particle size is between 0,005 µm and 5 µm.

Shades: A1, A2, A3, A3,5, B2, B3, O2, O4, O2B, SB, Baselliner

Applications:

- Fissurenversiegelung
- Erweiterte Fissurenversiegelung
- Class V fillings
- Minimal invasive fillings in areas not subjected to masticatory forces
- Minimally invasive Class I and II fillings in areas with moisture, blood or saliva.
- Smaller contour and shade adjustment on the enamel and dentin
- Small surface repairs to direct and indirect restorations in combination with a suitable bonding agent
- Occlusion
- Cavity lining
- Bracket retention

Kind of application: Prepares tooth structure for treatment by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Charisma shade using the Charisma shade guide.

1. Remove any caries and completely invasive preparation of the tooth structure according to the general guidelines for adhesive fillings. The use of a rubber dam is recommended.

2. A quality adhesive system (e.g. iBond® Self Etch, GLUMA® Comfort® Bond + Desensitizer, GLUMA® Solid Bond) should be applied to the enamel and dentin in accordance with the manufacturer's instructions before inserting the filling.

3. Apply Charisma flow in thin layers (max. 2 mm; Baselliner max. 1 mm) to the cavity. When applying Charisma flow ensure that it flows slowly and evenly to prevent air bubbles – apply a gentle stream of air to facilitate maximum marginal contact. If necessary, remove any air bubbles with a probe before polymerization.

4. Polymerize each layer according to the table below using a suitable light-curing unit (e.g. Translux® Power Blue®) for 20 or 40 sec. A dispersion layer forms during polymerization. This layer should not be touched or removed if other composite layers are to be applied.

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