
3M ESPE

Elipar™

- (en) **LED Curing Light**
- (ru) **Светодиодное устройство для фотополимеризации стоматологических материалов**
- (bg) **Светодиодна лампа за полимеризация**
- (hr) **LED svjetiljka za polimerizaciju**
- (hu) **Gyors polimerizációt biztosító LED-lámpa**
- (pl) **Diodowa lampa do polimeryzacji**
- (ro) **Lampă de fotopolimerizare cu LED**
- (sk) **LED Polymerizačná lampa**
- (sl) **Polimerizacijska LED luč**
- (cs) **Polymerační LED lampa**
- (tr) **LED Işık Cihazı**
- (et) **LED-valguslamp**
- (lv) **Gaismas diodes lampā**
- (lt) **Diodinis švitintuvas**

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PLEASE NOTE!

Prior to installation and start-up of the device, please read these instructions carefully!

As with all technical devices, the proper function and safe operation of this device depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.

1. Use of the device is restricted to trained personnel in accordance with the instructions below. The manufacturer assumes no liability for any damage arising from any other or improper use of this advice.
2. The charger must be accessible at all times. Do not use the charger for any use other than the charging of the Elipar™ handpiece. Disconnect the handpiece from the mains by unplugging the charger from the electrical outlet.

Treating patients using the handpiece while it is still connected to the charger is prohibited for safety reasons. Light-curing is possible only if the charger has been disconnected.
3. Use only the charger (AC adapter plug) which is provided with the device. The use of any other charger can result in damage to the battery.
4. The light must not be aimed at the eyes in order to avoid serious health consequences due to irradiation of the eyes. Exposure must be restricted to the area of the oral cavity in which clinical treatment is intended. Protect patient and user from reflection and intensive scattered light by taking the appropriate measures, e.g., glare shield or coverings.
5. CAUTION! The Elipar device generates high intensity light. High light intensity always involves the production of heat. The light emitted should be aimed directly above the material to be cured – exposure of the soft tissues (gingiva, oral mucosa, and skin) to high-intensity light should be avoided as such exposure may cause damage or irritation. If applicable, cover such areas. If exposure of soft tissues cannot be avoided, adjust the polymerization process to the light level, e.g., by shortening the polymerization times or increasing the distance between the light guide exit and the material to be cured. Longer exposure in the region of the pulp can lead to injury or irritation of the pulp, which is why the specified exposure times must not be exceeded.
6. Elipar may be operated only with the supplied light guide or original 3M ESPE Elipar replacement and accessory light guide. The light guide has to be seen as an applied part. The use of other light guides may

Table of Contents	Page
Safety Instructions	1
Glossary of Symbols	2
Product Description	3
Fields of Application	3
Technical Data	3
Charger	3
Handpiece	3
Charger and Handpiece	4
Transport and Storage Conditions	4
Installation of the Unit	4
Factory Settings	4
Initial Steps	4
Charger	4
Light Guide/Handpiece	4
Battery Charging	4
Battery Power Level Display on Handpiece	5
Operation	5
Selection of Exposure Time	5
Activating and Deactivating the Light	6
Inserting and Removing the Light Guide from/into the Handpiece	6
Positioning the Light Guide	6
Testing of Light Intensity	6
Sleep Mode	6
Audible Signals – Handpiece	7
Troubleshooting	7
Maintenance and Care	8
Care of the Handpiece	8
Cleaning the Light Guide	8
Clean Handpiece and Glare Shield	9
Storage of the Handpiece during Extended Periods of Non-Use	9
Return of Old Electric and Electronic Equipment for Disposal	9
Customer Information	10
Warranty	10
Limitation of Liability	10

- result in a reduction or increase in the light intensity. The product's warranty does not cover any damage resulting from the use of third-party light guides.
7. Condensation resulting from the device being transferred from a cold to a warm environment may be a potential risk. Never begin operating the device until it has reached the ambient temperature.
 8. In order to avoid electric shock do not introduce any objects into the device with the exception of replacement parts handled in accordance with the Operating Instructions.
 9. Use only genuine 3M ESPE parts when replacing defective components as directed in these Operating Instructions. The product's warranty does not cover any damage resulting from the use of third-party replacement parts.
 10. Should you have any reason to suspect the safety of the device to be compromised, the device must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective device. Safety may be compromised, e.g., if the device malfunctions or is noticeably damaged.
 11. Keep solvents, flammable liquids, and sources of intense heat away from the device as they may damage the plastic housing of the device, the seals, or the operating buttons.
 12. Do not operate the device in the proximity of flammable mixtures.
 13. Do not allow any cleaning agents to enter the device during cleaning as they could cause an electrical short or a dangerous malfunction.
 14. Only service centers authorized by 3M Deutschland GmbH may open the device housing and repair the device.
 15. Elipar must not be used in patients, or by users, with heart pacemaker implants who have been advised to be cautious with regard to their exposure to small electrical devices.
 16. Do not use Elipar in patients with a history of photo-biological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photosensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
 17. Individuals with a history of cataract surgery may be particularly sensitive to the exposure to light and should be discouraged from Elipar treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
 18. Individuals with a history of retinal disease should seek advice from their ophthalmologist prior to operating the device. In operating the Elipar device, this group of individuals must take extreme care and comply with any and all safety precautions (including the use of suitable light-filtering safety goggles).
 19. This device has been developed and tested in accordance with the relevant EMC regulations and standards. It is in conformity with legal requirements. Since various factors such as power supply, wiring, and the ambient conditions at the place of operation can affect the EMC properties of the device, the possibility that, under unfavorable conditions, there will be EMC disruptions cannot be completely excluded. If you should notice problems in the operation of this or other devices, move the device to a different location. The EMC manufacturer's declaration and the recommended separation distances between portable and mobile RF communications equipment and the Elipar unit are listed in the appendix.
 20. Prior to each use of the device ensure that the emitted light intensity is sufficient to safely guarantee polymerization. Check the light guide and the light guide mounting hole to make sure they are clean. If necessary, the light guide mounting hole and the light guide can be cleaned as described in the section "Maintenance and Care" (see also the section "Measurement of Light Intensity").

Glossary of Symbols

Follow instructions for use.

Attention, Consult Accompanying Documents

Type B Equipment –
Protection against electric shock

Protection Class II – double insulated

Use in closed spaces only


93/42/EEC

Battery power level

Icon to identify electric and electronic devices.
The device must be collected and disposed of separately

Product Description

Elipar is a high-performance light source for the polymerization of light-curing dental materials. The device has two components: a wireless handpiece with a built-in battery which can be replaced by 3M ESPE customer care and a charger. The device is a medical electrical device in accordance with IEC 60601-1 and is available as a tabletop device. Wall mounting is not possible.

The light source is a high-performance light-emitting diode (LED). In contrast to halogen lights, the emitted light specifically covers the light wavelength between 430 and 480 nm, the relevant range e.g. for camphor quinone-containing products.

The polymerization performance is so high that the exposure times can be reduced by 50% in comparison with a conventional halogen light (with light intensity typically ranging from 600 to 800 mW/cm²). Selected filling composites can be cured in as little as 5 seconds if the light guide can be used in their immediate proximity. Manufacturer's instructions for light-curing dental materials regarding layer thickness, exposure time, and light intensity must be followed at all times. Filling composites which can be cured in 5 seconds are shown in the list included with the device.

Settable exposure times:

- 5, 10, 15, 20 sec
- Continuous mode (120 sec)
- Tack-cure mode

Place the handpiece on a flat surface when not in use. The handpiece can be connected to the charger between applications so that the battery is charged.

The handpiece must be connected to the charger at the latest when the battery power level display glows

red steadily so that the battery is charged (cf. "Battery Power Level Display on Handpiece").

The device is shipped with a light guide with 10 mm diameter. It is not permissible to use the light guides of other devices.

The handpiece is equipped with a sleep mode to minimize the device's energy consumption. The handpiece switches to sleep mode if it is not used for a period of about 5 minutes or the incorrect charging voltage is detected.

The charger uses a maximum of 0.2 W when ready for operation.

- ☞ These Operating Instructions should not be discarded for the duration of use of the device.

Fields of Application

- Polymerization of light-curing dental materials with photo initiator for the wavelength range 430–480 nm.
 - Though the majority of light-curing dental materials are responsive in this range of wavelengths, you may wish to contact the manufacturer of the material in question.

Technical Data

Charger

Operating voltage: 100–240 V 50/60 Hz

Nominal consumption: 0.2 A max

Dimensions without country-specific adapter:

Length: 65 mm
Width: 40 mm
Depth: 31 mm

Weight:

75 g

Classification:

Protection class II, 

Manufacturer:

Click Technology Co., Ltd.

Model:

CPS 008050100

Handpiece

Power supply:

lithium-ion battery,
nominal voltage 3.7 V,
capacity 2100 mAh

Utilizable wavelength range:

430–480 nm

Wavelength peak:

455 nm +/- 10 nm

Light intensity
(between
400 and 515 nm): 1200 mW/cm² -10% / +20%
(independent of battery power level)

Light emission area: 60–65 mm² (optically active)

Intermittent operation: 1 min on, 15 min off
(ambient temperature of 40° C)
typically 7 min operating time
at ambient temperature (23° C)

Total exposure time
with new, fully
charged battery: typically 60 min

Dimensions: Diameter: 28 mm
Length: 270 mm

Weight: 180 g (incl. light guide)

Charger and Handpiece

Time to charge
empty battery: approx. 2 h

Operating
temperature: 10° C up to 40° C /
59° F up to 104° F

Relative humidity: 30% up to 75%

Atmospheric pressure: 700 hPa up to 1060 hPa

Transport and Storage Conditions:

Ambient temperature
range: -20° C up to +40° C /
-4° F up to +104° F

Relative humidity: 

Atmospheric pressure:  hPa

Subject to technical modification without prior notice.

Installation of the Unit

Factory Settings



The factory settings of the device are as follows:

- 10 sec exposure time

Initial Steps

Charger

- Select the plug adapter specific to the country and place it on the charger.

Light Guide / Handpiece

- Place the glare shield on the front of the device.
- Autoclave the light guide prior to first use.
- Then attach the light guide to the handpiece until it noticeably locks into place (see Section “Removing and Inserting the Light Guide from/into the Handpiece”).
- If the device malfunctions, insert the charger plug into the charging socket of the handpiece. The device will reset itself and can then be used again.


Battery Charging

- The device contains a powerful lithium-ion rechargeable battery. This type of battery does not have any

memory effect and can therefore be recharged at any time by inserting the charger plug into the charging socket of the handpiece (see the section “Battery Power Level Display on Handpiece”).

- Before using the handpiece for the first time, connect it to the charger for a period of about 2 hours so that the new battery is completely charged for the first time.

The green status light on the handpiece blinks while the battery is charging. The green status light glows steadily when the device is fully charged. **As a safety precaution, light-curing is not possible while the device is charging.**


Battery Power Level Display on Handpiece

Status LED	Operating status	
	Handpiece without charger	Charger connected
Steady green light	Handpiece ready for operation	Charging has been completed
Flashes green	—	Battery is being charged
Steady red light	Low battery charge	Problem during charging
Flashes red	Battery fully discharged, exposure cycle will be completed or, if in continuous mode, stopped	Problem in charging, battery is defective or cannot be charged

Operation

Selection of Exposure Time

Exposure time options: 5, 10, 15, 20 sec, continuous mode (120 sec), pulse curing function.

- Due to the high light output of the Elipar device, the times of 5, 10, 15, 20 sec are equivalent to 10, 20, 30 and 40 sec using a light polymerization unit of conventional performance (between approx. 600 and 800 mW/cm² for halogen technology or between 300 and 400 mW/cm² with LED). Thus, the normal exposure times for conventional units can be cut in half without compromising polymerization performance.
- Selected filling composites can be cured in as little as 5 seconds if they can easily be accessed using the light guide. Note: position the light guide at the minimum distance from the composite and expose for 5 sec. Filling composites which can be cured in 5 seconds are shown in the list included with the device.
- Reduced exposure time is permitted only if the full light intensity of the Elipar is possible (see section “Measurement of Light Intensity”).

Select the exposure time by pressing the TIME button.


- The selected exposure time is indicated by the 4 green LEDs.
- Each time the button is briefly pressed, the setting advances to the next (higher) value. All 4 green LEDs will be turned on for a setting of 20 sec. Pressing the button again will turn off all of the LEDs and enable the continuous mode.
- The display advances through the available settings, if the button is kept depressed.
- While exposure is ongoing, the button for selection of the exposure time is inactive.

Activating and Deactivating the Light

- ▶ Briefly press the START button; the light will turn on.
 - The LEDs first indicate the set exposure time; all 4 LEDs light up for 20 sec. Every 5 sec, as the time runs down, the LEDs will turn off one at a time; at 15 sec remaining time, 3 LEDs will still be on, at 10 sec remaining time 2 LEDs, etc.
 - The LEDs do not come on at all in continuous mode; an audible signal is emitted every 10 seconds.
- ▶ If desired, the light can be turned off by pressing the START button again before the exposure time is over.
- ▶ Holding down the START button activates the tack-cure function: the device emits a single short light pulse which enables the defined curing of Protemp™ Crown temporary restorations or a light-curing cement excess (e.g., RelyX™ Unicem) to enable easy removal.

Inserting and Removing the Light Guide from/into the Handpiece

- ▶ Attach the light guide to the handpiece by inserting it with a slight rotation until it clicks firmly into place and there are no gaps between the neck of the light guide and the handpiece (see picture).
- ▶ Remove the light guide from the handpiece by pulling it towards the front.

Positioning the Light Guide

- ▶ Rotate the light guide into the desired position for polymerization.
- ▶ To make full use of the light intensity provided, place the light guide as close to the filling as possible. Avoid directly contacting the filling material!
 - Keep the light guide clean at all times to obtain full light intensity.


- **Damaged light guides substantially reduce the light intensity and must be replaced immediately, sharp edges may cause serious injury!**

Testing of Light Intensity

Depending on the layer thickness of the filling material being cured, it is possible to check the function of the curing light using curing disks for composites:

- ▶ Place the curing disk on a mixing pad.
- ▶ Fill your preferred composite into a bore of the curing disk at least twice as deep as the recommended layer thickness of your composite.
- ▶ Cure the composite in the curing disk for the time recommended by the composite manufacturer.
- ▶ Scrape off soft material from the bottom of the cured material with a plastic spatula.
- ▶ The solid thickness of the cured material in the curing disk **divided by two** is the layer thickness which can be properly cured.
- ▶ If you have doubts about the correct function of your curing light, clean the light guide and the protection glass in the light guide mounting hole and ensure proper fit of the light guide in the handpiece and repeat the intensity test.
- ▶ If you have still doubts about the correct function of your curing light, contact 3M ESPE service.

The curing disk must only be used to test the function of the curing light! For clinical depth of cure, please refer to the Instructions for Use of the filling material.

Commonly available devices for measuring light intensity can also be used; their measurement values should not be regarded as absolute values. When such devices are used, we recommend recording the intensity of the curing light before its first use and measuring it again at regular intervals so that any

decrease in light intensity can be detected. The light intensity testing unit in an Elipar S10 base station has a device for measuring intensity on a percentage basis.

Sleep Mode

The handpiece begins charging automatically when connected to the charger (green status light blinks) if the battery requires charging. If the charging voltage is not correct (e.g., there is dirt on the contacts in the charger socket or the charger), the handpiece goes into sleep mode. If the handpiece is not connected to the charger and is not used for a period of about 5 minutes, it also goes into sleep mode. In this operation mode, all of the displays and signals of the handpiece are turned off so that power consumption is reduced to a minimum. To terminate the sleep mode press the START button.

- The sleep mode termination signal (two short audible signals) is emitted indicating that the handpiece is ready for operation; the handpiece displays the latest selected exposure mode and time settings.

Audible Signals – Handpiece

An audible signal is emitted

- every time a button is pressed,
 - every time the light is turned ON,
 - 1 time after 5 sec exposure time, 2 times after 10 sec, 3 times after 15 sec.
- Exception: in continuous mode; an audible signal is emitted every 10 seconds.

Two audible signals are emitted

- every time the sleep mode is terminated by pressing the START button,
- every time the light is turned OFF.

A 2 sec-error signal is emitted, if

- the handpiece over-heats,
- the battery lacks sufficient charge.

The audible signals from the handpiece can be turned off (except for the 2-sec error signal). Follow these instructions to turn them off. Plug the charger into a functioning power outlet, then press and hold down the TIME and START buttons at the same time with one hand. With the other hand, connect the charging cable from the power supply to the charging socket of the handpiece. An audible signal confirms that there has been a change from "Audible signals activated" to "Audible signals deactivated". Release the two buttons and disconnect the handpiece from the charger. Repeat the above procedure to activate the audible signals.

Troubleshooting

Error	Cause ► Solution
The status display on the handpiece glows red steadily.	The remaining battery charge is adequate only for typical 5 10-sec exposure cycles. ► Connect the handpiece to the charger and recharge the battery.
The status display of the handpiece flashes red. The ongoing exposure is interrupted (light off signal is emitted) followed by a 2 sec-error signal; the handpiece switches to "sleep" mode and resists further activation.	The battery lacks sufficient charge. ► Connect the handpiece to the charger and recharge the battery.
The status light on the handpiece blinks red whenever the handpiece is connected to the charger.	Charging problem. The battery is defective or at the end of its useful life. ► Contact 3M ESPE Service.
The handpiece has not been used for a long time and now it cannot be turned on.	There is not enough charge in the battery to turn on the handpiece. ► Connect the handpiece to the charger and recharge the battery.
The handpiece does not respond to the pressing of either button.	Software crash possible. ► Plug the included charger into an outlet and connect it to the handpiece. This causes the light-curing device to reset itself.
While the handpiece is connected to the charger, pressing the START button will not start the light emission.	The handpiece is connected to the charger. Light-curing is prevented as a safety precaution. ► Disconnect the charger from the handpiece and restart the light emission.

Error	Cause ► Solution
The light emission does not start when the START button is pressed; an error signal sounds for 2 sec.	The information about intermittent operation under the header “Technical data – handpiece” has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.
During light emission in continuous mode, an error signal sounds for 2 sec, the emission is stopped, and the handpiece changes to sleep mode.	The information about intermittent operation under the header “Technical data – handpiece” has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.
The light intensity is too low	► Clean the light guide and the protecting glass in the light guide mounting hole (please refer to “Cleaning the Light Guide”).

Maintenance and Care

The Elipar device is maintenance-free. No periodic maintenance is required. See the information contained in this chapter to secure problem-free operation.

Care of the Handpiece

- Use only the charger included with the product. The use of other chargers may damage the battery cells or result in inadequate charge.

Do not immerse the handpiece in water or incinerate! Please also observe the chapter on “Safety”.

Cleaning the Light Guide

Clean and disinfect the light guide before every use. The light guide is not sterile when delivered and must be autoclaved before being used for the first time.

Material Resistance

Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:

- Organic, mineral, and oxidizing acids (minimum acceptable pH value 5.5)
- Bases (maximum acceptable pH value 8.5)
- Oxidation agents (e.g., hydrogen peroxide)
- Halogens (chlorine, iodine, bromide)
- Aromatic/halogenized hydrocarbons

Please check the manufacturer’s information about the cleaning and disinfecting agents.

The light guide must not be exposed to temperatures higher than 134° C (273° F).

The light guide has been tested for up to 500 sterilization cycles.

Pre-Treatment

The pre-treatment must be carried out before either automatic or manual cleaning and disinfecting.

- Immediately after using (within a maximum of 2 hours), remove gross contaminations from the light guide.
- Rinse the light guide off thoroughly (at least 10 seconds) under running water or use a suitable disinfectant solution without any aldehyde (disinfectant should not contain any aldehyde to prevent blood from becoming fixed).

- Use a soft brush or a soft cloth to manually remove contaminations. Adhering polymerized composite should be removed with alcohol, a plastic spatula may help in removing the material. Do not use any sharp or pointed tools to protect the surface of the light guide from scratching.

Manual Cleaning and Disinfection of the Light Guide

- Place the light guide for the specified application time into the solution, making sure that it is completely covered (as needed, using ultrasonic support or careful brushing with a soft brush). A neutral enzymatic cleaning agent is recommended (e.g., Cidezyme/Enzol from Johnson & Johnson).
- Remove the light guide from the solution and rinse thoroughly (at least 10 seconds) in water with low germ count.
- To disinfect, place the cleaned light guide for the specified application time into the solution, making sure that it is completely covered. Disinfectants containing o-phthalaldehyde are recommended (e.g., Cidex OPA from Johnson & Johnson).
- Remove the light guide from the solution and rinse thoroughly (at least 10 seconds) in water with low germ count.
- Dry the light guide with a clean cloth.
- Check the light guide (see section "Check").

Automatic Cleaning/Disinfection (Disinfector/CDD (Cleaning and Disinfection Device))

Alternatively, cleaning and disinfecting can be conducted automatically. Information about validated procedures can be obtained from 3M Deutschland GmbH.

Sterilization

Effective cleaning and disinfection are absolutely essential requirements for effective sterilization.

Only steam sterilization is approved as a sterilization procedure:

- Maximum sterilization temperature 134° C (273° F)
- Sterilization time (exposure time at sterilization temperature) at least 20 min at 121° C (250° F) or at least 3 min at 132° C (270° F) / 134° C (273° F).

Check

Before using the light guide again, check it for damaged surfaces, discoloration, and contamination; do not use damaged light guides. If the light guide is still contaminated, repeat the cleaning and disinfection.

Clean Handpiece and Glare Shield

Clean all components with a soft cloth and, if necessary, a mild cleaning agent (e.g., dish-washing detergent). Solvents or abrasive cleaners can damage the components.

Cleaning agents must not enter the device.

- To disinfect all components, spray the disinfectant on a towel and use it to disinfect the device. Do not spray the disinfectant directly on the handpiece.
 - Disinfection agents must not enter the device!
- Dry residual disinfectants with a soft and fluff-free cloth, as they damage the plastic components.
- Make sure that disinfectants do not come into contact with the charging socket on the handpiece because this could impair proper charging.

If necessary, ask the manufacturer of the disinfectant if its constant use will damage plastic surfaces.

Clean the protection glass with a soft and fluff-free cloth. Beware of scratches!

Storage of the Handpiece during Extended Periods of Non-Use

- If the handpiece will not be used for a number of weeks – e.g., during vacation – charge the battery beforehand or connect the handpiece to the charger for this time. A safety switch within the battery prevents a total discharge.

Discharged or nearly discharged batteries must be recharged as soon as possible.

Return of Old Electric and Electronic Equipment for Disposal

1. Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources.

2. Return and Collection Systems

When your Elipar is no longer usable, do not dispose of the device with household waste. 3M Deutschland GmbH has set up special disposal facilities to handle the equipment. Details about the procedure for the specific country can be obtained from the pertinent 3M subsidiary.

3. Removing the Battery for Disposal

To dispose of the battery, disconnect the Elipar LED curing light from the charger, remove the two screws on each side of the charging socket and push the lower half of the housing towards the back, away from the upper half of the housing. Use a suitable tool to cut the connecting wires between the battery and the circuit board and remove the battery for disposal in accordance with 1. and 2.

4. Meaning of the Symbols

The EU Directive prohibits the disposal of any electric

or electronic devices marked with these symbols in combination with household waste.

Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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	Guidance and manufacturer's declaration – electromagnetic emission		
	The Elipar is intended for use in the electromagnetic environment specified below. The customer or the user of the Elipar should assure that it is used in such an environment.		
	Emissions test	Compliance	Electromagnetic environment – guidance
	RF emissions CISPR 11	Group 1	The Elipar uses RF energy only for its internal function. Therefor, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Elipar is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
	RF emissions CISPR 11	Class B	
	Harmonic emissions IEC 61000-3-2	Class A	
	Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The Elipar is intended for use in the electromagnetic environment specified below.
The customer or the user of the Elipar should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient /burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Elipar requires continued operation during power mains interruptions, it is recommended that the Elipar be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturer’s declaration – electromagnetic immunity

The Elipar is intended for use in the electromagnetic environment specified below.
The customer or the user of the Elipar should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 V	Portable and mobile RF communications equipment should be used no closer to any part of the Elipar, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [\frac{3,5}{V_1}] \sqrt{P}$
Conducted RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = [\frac{3,5}{E_1}] \sqrt{P}$ 80 MHz to 800 MHz $d = [\frac{7}{E_1}] \sqrt{P}$ 800 MHz to 2,5 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Elipar is used exceeds the applicable RF compliance level above, the Elipar should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Elipar.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Elipar

The Elipar is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Elipar can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Elipar as recommended below, according to the maximum output power of the communications equipment

	Separation distance according to frequency of transmitter m		
Rated maximum output of transmitter W	150 kHz to 80 MHz $d = [\frac{3,5}{V_1}] \sqrt{P}$	80 MHz to 800 MHz $d = [\frac{3,5}{E_1}] \sqrt{P}$	800 MHz to 2,5 GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0,01	0.12	0.40	0.77
0,1	0.37	1.26	2.42
1	1.17	4.00	7.67
10	3.69	12.65	24.24
100	11.67	40.00	76.67

For transmitters rated at a maximum output power not listed above the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

- NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

3M ESPE

3M Deutschland GmbH
Dental Products
Carl-Schurz-Str. 1
41453 Neuss – Germany

3M ESPE
Dental Products
2510 Conway Avenue
St. Paul, MN 55144-1000 U.S.A.

3M, ESPE, Elipar, Protemp and RelyX are trademarks of 3M Company or 3M Deutschland GmbH.
Used under license in Canada. © 2012, 3M. All rights reserved.

44000188035/02