CE
0051

97050551
Rev. 13
15.08

GIANO

EN
IT
FR
DE
ES
PT
RU
TR
ZH
NO
SV
FI
DA
PL
CS
HU

Cone Beam 3D Imaging
NewTom
what's next

Contents

1. INTRODUCTION AND INDICATIONS FOR USE	5
1.1. DESCRIPTION OF THE MANUAL	6
1.2. GENERAL WARNINGS	6
1.3. REQUIREMENTS (NOT PROVIDED WITH THE PRODUCT)	7
1.4. STANDARDS AND REGULATIONS	7
1.5. CLASSIFICATIONS	7
1.6. STYLISTIC CONVENTIONS	7
1.7. GENERAL SAFETY WARNINGS	8
1.7.1. INSTALLATION CONDITIONS	9
1.7.2. CONDITIONS OF USE	9
1.7.3. WARRANTY	10
1.7.4. MAINTENANCE AND DISPOSAL	10
1.7.5. CLEANING AND DISINFECTION	11
1.7.6. HYGIENE PROCEDURES FOR PATIENT PROTECTION	12
1.8. SAFETY WARNINGS	12
1.8.1. CONDITIONS OF USE	12
1.8.2. GENERAL SAFETY	12
1.8.3. SAFETY DURING X-RAY DEVICE MOVEMENTS	13
1.8.4. EMERGENCY BUTTON	13
1.8.5. CONDENSATE FORMATION	14
1.8.6. ELECTROSTATIC DISCHARGE	14
1.8.7. EXPOSURE TO LASER RADIATION	14
1.8.8. ELECTROMAGNETIC SAFETY	14
1.8.9. PROTECTION AGAINST RADIATION	16
1.8.10. APPLIED PARTS	16
1.8.11. LEAKED RADIATIONS	17
2. DESCRIPTION OF OPERATION	18
3. COMPONENTS	19
4. CONTROL PANEL	20
4.1. CONSOLE ONBOARD THE MACHINE	20
4.2. PUSHBUTTON PANEL ON TELE-X-RAY ARM	20
4.3. X-RAY EMISSION REMOTE CONTROL	21
4.4. PERFORM A SIMULATION (DUMMY RUN)	21
5. PERFORMING A 2D X-RAY EXAMINATION	22
5.1. STARTING THE SYSTEM	22
5.2. SELECTING THE EXAMINATION FROM THE CONTROL CONSOLE	22
5.2.1. 2D EXAMINATIONS AVAILABLE	22
5.2.2. SELECTING AN EXAMINATION	25
5.2.3. SETTING AN EXAMINATION FOR CHILDREN	26
5.2.4. SETTING A CURRENT EXAMINATION AS FAVOURITE	26
5.2.5. SETTING THE PROJECTION TYPE	27
5.2.6. SELECTING A REDUCED ANATOMIC REGION	28
5.2.7. CONFIGURATION OF THE X-RAY TECHNIQUE FACTORS	29
5.3. PREPARATION OF THE X-RAY EXAMINATION	30
5.3.1. DEVICES FOR PATIENT POSITIONING	30
5.3.2. SENSOR POSITIONING	31
5.3.3. PATIENT ACCESS STATUS – MINIMUM WAIT STATUS	32
5.3.4. EXAMINATION SUMMARY PAGE	33
5.3.5. EDENTOLOUS PATIENTS DEVICES	33
5.4. PATIENT POSITIONING	34
5.4.1. LASER TRACES	34
5.4.2. PATIENT POSITIONING DESCRIPTION (CRANIOSTAT)	35
5.4.3. CRANIOSTAT MOTOR-DRIVEN SUPPORT	36
5.4.4. PAN, DENT AND SIN EXAMINATIONS	36
5.4.5. TMJ EXAMINATION	38
5.4.5.1. TMJ LATERAL	38
5.4.5.2. TMJ FRONTAL	39
5.4.6. TELERADIOGRAPHIC (CEPH) EXAMINATIONS	40
5.4.7. POSITIONING FOR DTS EXAMINATION	41

5.5. ACQUISITION OF THE EXAM	42
6. 3D TOMOGRAPHIC EXAMINATION (CB3D)	43
6.1. SELECTING THE EXAMINATION FROM THE CONTROL CONSOLE	43
6.2. POSITIONING THE PATIENT FOR 3D EXAMINATIONS	46
6.3. PERFORMING THE EXAMINATION	46
7. VIEWING AND SAVING	47
7.1. VIEWING AND SAVING	47
8. CONTROL CONSOLE	48
8.1. LANGUAGE SELECTION	48
8.2. TOUCH SCREEN DISPLAY ICON	49
9. PERIODIC CHECKS AND MAINTENANCE	52
9.1. PERIODIC INSPECTIONS TO CHECK THE IMAGE	52
10. TECHNICAL DATA	52
10.1. ELECTRICAL CHARACTERISTICS	52
10.2. RADIOLOGICAL CHARACTERISTICS	53
10.3. RADIOLOGICAL CHARACTERISTICS IN CB3D MODE	54
10.4. ISODOSE CURVES FOR CB3D EXAMINATIONS	55
10.5. ISODOSE CURVES FOR 2D EXAMINATIONS	55
10.6. CTDI (COMPUTED TOMOGRAPHY DOSE INDEX) MEASUREMENTS	56
10.7. CB3D DETECTOR CHARACTERISTICS	58
10.8. PANORAMIC SENSOR CHARACTERISTICS(PAN)	58
10.9. TELERADIOGRAPHIC (CEPH) SENSOR CHARACTERISTICS (CEPH)	58
10.10. LASER CHARACTERISTICS	58
10.11. DIMENSIONAL CHARACTERISTICS	59
10.12. ENVIRONMENT CHARACTERISTICS	59
10.13. PC REQUIREMENTS	60
10.14. IDENTIFICATION LABEL POSITION	62
11. ERROR MESSAGES	63
12. USER'S LICENCE CONTRACT	65
12.1. GENERAL CONDITIONS OF THE IMAGE SOFTWARE LICENCE	65
12.1.1. USER'S LICENCE	65
12.1.2. COPYRIGHT	65
12.1.3. USE OF THE SOFTWARE PRODUCT AND EXPRESS TERMINATION CLAUSE	65
12.1.4. GUARANTEE AND EXCLUSIONS FROM THE GUARANTEE	66
12.1.5. LIMITS OF RESPONSIBILITY	66
12.1.6. APPLICABLE LAWS, JURISDICTION AND COURT OF JURISDICTION	66
12.1.7. WRITTEN FORM	66
12.1.8. TRANSLATION	66
12.1.9. SOFTWARE COMPLIANCE WITH THE REGULATIONS IN FORCE	66
13. INSPECTION AND MAINTENANCE	67
13.1. USER INSPECTION	67
13.2. TECHNICAL MAINTENANCE	68

1. INTRODUCTION AND INDICATIONS FOR USE

GIANO/VG3 is an extraoral X-ray system for digital panoramic X-rays, tele-X-rays and tomograms suitable for:

- (I) Production of orthopantomographic images of the maxillofacial region, diagnostic examination of the dentition (teeth), arches and other structures of the oral cavity.
- (II) Production of X-rays of the arches, parts of the cranium and the carpus for cephalometric examinations if equipped with a tele-X-ray arm (CEPH).
- (III) Production of tomographic images of the oral cavity and maxillofacial structures, diagnostic examination of the dentition (teeth), arches, oral cavity structures and some cranium bones if equipped with the CBCT option.

The system performs tomographic examinations with acquisition of a 360-degree rotary sequence of X-ray images and reconstruction of a three-dimensional matrix of the volume examined, producing two-dimensional views of this volume and displaying renderings of both two-dimensional and three-dimensional images. This technique is known as CBCT or CB3D.

GIANO/VG3 is a digital image reproducer suitable for expert professionals of the sector and allows easily and automatically obtaining highly detailed and accurate dental images with a low radiation dose. GIANO/VG3 is a digital system. The image is acquired by means of an X-ray detector and an X-ray source with constant potential powered by a high-frequency high-voltage generator. is a digital image reproducer suitable for expert professionals of the sector and allows easily and automatically obtaining highly detailed and accurate dental images with a low radiation dose. GIANO/VG3 is a digital system. The image is acquired by means of an X-ray detector and an X-ray source with constant potential powered by a high-frequency high-voltage generator. The image is then transferred to a computer in real-time (2D or 3D) or subsequently (2D) depending on the needs and choice of the operator.

GIANO/VG3 allows the following projections:

- Standard or paediatric panoramic views (PAN);
- Full or partial views of the dentition selected by the user (DENT);
- Frontal and lateral view of the maxillary sinuses (SIN);
- Lateral and posteroanterior view of the temporomandibular joints (TMJ) from several angles.

If the tele-X-ray arm (CEPH) is fitted, GIANO/VG3 allows the following projections:

- Cephalograms in standard or paediatric laterolateral view;
- Cephalograms in anteroposterior and posteroanterior view;
- X-ray of the hand (carpus).

If equipped with the CB3D option, GIANO/VG3 also allows tomographic images.

GIANO/VG3 is suitable for use in dentistry in the following sectors:

- Endodontics;
- Periodontology;
- Dental prostheses;
- Functional diagnostics and treatment of craniomandibular dysfunctions;
- Surgical dentistry;
- Implantation;
- Maxillofacial surgery;
- Orthodontics.

Contraindications:

- Viewing cartilaginous structures;
- The CBCT technique has a limited capability of detecting soft tissues.

1.1. DESCRIPTION OF THE MANUAL

This manual is an essential consultation tool and contains important information and instructions for use of the X-ray system and the relative controls.

These instructions describe how to properly and safely use the digital X-ray system.

Carefully read and familiarise yourself with the entire contents of the manual before attempting to use the system.

For use of the software, refer to the specific manual.

The manual is only provided in electronic format on a DVD and can be consulted directly on the PC screen during use.

It is advisable to keep a copy of this manual within reach with the aim of training the operators and as guide for consultation during use of the device. This manual also contains all the essential information for safety of the patient, the operator and the device.

It is therefore advisable to particularly carefully read the paragraphs on the safety rules.

The original text of this manual is in Italian.

1.2. GENERAL WARNINGS

The digital X-ray system and the relative software and drivers are developed and manufactured by Cefla S.C. - Imola (Italy), hereinafter referred to as the manufacturer, who is the manufacturer and distributor in compliance with the EC Medical Device Directive.

In order to use the system when it also includes the optional tomography function, a personal computer is required on which suitable software has been installed for image acquisition and storage. For software installation and use, refer to the specific software user manual. If the system is only used for the panoramic and cephalometric image acquisition functions, a personal computer can be used but is not essential. Please carefully read this manual as well as the PC and software manuals before using the device.

- No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior permission of the manufacturer.
- The manufacturer pursues a policy of continual improvement of its products, therefore, some specific instructions and images contained in this manual may differ from the product purchased.
- The manufacturer reserves the right to make any modification to this manual without prior notice.
- The information, technical specifications and illustrations contained in this publication are not binding. The Manufacturer reserves the right to make technical modifications and improvements without modifying these instructions.
- All the registered trademarks and the product names mentioned are the property of the respective owners.
- Carefully read the USER LICENSE AGREEMENT before using the product. At the time of program installation, you will explicitly be asked to accept the agreement; if you do not accept, the program cannot be installed.

Please pay particular attention to the sections in the manual where the following symbols appear:

Warnings regarding patient or operator safety.

Important information on product use.

Only for 3D devices.

In accordance with the Personal Data Protection Law in force in numerous countries, it is recommended to adequately protect sensitive data. In addition, before sending patient images or personal data via information systems, the patients' consent needs to be obtained. If required by the laws in force, dentists are obliged to protect the data using a protection password. Refer to the Microsoft® Windows operating system manual for the data access protection methods by means of password.

It is recommended to regularly (at least once a week) make a **backup copy of the databases**. This will allow restoring the data in the event of damage to the hard disk of the PC or the databases themselves.

1.3. REQUIREMENTS (NOT PROVIDED WITH THE PRODUCT)

For proper functioning, the device requires connection to a personal computer (PC) and the relative software. For the minimum requirements of the PC, refer to the paragraph TECHNICAL DATA.

The PC is not included with the device. It is recommended to only use a PC compliant to the IEC 60950-1:2007 standards for information technology devices.

1.4. STANDARDS AND REGULATIONS

The system has been designed to meet the following standards:

- Directive 93/42/EEC and s.c.
(dir. 2007/47/EC) - Medical Device Directive;
- Directive 2006/42/EEC - Machinery Directive.

Technical Standards:

IEC 60601-1:2005
IEC 60601-1-2:2007
IEC 60601-1-3:2008
IEC 60601-2-63:2012
IEC 60601-1-6:2010
IEC 62366:2007
IEC 60825-1:1993

The CE marking certifies compliance of the product as described herein with Medical Device Directive 93/42/EEC and subsequent amendments.

1.5. CLASSIFICATIONS

The system is classified as Class I and Type B as regards safety according to IEC 60601-1.

The system is classified as a Class IIB medical electrical X-ray device in accordance with Medical Device Directive 93/42/EEC and subsequent amendments.

1.6. STYLISTIC CONVENTIONS

The following symbols may be found on the X-ray device and in the manual:

Cefla S.C.– V. Bicocca 14/c – Imola (BO) ITALY	Name of the manufacturer and place of manufacture.
SN	Product serial number.
	Date of manufacture (month/year).
	Manufacturer.
	“Possible Danger: Read the user manual” symbol.

	Warnings regarding patient or operator safety.
	Consult the enclosed documentation before using the device relevant part.
	It is necessary to read the user manual before using the device.
	This symbol in the manual identifies paragraphs containing important information on use of the product.
	Applied part type B according to IEC 60601-1.
	Device in compliance with the requirements laid down in Directive 93/42/EEC and subsequent amendments.
	Ionising radiation warning symbol.
	Class 1 LASER radiation warning symbol.
	Crushing hazard.
	Symbol: dispose of the product in compliance with Directive 2002/96/EC - 2003/108/EC.
POWER	Power switch.
	Product/device identification code.
I	Unit on.
O	Unit off.

1.7. GENERAL SAFETY WARNINGS

These instructions describe how to use the system correctly. Please carefully read this manual before using the device.

The owner or the installation site manager is obliged to verify compliance with the local regulations in force and/or ask a Qualified Expert for advice. Pay particular attention to fulfilling the obligations of the law regarding the protection of workers, the population and patients against radiation.

The main REFERENCE STANDARDS are listed in this manual (1.4 - Standards and Regulations).

Do not use the system for uses different from those described in the indications for use (Introduction) and do not use it if you do not have the necessary knowledge in the dentistry and radiological field.

The law restricts and reserves the sale and use of this device exclusively to physicians, dentists or radiology specialists.

1.7.1. INSTALLATION CONDITIONS

- The system must not be used if it shows any electrical, mechanical or radiation defect. Like for all medical electrical systems, this device requires proper installation, use, maintenance and service with the aim of assuring safe and efficient operation.
- The entire system must be installed by a technician authorised by the Manufacturer under supervision of a Qualified Expert.
- The room where the system is installed must exclusively be for medical use and designed by an expert in protection against the risks associated with exposure to radiation in accordance with the regulations in force in the country of use.
- For Europe, the electric system in the room where the device is installed must be in conformity with the IEC 60364-7-710 standards (requirements for electric systems in rooms used for medical purposes).
- The X-ray system requires special precautions with respect to electromagnetic compliance and must be installed in accordance with the recommendations given in the paragraph "Electromagnetic safety" in this manual.
- The maximum dimensions reached while the unit is being handled must be taken into consideration to avoid banging into any objects present in the room. Refer to the dimensional diagram in the service manual.
- The installation must allow audio-visual communication between the operator and the patient during execution of the examination.
- The system can be installed in the following configurations:
 1. wall mounting;
 2. on the floor surface with a static baseplate (optional).
- Installation not in compliance with the instructions provided by the manufacturer might cause increased electromagnetic emission of the X-ray system and reduce its immunity to disturbances.
- In particular, use a screened cable for connection of the X-ray remote control and make the connection as specified in the technical manual.

For further details, refer to the installation template and the detailed instructions given in the service manual.

1.7.2. CONDITIONS OF USE

The device should exclusively be used by suitably trained and authorised persons (dentists and paramedics).

For installations in CANADA, make sure you comply with the following indications from Health Canada.

(a) 3D imaging should not be used for 'routine' screening (or 'screening') examinations. 3D imaging examinations must be indicated by the patient's clinical needs. Each exam must be justified by demonstrating that the benefits outweigh the risks.

(b) The use of lead aprons with thyroid collars for patients is required.

(c) Prior to imaging, women of childbearing age should be asked if they are pregnant or possibly pregnant. If so, the patient should not undergo the exam unless an accredited radiologist from a hospital setting has been consulted to discuss, with the patient and operator, the benefits and risks associated with this type of procedure along with other possible exam types.

(d) The operator should be protected by distance and shielding and should only be in the exam room close to the patient on those rare occasions when the patient requires assistance. If the operator is required to remain in the exam room, he or she must be protected with a lead apron with a thyroid collar.

In case of claims or need of technical assistance, users in Brazil are required to contact the following email address: servico.odontologico@cefla.it

1.7.3. WARRANTY

The manufacturer guarantees the safety, reliability and performance of the device.

The warranty is conditional on compliance with the following rules:

- the conditions indicated in the warranty certificate must strictly be observed;
- the device must be used exclusively according to the instructions contained in this manual;
- installation, technical service and updating of the devices must be carried out by persons authorised by the manufacturer;
- do not open the device covers. installation, repairs and in general all the operations that require opening the device must be carried out exclusively by technicians authorised by the manufacturer;
- the device must be installed exclusively in environments that comply with the restrictions given in this manual;
- the room where the X-ray unit is installed must be in compliance with the official directives that govern radiation protection in the country of use.

1.7.4. MAINTENANCE AND DISPOSAL

Never remove the device covers.

The device does not contain parts that can be repaired directly by the user. In the event of malfunctioning, do not attempt to carry out any type of maintenance operation. If you find or suspect any kind of system malfunctioning, do not attempt to carry out any type of maintenance operation and do not use the system on a patient, but directly contact your local distributor.

The user may not carry out maintenance on any mechanical or electronic part of the X-ray system.

Opening the cases to access the internal circuits may cause device breakage and failure of the electrical safety devices and will lead to forfeiture of the warranty.

Any maintenance, repairs and modifications to the device must be carried out only by personnel directly authorised by the manufacturer or by third parties expressly authorised by the manufacturer and must be carried out according to the laws in force and the generally accepted technical standards.

All the system components must be checked and replaced if necessary by qualified personnel.

For any maintenance operation, please contact the manufacturer via the website indicated on the cover of this manual and fill in the Information Request form.

Further information about device regular maintenance and inspection is provided in the document “GIANO/VG3 - Inspection and Maintenance”.

Should you for any reason need to return the device or its parts to the manufacturer or a Technical Service centre, disinfect all the external parts of the device using a specific product (see the paragraph “Cleaning and disinfection”) and preferably return it in its original packaging.

At the end of its lifetime, dispose of the device in accordance with the regulations in force. It is also advisable to disinfect all the external parts of the device before disposal and to separate the materials for differentiated waste collection. In accordance with Directives 2002/95/EC, 2002/96/EC and 2003/108/EC regarding reduced use of dangerous substances in electrical and electronic devices as well as waste disposal, do not dispose of the devices in the household waste but collect them separately. When purchasing a new device of an equivalent type, one for one, the device that has come to the end of its lifetime should be returned to the dealer for disposal. As regards reuse, recycling and other forms of recovery of the above mentioned waste, the manufacturer carries out the functions defined in the individual national legislations. Appropriate differentiated waste collection for subsequent recycling, treatment and disposal respectful of the environment contributes to preventing possible negative effects on the environment and health and encourages recycling of the materials of which the device is made up. The crossed-out bin symbol on the device indicates that the product must be collected separately from other waste at the end of its useful life. Abusive disposal of the product is liable to a fine as laid down in the individual national legislations.

1.7.5. CLEANING AND DISINFECTION

Cleaning is the first step necessary in any disinfection process. The physical action of rubbing a surface with detergents and surface-active agents and rinsing with water removes a large number of microorganisms. If a surface is not cleaned first, the disinfection process cannot be successful.

When a surface cannot be adequately cleaned, it should be protected with barriers.

The external parts of the device must be cleaned and disinfected using a product for hospital use with indications for HIV, HBV and tuberculosis (medium-level disinfectant) specifically for small surfaces.

The various drugs and chemical products used in dentist's surgeries may damage the painted surfaces and the plastic parts. The tests and research conducted have demonstrated that the surfaces cannot be fully protected against attack of all the products found on the market. It is therefore recommended to use barrier protections whenever possible.


The aggressive effects of chemical products also depend on the time they stay in contact with the surfaces. It is therefore important not to leave the product on the surfaces of the device for longer than the time indicated by the manufacturer.

It is recommended to use the specific medium-level disinfectant, STER 1 PLUS (CEFLA S.C.), which is compatible with the painted surfaces, plastic parts and unpainted metal surfaces. Alternatively, it is recommended to use products that contain:

- 96% ethanol Concentration: maximum 30 g for every 100 g of disinfectant.
- 1-propanol (N-propanol, propyl alcohol, N-propyl alcohol). Concentration: maximum 20 g for every 100 g of disinfectant.
- Combination of ethanol and propanol. Concentration: the combination of the two must be at maximum 40 g for every 100 g of disinfectant.
- **Painted surfaces and plastic parts.**
 - Incidin Spezial (Henkel Ecolab);
 - Omnidid (Omnident);
 - Plastisept (Alpro) (not tuberculocide as not an alcohol-based disinfectant);
 - RelyOn Virkosept (DuPont);
 - Green & Clean SK (Metasys) (not tuberculocide as not an alcohol-based disinfectant).

- Do not use products containing isopropyl alcohol (2-propanol, isopropanol).
- Do not use products containing sodium hypochlorite (bleach).
- Do not use products containing phenols.
- Whatever product you use, follow the instructions given by the manufacturer.
- Do not mix the STER 1 PLUS disinfectant with other products.
- Do not spray the product directly on the device surfaces.

For cleaning and disinfection use disposable, soft, non-abrasive paper (do not use recycled paper) or sterile gauze.

- It is recommended to turn off the device before cleaning and disinfecting the external parts.
- Everything used for cleaning and disinfection must be thrown away after use.

1.7.6. HYGIENE PROCEDURES FOR PATIENT PROTECTION

Single-use hygienic protections are the main means of protection against transmission of cross-infections between patients. **In order to prevent transmission of infectious diseases from patient to patient, it is essential to always use the single-use protections provided. The single-use protections are classified as Class I medical device and may not be replaced with others in a lower class.**

The single-use protections must be in compliance with the ISO 10993-1 standards on biocompatibility and approved by the control bodies where required (e.g. FDA, CE).

Always replace the single-use hygienic protections of the bite piece before positioning a new patient.

The single-use hygienic protections (Cefla code 97901337) must be stored in a dry and clean place without direct exposure to sunlight or UV rays.

The bite piece and the chinrest can be disinfected by immersing them in a cold sterilizing liquid. For sterilization of these parts, follow the instructions of the supplier of the sterilizing product.

Cover all the components that will come into contact with the hands of the dental staff with single-use protections, as they might be contaminated by indirect contact with the patient's mouth. In particular, be careful how you handle the control console of the device and the touch screen, mouse and keyboard of the PC.

Prior to positioning patients for any x-ray exposure, always cover the biteblock with a new clean non-sterile plastic barrier, to avoid cross contamination.

Note to users in Canada: ask your trusted dental material distributor for any plastic barrier that is suitable in size and is legally marketed in Canada.

According to Health Canada, bite block covers are Class I devices and are distributed by authorized establishments only, as listed in the MDEL database.

1.8. SAFETY WARNINGS

1.8.1. CONDITIONS OF USE

In order to use the device in safe conditions, refer to the following paragraphs in the manual.

1.8.2. GENERAL SAFETY

- Before leaving the surgery, deactivate the main switch of the device.
- The device is not protected against liquid penetration (Class IPX0 – common protection).
- The device is not suitable for use in the presence of a mixture of flammable anaesthetic gas with oxygen or nitrous oxide.
- Portable telecommunications devices (RF) may interfere with the X-ray device; use in the vicinity of the X-ray device should therefore be prohibited.
- The device must be looked after and kept in perfect efficiency.
- The device must never be left unattended when it is on or ready to start, in particular in the presence of minors or in general persons not authorised for its use.
- The manufacturer declines all responsibility (civil and penal) for any abuse, negligence or improper use of the device.
- Unauthorised technicians that modify the device replacing parts or components with others different from those used by the manufacturer assume the same responsibility as the manufacturer.
- Computers, monitors, printers, mice, keyboards and other devices connected to the X-ray device must be in compliance with the ISO, IEC and EN standards or the respective local regulations.
- The manufacturer is not responsible for problems or malfunctions of parts or components not approved by the manufacturer and not in compliance with the standards installed by qualified technicians acknowledged by the manufacturer.
- Do not use electronic devices that **do not comply with IEC 60601-1-2:2007** in proximity of life-support equipment (e.g. pacemakers or heart stimulators) and hearing aids. Before using any electronic device in health facilities, always check that it is compatible with the other equipment present.
- The X-ray tube contains insulating mineral oil. This oil is potentially hazardous if ingested or if it comes into contact with the skin or mucous membranes. In the event of a defect or fault, the oil may leak out. Avoid direct contact with the oil and do not inhale the vapours.
- Do not eat, drink or smoke near the device.

1.8.3. SAFETY DURING X-RAY DEVICE MOVEMENTS

The X-ray device carries out movements in the vicinity of the patient and the operator.

During the execution of the X-rays, the operator controls the movements by holding down dedicated buttons.

The reset procedure must be carried out before the patient accesses the device.

The operator must stand at a suitable distance from the moving parts. Movements can be stopped at any time by pressing the emergency button.


During all the movements of the X-ray device, the operator must:

- closely watch over the patient, and if there is a risk of collision between the X-ray device and the patient, immediately interrupt the movement by releasing the control button;
- not allow the patient to assume incorrect positions (resting the hands or other parts of the body in inappropriate areas) or to move away from the examination area.

The operator must be careful not to interfere with the sensor movement during its servo-controlled activation.

1.8.4. EMERGENCY BUTTON

The system is equipped with an emergency button to stop X-ray device operation, located under the patient support arm in proximity of the telescopic column. A remote emergency button is connected using the dedicated connector on the board located at the feet of the lifting column.

1 - Emergency button

These buttons must be activated in the event of danger and emergency, for example, failed interruption of radiation from the source, in situations of evident danger to persons or when an emergency is signalled.

If the emergency button is pressed, radiation emission is immediately interrupted and all the movements of the patient support and the rotary arm are stopped, all the linear motors run a distance of less than 10 mm and the rotary movements a distance of less than 2 degrees, and the button remains locked in safe conditions.

After the emergency, turn the button in the direction indicated by the arrow to resume normal operation.

1.8.5. CONDENSATE FORMATION

Following strong temperature oscillations, condensate may form in the X-ray device. Activate the X-ray device only after an adequate ambient temperature has been reached. See the chapter “Ambient characteristics”.

1.8.6. ELECTROSTATIC DISCHARGE

Electrostatic discharge (ESD).

Electrostatic discharge from persons may cause damage to the electronic components in the event of contact. Generally, damaged components must be replaced. The repair must be carried out by qualified technicians.

Do not touch the parts at risk marked with the symbol

1.8.7. EXPOSURE TO LASER RADIATION

The system contains some Class 1 laser diodes in compliance with IEC 60825-1. Three of them are positioned on the X-ray generator and one on the chin support base.

Both the patient and the operator may be dazzled by laser traces.

- Do not look directly into the laser beam. Be careful that the laser beam does not hit the patient's eye.
- Keep a distance of at least 10 cm between the eye and the laser.

The position of the laser sources is indicated with the following symbol

1.8.8. ELECTROMAGNETIC SAFETY

Do not use electronic devices that **do not comply with IEC 60601-1-2: 2007** in proximity of life-support equipment (e.g. pacemakers or heart stimulators) and hearing aids.

In healthcare facilities, before using any electronic device, always make sure that it is compatible with the other devices present, including the X-ray device.

Guidance and manufacturer's declaration – Electromagnetic emissions		
The GIANO/VG3 is suitable for use in the specified electromagnetic environment. The purchaser or user of the GIANO/VG3 should assure that it is used in an electromagnetic environment as described below:		
Emissions test	Compliance	Electromagnetic Environment
RF emissions CISPR 11	Group 1	This GIANO/VG3 uses RF energy only for its internal function. Therefore, the RF emission is very low and not likely to cause any interference in nearby electronic equipment. This GIANO/VG3 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic immunity			
The GIANO/VG3 is suitable for use in the specified electromagnetic environment. The purchaser or user of the GIANO/VG3 should assure that it is used in an electromagnetic environment as described below:			
Immunity Test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	IEC 60601-1-2 Test level	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst 61000-4-4	2 kV for power supply lines 1 kV for input/output lines > 3 m	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_n for 0.5 cycles 40 % U_n for 5 cycles 70 % U_n for 25 cycles 0 % U_n for 5 s	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment. If the user of the GIANO/VG3 requires continued operation during power mains interruptions, it is recommended that the GIANO/VG3 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	IEC 60601-1-2 Test level	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – Electromagnetic immunity			
The GIANO/VG3 is intended for use in the electromagnetic environment specified below. The customer or the user of the GIANO/VG3 should assure that it is used in such an environment.			
Immunity Test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic Environment
			Portable and mobile RF communications equipment should be used no closer to any part of the GIANO/VG3, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance.
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \times \sqrt{P}$ 80 MHz to 800MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5GHz
Conducted RF EN 61000-4-6	3 V 150 kHz to 80 MHz	3 V	$d = 1.2 \times \sqrt{P}$
			Where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths for fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range interference may occur in the vicinity of equipment marked with the following symbol: 

Recommended separation distances between portable and mobile RF communications equipment and the GIANO/VG3.			
The GIANO/VG3 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the GIANO/VG3 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the GIANO/VG3 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150KHz to 80MHz $d = 1.2 \times \sqrt{P}$	80KHz to 800MHz $d = 1.2 \times \sqrt{P}$	800KHz to 2.5MHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at the maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: (1) At 80MHz and 800MHz, the separation distance for the higher frequency range applies. (2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

1.8.9. PROTECTION AGAINST RADIATION

The system is an X-ray device. As such, it exposes the patient and the operators to the risks deriving from radiation. It must be used in compliance with the safety regulations set out in the radiation protection standards in force in the country of use. Some requirements are listed below:

- Start X-ray emission only from the control room. The radiation room must be adequately shielded (if required by regulations currently in force in the country of use).
- Make sure the radiation room's doors are closed before starting the examination.
- Only the patient shall be present in the radiation room during X-ray emission. If the presence of a person is necessary during the examination (for example to help patients who are not self-sufficient), personal equipment must be used to protect the individual against scattered radiation. In any case, no body parts should be exposed directly to the X-rays. Patients may not be assisted by pregnant women or minors.
- During the system initialisation procedure, no one may stay in the examination room and the doors must be closed. The room must be manned on the outside by authorised personnel until projection is complete.
- Make sure that you do not release the X-ray emission button prematurely. Note that radiation may be emitted numerous times during an exposure cycle. Wait until the exposure cycle is complete.
- The following points must always be observed:
 - During exposure, keep a distance of at least 2 metres from the X-ray source. For installations in Canada, the required distance is 3 metres.
 - Anyone not directly involved with the patient should be outside the room where the examination is carried out or stand behind a lead shield or lead glass panel during exposure.
 - Make sure that the operator can communicate verbally and visually with the patient.
 - If required, use a dosimeter for personal monitoring.

1.8.10. APPLIED PARTS

The parts that, during standard use, necessarily come into contact with the patient in order for the device to carry out its functions correctly, are: chinrest, bite piece and hygienic protections, headrest, handles, cephalon and earpiece protections.

The non-applied parts that may come into contact with patient are the external covers and the patient arm.

1.8.11. LEAKED RADIATIONS

Leaked radiation measurement is highly influenced by ambient conditions, such as wall composition and positioning, therefore, under certain circumstances, the detected values can be significantly different.

The measurement points used are at 0.5 m, 1.0 m and 2.0 m, respectively, from a central rotation axis.

The circular measurement points are determined based on the patient bite position on the machine.

Leaked radiations are measured in full panoramic mode or in 3D mode, with maximum field (13x13) for well-built patients, selected by means of a cylindrical PMMA phantom (16 cm Φ x 17 cm height), simulating the patient head.

1) Standard panoramic examination: leaked radiations measured at the maximum usage percentage allowed by the X-ray generator (corresponding to an average anode power of 42 W).

Distance between rotation axis and measurement point (Circular measurement points)	Leaked radiations *
0.5 m	40 $\mu\text{Gy} / \text{h}$
1.0 m	13 $\mu\text{Gy} / \text{h}$
2.0 m	3 $\mu\text{Gy} / \text{h}$
3.0 m	1.5 $\mu\text{Gy} / \text{h}$

2) Standard panoramic examination: leaked radiations measured at the current average usage percentage or based on 4 examinations per hour.

Distance between rotation axis and measurement point (Circular measurement points)	Leaked radiations *
0.5 m	6 $\mu\text{Gy} / \text{h}$
1.0 m	1.5 $\mu\text{Gy} / \text{h}$
2.0 m	0.5 $\mu\text{Gy} / \text{h}$
3.0 m	0.2 $\mu\text{Gy} / \text{h}$

3) 3D High Resolution 13x13: leaked radiations measured at the maximum usage percentage allowed by the X-ray generator (corresponding to an average anode power of 42 W).

Distance between rotation axis and measurement point (Circular measurement points)	Leaked radiations *
0.5 m	640 $\mu\text{Gy} / \text{h}$
1.0 m	160 $\mu\text{Gy} / \text{h}$
2.0 m	40 $\mu\text{Gy} / \text{h}$
3.0 m	20 $\mu\text{Gy} / \text{h}$

4) 3D High Resolution 13x13: leaked radiations measured at the current average usage percentage or based on 4 examinations per hour.

Distance between rotation axis and measurement point (Circular measurement points)	Leaked radiations *
0.5 m	45 $\mu\text{Gy} / \text{h}$
1.0 m	11 $\mu\text{Gy} / \text{h}$
2.0 m	3 $\mu\text{Gy} / \text{h}$
3.0 m	1.5 $\mu\text{Gy} / \text{h}$


* It is the maximum value at 15 cm above the horizontal plane of intersection with patient bite. Other values on the vertical axis are lower than these values.

2. DESCRIPTION OF OPERATION



For execution of panoramic X-rays or tomographic examinations, the X-ray device has a rotating arm coupled to a column support. The rotating arm performs motor-driven rototranslatory movements that allow moving the X-ray emission system and the image detector around the patient according to complex orbits that follow the morphologic profile. The rotary arm is fitted on a column support that can slide vertically by means of a motor-driven movement. The X-ray device position shown in the figures is the Patient Access Position.

The X-ray device can be equipped with a tele- X-ray arm coupled to the column support. The arm hosts a cephalostat to hold the patient in position during the examination and an image detector that translates synchronised with the X-ray source movement.


By choice of the user, the X-ray device can be equipped with only one image detector (the operator must then position it on the rotary arm for the panoramic examinations or on the tele- X-ray arm for teleradiographic examinations (CEPH)) or with two distinct image detectors (which are not removable, one on the rotary arm and the other on the tele-X-ray arm).

3. COMPONENTS

	Basic machine
	Panoramic image sensor
	Optional tele-X-ray arm
	Craniostat
	Optional tele-X-ray sensor (also usable for panoramic X-rays)
	CB3D detector
	Optional standard stand
	Optional "easy access" stand
	Remote X-ray button
	DVD containing the instruction manual, drivers and software for image display.
	Declaration of conformity
	Warranty certificate

4. CONTROL PANEL
4.1. CONSOLE ONBOARD THE MACHINE

1 –Touch screen


Control panel area	
Button	Description and use
	Confirm button
	Column up movement
	Column down movement
	Craniostat movement (up) + (only systems equipped with CB3D)
	Craniostat movement (down) – (only systems equipped with CB3D)
	Vertical Laser Positioning - canine or condyle cusp for patient framing (+)
	Vertical Laser Positioning - canine or condyle cusp for patient framing (-)

The control panel (console onboard the machine) allows full control of the X-ray device. The layout of the controls has been designed dividing the control panel into two areas.

The top area consists of a touch screen display that allows the user to select the settings required to start the examination, select the examination to be performed and view the status information.


The bottom area consists of 7 buttons with functions to move parts of the X-ray device (column, chinrest, laser traces) used to position the patient.

3 LEDs on the control panel indicate the system status:

	If the X-ray device is on, the green LED on the left is on. If the X-ray device is in standby, the green LED on the left flashes.
	If the central green LED is on fixed, it means that the X-ray device is ready for X-ray emission. If the remote X-ray button is pressed, the X-ray examination starts. If the central green LED flashes, it means that the device is connected to an interlock switch that signals to the system that the access door to the X-ray room is open. In this case, the device is only potentially ready, as X-ray emission will be enabled only when the door is closed. If the remote X-ray button were pressed, an error message would be shown on the display (see the chapter Error messages).
	The yellow LED on the right comes on when X-ray emission is in progress.

4.2. PUSHBUTTON PANEL ON TELE-X-RAY ARM

For units equipped with a tele-X-ray arm there is a pushbutton panel that allows vertically moving the system and pressing the CONFIRM button.

A) Upward column movement
B) Confirm button
C) Downward column movement

4.3. X-RAY EMISSION REMOTE CONTROL

The system is equipped with a remote control to enable X-ray emission.

The remote control has:

- 1) A button to confirm X-ray emission
- 2) Two LEDs
 - a. Green (device ready for emission)
 - b. Yellow (emission in progress)

When the system is ready (green LED on fixed), you can start X-ray emission by pressing the button on the remote control and holding it down for the entire duration of the examination.

When the examination is started, the yellow LED on the remote control comes on and simultaneously an acoustic signal is emitted.

Button pressed

Button released

If the button is released before the examination is complete, image acquisition will be interrupted.

If acquisition is interrupted before it is complete, the entire portion of tissue analysed will however be shown on the PC together with an error message.

If the image is exhaustive, even if partial, the examination does not need to be repeated, thus avoiding subjecting the patient to a further dose of radiation.

In order to comply with the operator safety regulations on ionising radiation (for clarifications, see the paragraph on radiation safety), the operator must stand at least 2 metres away from the X-ray device during emission.

4.4. PERFORM A SIMULATION (DUMMY RUN)

Act as follows for the device movements. The operator must:

- Closely watch over the patient, and if there is a risk of collision between the X-ray device and the patient, immediately interrupt the movement by releasing the control button.
- Not allow the patient to assume incorrect positions (resting the hands or other parts of the body in inappropriate areas) or to move away from the examination area.

To perform a simulated cycle, a simulation of the projection mode or program selected, but without emission of radiation

(dummy run), hold down the button  (CONFIRMATION).

This button must be held down for the entire cycle, otherwise the movement will stop.

This simulation can be useful for particularly emotional patients, especially children, showing them what the examination consists of without having to leave them alone in the room.

5. PERFORMING A 2D X-RAY EXAMINATION

The steps to follow to properly perform a 2D X-ray examination are:

- 1 – Turn on the system
- 2 – Select the X-ray examination on the touch screen display
- 3 – Prepare the X-ray examination
- 4 – Position the patient
- 5 – Perform the examination

5.1. STARTING THE SYSTEM

Turn on the X-ray device by pressing the on/off button  on the rear near the column base: the display lights up and an acoustic signal is emitted. Wait until the progress bar indicates that system start is complete.

If the X-ray device is in standby (green LED  flashing) press the CONFIRM button to reset the functions.

Once properly started, the green LED  stays on fixed.


5.2. SELECTING THE EXAMINATION FROM THE CONTROL CONSOLE

5.2.1. 2D EXAMINATIONS AVAILABLE


Panoramic examinations (PAN):

Preview	Name	Description
	STANDARD Scanning time: 9.3 s – 8.5 s	Projection that shows the view of the entire dental arch and the temporomandibular joints (TMJ); in the paediatric version it adapts to the morphology of the child and lasts 8.5 s. Magnification 1.23
	RH HALF-PAN Scanning time: 6.5 s	Like standard, limited to the right-hand side only. Magnification 1.23
	LH HALF-PAN Scanning time: 6.5 s	Like standard, limited to the left-hand side only. Magnification 1.23
	WIDE FOCUS Scanning time: 9.12 s	Like standard, with the focal trough enhanced in the front areas for better root exposure. Magnification 1.23

Dentition examinations (DENT):

Preview	Name	Description
	BITEWING Scanning time: 12.9 s	Set of 4 optimised images for representation of the crowns of the entire dentition. Magnification 1.43
	RH BITEWING Scanning time: 6.4 s	Set of 2 optimised images for representation of the crowns of the right dentition. Magnification 1.43
	LH BITEWING Scanning time: 6.5 s	Set of 2 optimised images for representation of the crowns of the left dentition. Magnification 1.43
	DENTITION Scanning time: 7.3 s	Projection of the entire dental arch excluding the temporomandibular joints (TMJ) and with improved orthogonality to reduce crown overlapping. Magnification 1.23
	FRONTAL Scanning time: 2.3 s	Projection of the front dentition only with enhanced focal trough to improve root exposure. Magnification 1.23
	RH HALF-ARCH Scanning time: 5.5 s	Like DENTITION, limited to the right half-arch only and with improved orthogonality to reduce crown overlapping. Magnification 1.23
	LH HALF-ARCH Scanning time: 9 s	Like DENTITION, limited to the left half-arch only and with improved orthogonality to reduce crown overlapping. Magnification 1.23

Maxillary sinus examinations (SENI):

Preview	Name	Description
	MAXILLARIES Scanning time: 6.2 s	Linear projection of the cranium in posteroanterior view on the level of the maxillary sinuses. Magnification 1.24
	LH MAXILLARIES Scanning time: 2.6 s	Linear projection of the cranium in lateral view on the level of the left maxillary sinuses only. Magnification 1.24
	RH MAXILLARIES Scanning time: 2.6 s	Linear projection of the cranium in lateral view on the level of the right maxillary sinuses only. Magnification 1.24

Temporomandibular joint examinations (ATM):


Preview	Name	Description
	DOUBLE R-L Scanning time: 11.8 s	Executes both the lateral examination (Magnification 1.30) and the frontal examination (Magnification 1.50) of both the left and the right joints; in total four images.
	FRONTAL R-L Scanning time: 5.9 s	Posteroanterior projection of both joints. Magnification 1.50
	LATERAL R-L Scanning time: 5.8 s	Lateral projection (along the major axis of the mandibular condyles) of both temporomandibular joints. Magnification 1.30
	FRONTAL R3 Scanning time: 8.8 s	Posteroanterior projection of the right joint only; executes three sections at a different angle to achieve better anatomic visibility. Magnification 1.50
	FRONTAL L3 Scanning time: 8.8 s	Posteroanterior projection of the left joint only; executes three sections at a different angle to achieve better anatomic visibility. Magnification 1.50
	LATERAL R3 Scanning time: 8.8 s	Lateral projection of the right joint only; executes three sections at a different angle to achieve better anatomic visibility. Magnification 1.30
	LATERAL L3 Scanning time: 8.8 s	Lateral projection of the left joint only; executes three sections at a different angle to achieve better anatomic visibility. Magnification 1.30

For each of the TMJ diagnostic programs, you can perform the same examination with the mouth closed, in rest conditions and with the mouth open.

Teleradiographic examinations (CEPH):

Preview	Name	Description
	ANTERO-POST Scanning time 6.6 s	Anteroposterior (AP) or posteroanterior (PA) examination of the cranium (Length 20 cm). Magnification 1.11
	SUB-CHIN VX Scanning time 10.0 s	Submentovertex examination of the cranium (length 20cm), including reverse Waters e Towne projections. Magnification 1.11
	LATERO-LAT SHORT/STANDARD/LONG Scanning time 4.6 s	Laterolateral examination of the cranium (length 18, 24 or 30cm) with automatic soft tissue improvement function, which highlights the aesthetic profile of the face. Magnification 1.11
	LATERO-LAT CHILD Scanning time 3.4 s	Laterolateral examination of the cranium in a reduced format of 19cm in height and 18cm in length, particularly suitable for paediatric examinations. Magnification 1.11
	CARPUS Scanning time 3.6 s	Scans the carpus (length 18 cm) with the aid of a special support for the hand. Magnification 1.03


5.2.2. SELECTING AN EXAMINATION

After turning on, the HOME page appears listing the examinations selected as favourites by the user (maximum six examinations).

The examinations set as favourites at the factory are:
Standard PAN, Bitewing,
TMJ Lateral R-L, TMJ Frontal R-L,
Standard CEPH Lateral, CB3D.

If the examination to be performed is not among the favourites, tap the EXAM icon to access the complete list of examinations available.


Selecting one of the favourite examinations by tapping the relative icon (if present), sets that examination as CURRENT examination and displays the dedicated EXAMINATION page.

Select the category of the examination to be performed from the various types: PAN - DENT - SIN - TMJ - CEPH - 3D.

The category selected is highlighted (in the example on the side it is PAN).


After this operation, the display automatically shows the EXAMINATION page of the current examination category selected.

5.2.3. SETTING AN EXAMINATION FOR CHILDREN

To activate the examination parameters adapted to the build of a child, tap the CHILD icon , if available. The icon is highlighted  if the settings for children are active and goes back to the initial state when they are deactivated by tapping on the icon again.

5.2.4. SETTING A CURRENT EXAMINATION AS FAVOURITE

To set the CURRENT examination as favourite, tap the icon .

The FAVOURITES home page will be displayed where you can insert the selected examination by tapping on a free position. If you select an already occupied position, the examination previously set as favourite is overwritten.

To go back to the EXAMINATION page tap the icon . If the CURRENT examination is already among the favourites, the relative icon is red .


5.2.5. SETTING THE PROJECTION TYPE

In some examination categories, different projection types are available.

In the PAN and TMJ categories, tapping the icon , you can switch between one type of projection and the next.


PAN

For the panoramic examinations you can choose between STANDARD and WIDE FOCUS.


DENT

For the dentition examinations you can choose between STANDARD and BITEWING.

TMJ

For the TMJ examinations you can choose from 3 projections: LATERAL, FRONTAL, DOUBLE R-L.


For the TMJ examinations you can choose between Mouth Closed  and Mouth Open  for each projection.


CEPH

For the teleradiographic examinations you can choose between LATERAL and FRONTAL projections.

The icon changes status as shown .

CEPH FRONTAL

CEPH LATERAL

For ceph lateral examinations, it is possible to enable the shoulder anti-collision option if the patient anatomy requires it.
The icon changes status as shown:

(option enabled)


(option disabled)

 For a CEPH examination, turning the cephalostat from frontal to lateral position and vice versa also determines the setting shown on the display of the control console.


5.2.6. SELECTING A REDUCED ANATOMIC REGION


For the X-ray examinations that permit it, you can choose from the anatomic regions available by tapping the boxes that appear on the graphic representation of the anatomic region involved in the examination.

STANDARD PAN examination


HALF PAN RIGHT examination


HALF PAN LEFT examination

5.2.7. CONFIGURATION OF THE X-RAY TECHNIQUE FACTORS


On the EXAMINATION page, tapping the icon  you can access CONFIGURATION of the X-ray technique factors (kV, mA, and seconds of exposure) according to the patient's build.

Configuration page.


AUTOMATIC: the X-ray technique factor configuration ensures the best quality of the final image, while minimising the amount of X-rays received by the patient.

PRESET: sets the parameters following the default values for man, woman and child. For each of the adults presets, three different builds are available:

big  *average*  *small* 

Tap the icon corresponding to the adult/woman/child default setting according to the patient's build.

CUSTOM: sets the parameters manually with an accuracy of 1 kV or 1 mA.

Select CUSTOM only if you decide to manually set the exposure parameters.

To increase or decrease the kV and mA values, tap on the icons  .

5.3. PREPARATION OF THE X-RAY EXAMINATION

5.3.1. DEVICES FOR PATIENT POSITIONING

Examination type	Movable support	Image
PAN DENT SIN	Craniostat, chinrest and bite piece	
TMJ	Subnasal support and craniostat	
CEPH	Tele-X-ray cephalostat	

Remember to change the single-use protections before positioning each new patient.

Before each X-ray examination, **make sure that the patient has taken all metal objects off**, such as glasses, removable prostheses, earrings and other removable metal objects at the level of the head or the neck. If a radiation protection apron is used, make sure that the patient's neck is not covered otherwise there would be an area not exposed to the X-ray.

5.3.2. SENSOR POSITIONING

Check that the sensor to be used is fitted in a suitable position for the examination to be performed, otherwise relocate it.

If the sensor is not in a suitable position for the examination to be performed, a warning will appear on the control console onboard the machine and you will not be able to continue with the examination selected.

Normally, the sensors for the PAN and CEPH examinations can not be removed by the user.

Only if the X-ray device is equipped with the CEPH arm for tele-X-rays, but with only one sensor, the sensor needs to be moved from the CEPH to the PAN position and vice versa according to the type of examination you want to perform. The system recognises automatically if the sensor is attached to the correct position according to the type of examination planned. If the sensor is not in the correct position for the examination, the system allows you to move it to the correct position.

The sensor locking system contains both electronic parts and a mechanical lock.


During operation, the sensor is mechanically locked and cannot be removed.

Do not attempt to remove the sensor if this function is not requested and set up by the X-ray device. The sensor is a delicate electronic part. If you try removing the sensor by forcing the holder, there is a risk of damaging both the sensor and the locking system.

The X-ray device automatically sets up for removal or fitting when you select an examination type different from the current one and press the CONFIRM button on the control console.

If the sensor is not enabled (because it is the wrong one or absent), an error message appears on the console display and you cannot perform the examination.

The sensor automatically moves to the release position when necessary. The position can easily be detected by looking if the sensor is free in the upper area (A).

On units equipped with two sensors, be particularly careful to fit each sensor on the correct holder. The two sensors can be recognised by the height of the black band that identifies the area sensitive to X-rays (see images).

PAN sensor


CEPH sensor

5.3.3. PATIENT ACCESS STATUS – MINIMUM WAIT STATUS

When the X-ray device is ready, it may be in two distinct and consecutive conditions which you access by pressing the CONFIRM button  :

Patient Access Status = the X-ray device sets up to allow patient access and correct cranium positioning.

Minimum Wait Status = this is obtained when pressing the CONFIRM button again after positioning the patient; this condition anticipates some preparatory movements of the device, reducing to a minimum the time the patient remains alone during X-ray emission and hence not assisted by the operator.

Both conditions are indicated on the control console by the green LED being on (or flashing if an INTERLOCK connection is open).

When pressing the CONFIRM button, wait for a few seconds for the X-ray device to complete its movements; in this period you will not be able to use the control console, the only exception being the CONFIRM button to interrupt the process.

In both the Patient Access Status and the Minimum Wait Status condition, the Examination Summary page will be shown on the control console display.

In “Patient Access and Minimum Wait” statuses, during patient access and, in any case, before starting the

examination, ensure that the X-ray device is not hit unintentionally: if this is the case, press  in the

Examination Summary page and position the X-ray device again with .

5.3.4. EXAMINATION SUMMARY PAGE

This page can be viewed only if the device is ready. This page shows:

- the various X-ray technique factors correctly set and an icon for the relevant setting mode  (AUTOMATIC),
 (PRESET) or  (CUSTOM);
- the examination type selected;
- icons at the bottom of the page with which you can redefine the projections and the anatomic region involved in the examination.

5.3.5. EDENTULOUS PATIENTS DEVICES

In case of edentulous patients, for the PAN, DENT, SIN, DTS or CB3D examinations requiring the use of the bite piece, apply the single-use soft insert supplied, as follows:

- 1 - remove the soft bite from the mould (Fig. 1);
- 2 - apply the single-use protection on the bite piece (Fig. 2);
- 3 - insert the soft bite in the bite piece as indicated in the figure (Fig. 3);
- 4 - push the bite until the material block is ejected (Fig. 4).

Always apply the single-use protection on the bite piece before inserting the soft bite.

Then proceed to standard positioning for the exam, asking the patient to tighten the mucous membranes at the centre of the soft insert.


Fig. 1

Fig. 2

Fig. 3

Fig. 4

5.4. PATIENT POSITIONING


Set the X-ray device to Patient Access Status before having the patient access and before starting any cranium positioning operation.

If the X-ray device is not in Patient Access Status, press the CONFIRM button once, wait for the X-ray device to complete its movements and for the laser traces to turn on.

If you want to interrupt device movement, press the CONFIRM button again.

5.4.1. LASER TRACES

The X-ray device is equipped with four illuminated traces to help position the patient:

1)Upper Horizontal Laser Trace	used in all PAN, DENT, SIN, ATM, 3D examinations
2)Sagittal Vertical Laser Trace	used in all examinations
3)Focusing Vertical Laser Trace	used in PAN, DENT, SIN, TMJ examinations
4)Lower Horizontal Laser Trace	used in CB3D, CEPH examinations

Upper horizontal trace [1]

This is generated by a laser projector located on the side of the X-ray generator; it can be adjusted up or down to adapt it to different-sized heads by acting on the lever situated next to the light beam opening. This trace is used in PAN and DENT examinations to ensure that the Frankfurt plane of the patient is horizontal and, more in general, correct positioning.

The Frankfurt plane is represented by an imaginary line that runs from the upper edge of the auditory meatus to the lower limit of the orbit.

Sagittal vertical trace [2]

This ensures symmetry of the patient's head with respect to the sagittal median line. Observing this trace, make sure that the patient looks straight ahead of himself or herself to avoid that his or her head is tilted sideways or slightly turned.

Focusing vertical trace [3]

This indicates the exact position of the focal trough; for good focusing in PAN and DENT examinations, this trace must coincide with the upper canine cusp. The same trace is used to aim at the exact position of the condyle head for TMJ examinations.

Lower horizontal trace [4]

This has the dual purpose of tracing the lower limit of the field of view during CB3D examinations or is used in teleradiographic (CEPH) examinations to ensure that the Frankfurt plane of the patient is horizontal.


The illuminated traces appear when you press the CONFIRM button once and the X-ray device has completed its movements to go into the patient access position. After 30 seconds of inactivity, the traces disappear; to make them

reappear, just press one of the buttons , **BUT NOT THE CONFIRM BUTTON.**

Each time you press the CONFIRM button, all the movable parts of the device move. Be careful NOT to press the button during patient positioning and make sure that the device can move in safe conditions.

5.4.2. PATIENT POSITIONING DESCRIPTION (CRANIOSTAT)


Craniostat components:

- 1 – Base
- 2 – Rods
- 3 – Chinrest
- 4 – Bite piece
- 5 – Bite piece block lever
- 6 – Cross member
- 7 – Front support
- 8 – Arms
- 9 – Anatomic arches
- 10 – Arms block levers
- 11 – Undenose support


The craniostat consists of a lower part and an upper one, connected by two carbon rods (2).


The lower part consists of a base (1) which is fitted in the chinrest seat with metal pins and that can be removed by simply pulling upwards.

The chinrest (3) is fitted in the dedicated inserts in the base (1) by means of pins and can be removed simply by sliding it upwards.


The bite piece (4) is fitted in the suitable hole of seat (1): once it is set at the desired height, pull the central lever (5) from the left (unlocking position ) to the right (locking position ) in order to centre it and lock it in place.

Please remember that, in order to remove it, it is necessary to take lever (5) to the unlocking position and to slide it out of its seat.

The upper part consists of a cross member (6) that can slide vertically along the carbon rods (2).

The cross member also features two arms (8), installed on its sides, that can be extended or closed depending on the size of patient's cranium. The anatomic arches (9) are secured at the arm ends: they can rotate freely on their pin to adjust to patient's anatomy.

After having correctly positioned the patient's head,:

- a - push front support (7), making the pad rest against the patient's forehead,.
- b - rotate the arms (8) towards the patient's temples, ensuring that the rubber blocks at the ends of arches (9) adhere to the head anatomy.
- c - rotate levers (10) downwards to ensure a correct locking.

At the end of the exposure, rotate the levers upwards to unlock arms, thus allowing the patient to get out easily

In examinations where the undernose support (11) is used, it should be fitted **INSTEAD OF THE CHINREST** in the inserts in the base and pushed **FULLY DOWN INTO PLACE**.

5.4.3. CRANIOSTAT MOTOR-DRIVEN SUPPORT

The chinrest is movable only for 3D examinations. Use the   buttons located on the left-hand side of the control console to raise or lower the chinrest to suit the patient.

Briefly pressing the button, moves the mechanism 1 mm up or down, while holding it down moves it continuously at constant speed.

5.4.4. PAN, DENT AND SIN EXAMINATIONS


- 1) Adjust the height of the unit to facilitate patient access using the column up or down movement buttons. The telescopic column moves slowly at first and then picks up speed. Adjust the height until the bite piece is slightly higher than the occlusal plane of the patient. That way you induce the patient to stretch to reach the bite piece helping him or her to stretch and straighten the neck.

- 2) Make sure that the bite piece is turned sideways to leave room for the patient. Fit the single-use bite piece protection.
- 3) Guide the patient towards the unit so that he or she is in front of the bite piece and can grip the large handles. The operator and the patient face each other. Ask the patient to rest his or her chin on the chinrest.




- 4) Ask the patient to move forward one step still gripping the handles until reaching the position shown in the figure.

- 5) Adjust the height of the bite piece and turn it inside the mouth letting the patient bite on it as shown in the figure. The tip of the upper and lower incisors must sit in the groove of the bite piece. The interproximal space of the incisors must be on the median line of the bite piece.

The support abutment can slide up and down facilitating proper positioning of the bite piece. Then tighten the knob to lock the bite piece in the correct position.

If you are performing a BITEWINGS examination, you **NEED TO ADJUST THE BITE PIECE TO A PRECISE HEIGHT**. It is advisable to remove the chinrest.

Then adjust the bite piece vertically, in the position where a clack is heard and a slight resistance to vertical sliding is perceived. After having reached this position, lock the bite piece as described in par. 5.4.2, by turning the lever to the locking position .

6) Check the symmetry of the patient's head using the sagittal vertical trace as guide.

7) Check proper positioning of the Frankfurt plane superimposing the upper horizontal light trace. To adjust the patient's head inclination, act on the column up and down movement buttons. Make sure that the patient keeps his or her back straight and extended.

8) After having reached the correct positioning, lock the craniostat as explained in paragraph 5.4.2:

- A- Push the front support forward until resting it against the patient's forehead. The forehead pressure on the support stabilizes it automatically.
- B- Close the arms until they adhere to the patient's cranium.
- C- Rotate the upper levers downwards until locking the patient laterally.

9) Ask the patient to smile to uncover the upper teeth. Usually the vertical light trace falls between the canine cusp and the distal surface of the upper premolar*.

In case of patient dysmorphia, approach the light trace to the canine, by using the console buttons  in order to optimise dentition focusing.

* the canine can be used as a useful reference to optimize the patient alignment, but this is not strictly necessary.

10) Press the CONFIRM button and just before leaving the room to press the X-ray emission button, ask the patient to close his or her eyes and to swallow.



5.4.5. TMJ EXAMINATION

5.4.5.1. TMJ LATERAL

1) Remove the chinrest and the bite piece and fit the subnasal support.

2) Adjust the height of the unit to facilitate patient access using the column up or down movement buttons

until the subnasal support is at the height of the base of the nose. The telescopic column moves slowly at first and then picks up speed.

3) Guide the patient towards the unit so that he or she is in front of the subnasal support and can grip the large handles. The operator and the patient face each other. Let the patient rest the base of the nose on the subnasal support as shown in the figure..

4) Check the symmetry of the patient's head using the sagittal vertical light trace as guide. Check proper positioning of the Frankfurt plane superimposing the upper horizontal light trace as shown in the figure. If the examination requires it and it is necessary, slightly tilt the patient's head forward to help him/her open the mouth as wide as possible.

5) After having reached the correct positioning, lock the craniostat as explained in paragraph 5.4.2:

D- *Push the front support forward until resting it against the patient's forehead. The forehead pressure on the support stabilizes it automatically.

E- Close the arms until they adhere to the patient's cranium.

F- Rotate the upper levers downwards until locking the patient laterally.

*not essential

6) Check that the required examination is correctly

selected by observing the MOUTH CLOSED

or

MOUTH OPEN icon


Note: in the case of an examination with the mouth open, the light trace moves forward: opening the mouth wide, the patient's condyle moves out of the mandibular fossa and forward.

Then act on the buttons


to position the focusing vertical light trace exactly on the head of the condyle as shown in the figure.

- 7) Press the CONFIRM button and just before leaving the room to press the X-ray emission button, ask the patient to close his or her eyes and remain immobile.

5.4.5.2. TMJ FRONTAL

Follow the steps as for the TMJ LATERAL examination but acting differently in steps 4 and 6 as follows:

Position the patient's head not according to the Frankfurt plane but so that the horizontal light trace passes from the upper margin of the orbit to the upper margin of the auditory meatus in such a way that these two points are on the same horizontal plane, as shown in the figure.

To adjust the head inclination, just raise or lower the unit using the height adjustment buttons.

5.4.6. TELERADIOGRAPHIC (CEPH) EXAMINATIONS

Teleradiographic examinations can be performed only if the system is equipped with a tele-X-ray arm with relative cephalostat. For these examinations the patient generally remains standing. For very tall or very short patients or in a wheelchair, the examinations may be performed with the patient sitting down.

If a chair is used, make sure that the backrest or armrests do not obstruct proper movement of the device.

- 1) Open the earpiece supports wide apart acting on the upper white sections and not on the transparent rods. Fit the single-use earpiece protections.
- 2) Turn the NASION support upward.
- 3) Adjust the height of the motor-driven column using the dedicated buttons until the earpiece rings are at the height of the external auditory canal of the patient.
- 4) Have the patient access the cephalostat. The patient must stand with his or her back straight and staring at the horizon.
- 5) Close the earpieces in such a way that they slip slightly into the auditory meatuses being careful not to cause discomfort.
- 6) For LATEROLATERAL projections, position the patient's head so that the Frankfurt plane coincides with the horizontal light trace. For ANTEROPOSTERIOR, POSTEROANTERIOR, SUBCHIN-VERTEX, REVERSE WATERS AND TOWNE projections, position the patient's head according to the shooting angle required for the method selected.
- 7) Turn the NASION support downward and adjust the depth and height so that it rests in correspondence to the Nasion point of the patient without pushing and without changing the position previously set.
- 8) For ceph lateral examinations, it is possible to enable the shoulder anti-collision option if the patient anatomy requires it.
The icon changes status as shown:
 (option enabled)
 (option disabled)
- 9) Press the CONFIRM button and just before leaving the room to press the X-ray emission button, ask the patient to close his or her eyes, clench the teeth and keep the lips relaxed.

5.4.7. POSITIONING FOR DTS EXAMINATION

Carefully read the instructions provided at the beginning of paragraphs 5.4 and 5.4.1.

The positioning devices to be used for this family of examinations are those specific for PAN/DENT/SIN, as indicated in paragraph 5.3.1.

For patient positioning for PAN examinations, follow the instructions provided in paragraph 5.4.4.

In case the volume reconstruction involves the rear side of the mandible (see figure 1), the patient must be positioned as shown in figure 3: The horizontal laser must pass through the auditory meatus and the base of the nose so that it is parallel with the mandible plane, whereas the vertical laser must always be positioned close to the upper canine.

For the reconstruction of front areas and all maxillary sinuses (see figure 2), the positioning is the same as for the standard panoramic examination (figure 4), therefore the horizontal light trace must pass through the auditory meatus and the lower margin of the orbit, and the vertical laser must be close to the upper canine.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

For further information on how to carry out exposure correctly, refer to document "DOT User Manual".

5.5. ACQUISITION OF THE EXAM

- Visually check that the patient is correctly positioned and that the central green LED on the X-ray remote control is on fixed.
- Optionally, press the CONFIRM button to go to the Minimum Wait Status position.
- Ask the patient not to move during the examination and to breathe slowly and regularly. For panoramic (PAN) projections, ask the patient to swallow just before taking the X-ray so that the tongue is up against palate.
- Have anyone not necessary move away from the area exposed to radiation and where required stand behind the special protection.

- Press the X-ray button on the remote control and hold it down for the entire duration of the examination. The examination time is determined by the intermittent yellow LED on the X-ray remote control device. X-ray emission is signalled with an acoustic signal.

The X-ray device emits rays to acquire the images only if it is ready, i.e. when the green LED on the control console and on the remote control is on. It may occur that as a result of a user or X-ray device error, the device does not indicate that it is ready and X-rays can hence not be emitted. Correct the error (see the chapter Error Messages) and press the CONFIRM button.

6. 3D TOMOGRAPHIC EXAMINATION (CB3D)

ONLY FOR 3D DEVICES

The CB3D examination is obtained by three-dimensional reconstruction of the anatomic region X-rayed and can be consulted both through two-dimensional views and three-dimensional representations generated by a program run on a workstation (PC).

Read the NNT software user manual for the instructions for image processing.

Remember to change the single-use protections before positioning each new patient.

Before positioning the patient, make sure that the patient has taken off all metal objects, such as glasses, removable prostheses, earrings and other removable metal objects at the height of the head. If a radiation protection apron is used, make sure that the patient's neck is not covered, otherwise there would be an area not exposed to the X-ray.

6.1. SELECTING THE EXAMINATION FROM THE CONTROL CONSOLE

The 3D examination can only be performed if the X-ray device is connected to a PC with the NNT software installed.


- 1) Turn on the PC and run the NNT program .
- 2) At the first start of the day, you are asked to complete the daily check procedure before being able to perform any 3D examination. The daily check procedure is a service procedure with emission of X-rays to be carried out WITHOUT A PATIENT.
- 3) For the instructions, see "Acquisition operations" in the NNT manual.


The patient does NOT need to access the X-ray device during the examination selection operations, but only at the end of the procedures described in this paragraph.


View the FAVOURITES home page by tapping the HOME icon .

If the 3D examination is among the favourites, you can select it by tapping the relative icon, otherwise access the 3D section after pressing EXAM.

In both cases, the page shown here on the side will be displayed showing the correctly selected CB3D program.


1) Selecting the field of view (FOV)

Set the field of view by tapping the icon .

The SET FOV page appears: select the field of view you want to use. The fields of view available vary depending on the model of the X-ray device installed.


The fields of view available vary depending on the model of the X-ray device installed and on the software licenses purchased.

In particular, with the Extra FOV software licence more fields of view are available with respect to the basic licence. A field of the Extra FOV type envisages the automatic execution of a double scan that provides higher volumes with respect to those that can be achieved with standard fields of view. On the control console, the additional fields of view feature the letter "e" next to the indication of their size.



2) Setting the protocol

Select the examination protocol by tapping the icon .

The SET ARCH page appears. Depending on the machine equipment the following choices will be available or not:

3) Selecting the mode

Select the mode with which you want to perform the examination:


High resolution; total scanning time ~18s, minimum voxel size of 150µm.

Maximum resolution; total scanning time ~36s, with minimum voxel size of 75µm.


5x5


8x5


8x8


11x5

11x8

11x13

4) Selecting the centre of volume reconstruction

Select the centre of CB3D reconstruction from the EXAMINATION page by tapping on the anatomic area you want to reconstruct. A red circle highlights the area selected.

The selection of the centre of reconstruction is only preliminary; after positioning the patient, fine adjustment of the anatomic region involved in the 3D examination will be made from the PC workstation through a guided procedure. Note: it is not required nor permitted to modify the X-ray exposure technique factors (kV, mA, seconds of exposure) as they are automatically optimised by the system during execution of the examination.


You can now have the patient access as described in the next paragraph.

6.2. POSITIONING THE PATIENT FOR 3D EXAMINATIONS

1) Once you have completed preparation of the X-ray device, have the patient access.

2) Adjust the height of the motor-driven column using the buttons   to facilitate patient access. Bring the column up to the height of the patient.

3) The patient should grip the handles with both hands and be in an erect position.


4) If necessary, adjust the height of the chinrest using the

buttons at the bottom of the control panel   and centre the acquisition area marked by the lower horizontal laser trace, which determines the lower limit of the area radiated.

5) Arrange the patient's head using as guide the vertical laser trace that identifies the sagittal plane.

6) After having reached the correct positioning, lock the craniostat as explained in paragraph 5.4.2:

- G- Push the front support forward until resting it against the patient's forehead. The forehead pressure on the support stabilizes it automatically.
- H- Close the arms until they adhere to the patient's cranium.
- I- Rotate the upper levers downwards until locking the patient laterally.

7) Go to the PC workstation.

6.3. PERFORMING THE EXAMINATION

Follow the instructions for the acquisition operations described in the NNT manual to complete image acquisition. Refer to "eFOV" mode for the acquisition with fields of the Extra FOV type.

7. VIEWING AND SAVING

In order to view and save the examination, you need a PC and dedicated software.

The X-ray system is supplied with the NNT program for viewing and saving the examinations. If you have this software, refer to the NNT user manual.

If you use third-party programs for viewing and saving the examinations, refer to the instructions provided with the software application used.

Use of the NNT software is optional in the case of 2D examinations (e.g. panoramic and cephalometric).

However, NNT must be used to acquire tomographic examinations as it contains the image volume reconstruction technology.

If you need to hand the X-ray examination to the patient or another operator, NNT automatically guides you in creating a DVD that will include a redistributable copy of NNT for image viewing (NNT Viewer).

Alternatively, you can export only the X-ray images in a standard format (DICOM 3.0) so that they can be viewed with third-party programs.

7.1. VIEWING AND SAVING

The last 2D projection acquired with the device remains stored in its internal memory until the device is turned off or the projection is overwritten when performing a new examination.

If an image is available in the internal memory, the icon  is active on the home page  and the image can be saved to a USB key.

To save an image:

- 1) check that the icon  is active (if not, it means that no image is available in the internal memory);
- 2) insert the USB key in the connector at the base of the machine console;
- 3) press the button  on the home page.

The saving procedure is started and a progress bar indicates the status.

NOTE: DO NOT REMOVE THE KEY UNTIL WRITING IS COMPLETE.

When complete, the console automatically returns to the home page.

The following errors may occur while saving the image:

Err.0.31: "USB key not found" => the key has not been inserted properly or has not been recognised by the device. Try reinserting the key and repeat the procedure. If the error persists, change the key with a compatible one (see the list at the end of this chapter).

Err.0.32: "USB error on writing" => the key is corrupted, incompatible, write-protected, or there is insufficient free space. Retry the saving procedure and if the error persists, check that the key is not write-protected and that there is sufficient space, then retry. If necessary, replace the key with a compatible one.

Err.0.31: “USB key not found” => the key has not been inserted properly or has not been recognised by the device. Try reinserting the key and repeat the procedure. If the error persists, change the key with a compatible one (see the list at the end of this chapter).

Err.0.32: “USB error on writing” => the key is corrupted, incompatible, write-protected, or there is insufficient free space. Retry the saving procedure and if the error persists, check that the key is not write-protected and that there is sufficient space, then retry. If necessary, replace the key with a compatible one.

A PNG image will be saved on the USB key, which needs to be dragged to NNT in order to be viewed correctly.

Listed below are the USB keys that have been tested for compatibility with the device:

- Sandisk Cruzer 4GB,
- Sandisk Cruzer 8GB,
- Sandisk Cruzer 16GB,
- Kingston Traveler 16GB,
- Kingston Traveler 32GB.

8. CONTROL CONSOLE

The upper part of the console onboard the machine consists of a colour touch screen. The user interacts with the display by tapping graphic icons with selection button functions. The intuitive operator interface makes it easy for the operator to select the settings necessary to start the examination.

8.1. LANGUAGE SELECTION

When the system is turned on, the FAVOURITES home page is displayed.

If a different page is displayed, go to the HOME page by tapping the button





Tap the icon to access the CONSOLE CONFIGURATION menu.

Tap the icon to access LANGUAGE selection where you can select a language from those available by tapping the icon depicting the country.

8.2. TOUCH SCREEN DISPLAY ICON

Button	Description and use	Button	Description and use
FAVOURITES home page			
	HOME button		
CONSOLE CONFIGURATION menu			
	Tapping this icon you access the CONSOLE CONFIGURATION menu		Selects the current language
	Resets the list of favourites		Selects the touch screen configuration menu
	Calibration button With a fingertip tap in the centre of the blue circles displayed until calibration is complete		Touch screen button sensitivity Increase (+): higher sensitivity Decrease (-): lower sensitivity (you need to tap harder to activate the button)
	General information on the X-ray device		CANCEL: Cancels the action performed The display shows the previous page.
CURRENT EXAMINATION page icon			
	This icon activates manual examination selection mode where you scroll a list using two arrow buttons instead of the PAN, DENT, TMJ, SIN, CEPH and 3D icons, and it allows you to select the diagnostic programs without acting on the graphic icons You can also use it as alternative method to select from the possible fields of view		NB: Use of the arrows is always optional You can access all the machine functions by tapping the relative icons with a fingertip
	CHILD BUILD button Tapping this button, you can switch between a generic adult build and a child build 	 	FAVOURITES button Allows inserting the current examination in the favourites list; after tapping this icon, you will be asked where in the list you want to insert the examination If the current examination is already listed among the favourites, the icon is highlighted and you can remove the examination from the favourites by tapping the highlighted icon

Button	Description and use	Button	Description and use
	<p>X-RAY TECHNIQUE FACTOR CONFIGURATION page</p> <p>This icon allows access to configuration of the X-ray technique factors: the parameters KV, mA and seconds of exposure</p> <p>PRESET: sets the parameters following the default values for man, woman and child. For adults you can select three different builds</p> <p>high  medium  low </p> <p>CUSTOM: manually sets the parameters with an accuracy of 1 kV or 1 mA</p>		
Technique Factor Setting			
	Factor setting from the PRESET page		PRESET factors icon on the EXAMINATION SUMMARY page tap to modify the factors
	Increase or decrease of the technique factors		CUSTOM factors icon on the EXAMINATION SUMMARY page tap to modify the factors
	BACK button to go back to the previous page		
Only for TMJ examinations			
	<p>Mouth Open</p> <p>Selects the examination mode with the patient's mouth open</p>		<p>Mouth Closed</p> <p>Selects the examination mode with the patient's mouth closed or in rest position</p>
	<p>Mouth Open action</p> <p>This icon appears on the Examination Summary page and allows redefining Mouth Open or Mouth Closed mode.</p>		<p>Mouth Closed action</p> <p>This icon appears on the Examination Summary page and allows redefining Mouth Open or Mouth Closed mode</p>
Only for CEPH examinations			
	<p>Ceph Frontal</p> <p>Tap to select frontal teleradiographic examination mode</p>		<p>Ceph Lateral</p> <p>Tap to select lateral teleradiographic examination mode</p>
	<p>Shoulder anti-collision option disabled (default), for ceph lateral examinations</p> <p>Press to enable the patient shoulder anti-collision option</p>		<p>Shoulder anti-collision option enabled for ceph lateral examinations</p> <p>Press to disable the patient shoulder anti-collision option</p>

Button	Description and use	Button	Description and use
Only for 3D examinations			
	High-resolution acquisition mode CB3D		Maximum resolution CB3D acquisition mode
	FOV Tapping this button you access the SET-FOV menu where you can select the various dimensions of the scanning area		Back to FOV Icon on the Examination Summary page that allows going back to the SET-FOV menu
	Selecting both arches		Selection of the temporomandibular joint (TMJ) in CB3D mode
	Selection of the lower arch		Selection of the maxillary sinuses in CB3D mode
	Selection of the upper arch		

9. PERIODIC CHECKS AND MAINTENANCE

In the interest of safety and health of the patient, the staff or third parties, inspections and maintenance need to be carried out at set intervals.

Period	Operator	Object	Description
Yearly	Specialised technician of the dealer that initially installed the device or another technician authorised by the manufacturer	The X-ray device as a whole	In order to assure the operating safety of the device, it is advisable to inspect the x-ray device in all its parts, in order to prevent or repair any faults

For installations in the U.S.: please refer to the document “Addendum: Inspection and Maintenance”.

9.1. PERIODIC INSPECTIONS TO CHECK THE IMAGE

- Blank image acquisition

Acquisition of the blank image allows optimising scanning performance. This brief procedure is obligatory and requested by the NNT software every two weeks. To carry out the procedure, see “Acquisition operations” in the NNT manual.

- Quality control by means of a phantom (QA Phantom)

Image resolution (sets of lines/mm)	2,5 lp/mm
Low contrast resolution	2 low contrast points at least can be distinguished
Artefacts	The image must not feature artefacts such as visible horizontal lines or halftoning
Frequency of check	Yearly

Use the optional phantom, available upon request, with the relevant user instructions, or, as an alternative, use a dental phantom compliant with IEC 61223-3-4:2000 standard.

The quality control consists of performing a normal examination on QA Phantom through an automatic procedure guided by the NNT software. This periodic quality control guarantees proper functioning of the device and the validity of the results obtained.

The control is carried out by a qualified expert or by the dealer who carried out the first installation of the X-ray device or by another technician authorised by the manufacturer.



To correctly carry out this procedure, see “Acquisition operations” in the NNT manual and the technical manual.

10. TECHNICAL DATA


10.1.ELECTRICAL CHARACTERISTICS

Power supply voltage	115 – 240 VAC +/-10% Single phase
Power supply frequency	50 / 60 Hz
Current absorbed in rest conditions	1A @ 115V; 0.5A @ 240V
Maximum current absorbed in operating conditions	14A @ 115V; 7A @ 240V
Column movement operating cycle	25 s ON, 400 s OFF
Maximum apparent line resistance	0.5 Ω @ 240 V - 0.25 Ω @ 115 V
Overvoltage type	II

10.2.RADIOLOGICAL CHARACTERISTICS

Generator voltage	60 – 85 kV continuous mode (2D MODE) 90kV pulsed mode (3D CB3D MODE) (automatically or manually selectable in steps of 1 kV)
Maximum continuous input anode power	42 W
Anode current	4 – 10 mA automatically or manually selectable in steps of 1 mA
Maximum deviations from declared values	kV: < 5% (10% IEC) mA: < 10% (20% IEC) Linearity error < 0.2 Variation coefficient < 0.05
X-ray tube	CEI OPX/105 TOSHIBA D-054SB
Focal spot size	0.5 mm according to IEC 336
Anode structural material	Tungsten (W)
Anode inclination	5° (CEI)
Anode thermal capacity	30 kJ (CEI)
Maximum continuous thermal dissipation	270 W (CEI)
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  <p>CURVE TERMICHE DI RISCALDAMENTO E DI RAFFREDDAMENTO</p> </div> <div style="text-align: center;">  <p>hyperion Tube housing heating and cooling curve - 85 kV 10 mA (duty cycle 1:20)</p> </div> </div>	
Inherent filtration	3.2 mm Al @ 85 kV (2D mode) - (+/-0,2mm) 6.2 mm Al @ 85kV (CB3D mode) - (+/-0,2mm)
Hemivalent layer	2D (PAN,CEPH): HLV=3.6mm @85kV 3D (CBCT mode): HLV=5.0mm @90kV (with automatic added filter)
Leaked radiation	<0.88mGy / h at 1 metre from focusing at 90 kV 10 mA, duty cycle 1:20
Maximum current	7A for the 230Vac version 14A for the 115Vac version at 85KV 10mA
Exposure time (PAN/CEPH/CB3D)	2D: 0.5s - 14s continuous radiation (depending on the 2D examination type selected) 3D: 18s-36s pulsed mode with actual ray total emission time of 3.6s- 9s (3D mode)
X-ray unit load factor	Max 20 s ON, 400 sec. OFF
Distance between source and primary collimator	100 mm
Distance between source and skin (PAN)	> 150 mm
Distance between source and skin (CEPH)	> 1000 mm
Distance between source and skin (CB3D)	> 150 mm
Output radiation intensity (output dose)	Typical 0.30 mGy/s @ 70 kV, 10 mA, 1 m SDD. (+/-30%) (1.0 mGy/s 70kV 550 mm 10mA)
Source - detector distance (SID) - PAN	550 mm
Source - detector distance (SID) - CEPH	1512 mm
Source - detector distance (SID) - CB3D	660 mm
Sensor life	at least 30,000 exposures


Generator reference axis:


10.3.RADIOLOGICAL CHARACTERISTICS IN CB3D MODE

Operating voltage	90kV with all projections
Generator voltage at maximum current	90kV - 10mA
Load factors for maximum output power	90kV - 10mA
Maximum rated output power (value averaged with N>4s)	810VA (90kV 10mA - 20/25% duty cycle)
Exposure time (CB3D)	18s-36s pulsed mode with actual ray total emission time of 3.6s- 9s
Pulsed mode operating cycle for CB3D exposure	20% - 90%

10.4.ISODOSE CURVES FOR CB3D EXAMINATIONS

10.5.ISODOSE CURVES FOR 2D EXAMINATIONS

10.6.CTDI (COMPUTED TOMOGRAPHY DOSE INDEX) MEASUREMENTS

The CTDI was measured using a cylindrical PMMA phantom with a diameter of 16 cm with holes at 12, 3, 6 and 9 o'clock and in the centre.

Diagram of PMMA cylinder for CTDI measurements

During the measurements, the holes that remained empty were filled with PMMA cylinders and the ionisation chamber was inserted in the hollow reference cylinder.

The possible operating conditions of the X-ray device consist of the following 11 scanning modes, distinguished by size of the field of view selected and by choice between high resolution and maximum resolution (HiRes):

[11x13e], [11x8], [11x8] HiRes, [11x5], [11x5] HiRes, [8x8], [8x8] HiRes, [8x5], [8x5] HiRes, [5x5], [5x5] HiRes.

he version with the smaller panel, 1308 DXT instead of 1313 DXT, has only 7 of these modes:

[11x8e], [11x5], [11x5] HiRes, [8x5], [8x5] HiRes, [5x5], [5x5] HiRes.

The load factors (mA, exposure time) are automatically determined by the system. The X-ray load automatically set by the machine for this dosimetric phantom for CTDI is reasonably similar to that of an average patient. The table indicates the X-ray load used, typically suggested by the X-ray device.

Type of examination	Technical factors	Central	Lateral (left of patient)	Lateral (right of patient)	Lateral (back of patient)	Lateral (front of patient)
[11x13e] ***	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	7.2	7.2	7.2	7.2	7.2
	mAs	28.80	28.80	28.80	28.80	28.80
	CTDI100 [mGy]	2.946	2.568	2.629	2.492	2.551
[11x8e] *	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	7.2	7.2	7.2	7.2	7.2
	mAs	28.80	28.80	28.80	28.80	28.80
	CTDI100 [mGy]	3.202	2.808	2.891	2.672	2.761
[11x8] **	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	3.6	3.6	3.6	3.6	3.6
	mAs	14.40	14.40	14.40	14.40	14.40
	CTDI100 [mGy]	2.104	1.838	1.878	1.780	1.822
[11x8] HiRes **	kV	90	90	90	90	90
	mA	6.00	6.00	6.00	6.00	6.00
	exposure time (s)	9.0	9.0	9.0	9.0	9.0
	mAs	54.00	54.00	54.00	54.00	54.00
	CTDI100 [mGy]	7.271	6.182	6.206	6.144	6.182
[11x5]	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	3.6	3.6	3.6	3.6	3.6
	mAs	14.40	14.40	14.40	14.40	14.40
	CTDI100 [mGy]	2.287	2.006	2.065	1.909	1.972

[11x5] HiRes	kV	90	90	90	90	90
	mA	6.00	6.00	6.00	6.00	6.00
	exposure time (s)	9.0	9.0	9.0	9.0	9.0
	mAs	54.00	54.00	54.00	54.00	54.00
	CTDI100 [mGy]	7.748	6.647	6.944	6.468	6.607
[8x8] **	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	3.6	3.6	3.6	3.6	3.6
	mAs	14.40	14.40	14.40	14.40	14.40
	CTDI100 [mGy]	2.080	1.519	1.664	1.448	1.600
[8x8] HiRes **	kV	90	90	90	90	90
	mA	6.00	6.00	6.00	6.00	6.00
	exposure time (s)	9.0	9.0	9.0	9.0	9.0
	mAs	54.00	54.00	54.00	54.00	54.00
	CTDI100 [mGy]	7.071	5.137	5.622	5.224	5.373
[8x5]	kV	90	90	90	90	90
	mA	4.00	4.00	4.00	4.00	4.00
	exposure time (s)	3.6	3.6	3.6	3.6	3.6
	mAs	14.40	14.40	14.40	14.40	14.40
	CTDI100 [mGy]	1.780	1.266	1.266	1.266	1.266
[8x5] HiRes	kV	90	90	90	90	90
	mA	6.00	6.00	6.00	6.00	6.00
	exposure time (s)	9.0	9.0	9.0	9.0	9.0
	mAs	54.00	54.00	54.00	54.00	54.00
	CTDI100 [mGy]	5.902	4.204	4.204	4.204	4.204
[5x5]	kV	90	90	90	90	90
	mA	5.00	5.00	5.00	5.00	5.00
	exposure time (s)	3.6	3.6	3.6	3.6	3.6
	mAs	18.00	18.00	18.00	18.00	18.00
	CTDI100 [mGy]	2.421	1.520	1.550	1.437	1.502
[5x5] HiRes	kV	90	90	90	90	90
	mA	7.00	7.00	7.00	7.00	7.00
	exposure time (s)	9.0	9.0	9.0	9.0	9.0
	mAs	63.00	63.00	63.00	63.00	63.00
	CTDI100 [mGy]	7.609	4.760	4.812	4.476	4.724

(*=available only with CB3D 1308 DXT detector, with Extra FOV software licence)

(**=available only with CB3D 1313 DXT detector)

(***=available only with CB3D 1313 DXT detector, with Extra FOV software licence)

Accuracy of declared values	+/-30% according to the standards
-----------------------------	-----------------------------------

10.7.CB3D DETECTOR CHARACTERISTICS

Pixel size	127 x 127 μm
Sensitive area dimensions (depending on the detector installed)	81x130, 130x130, 80x146, 146x146 mm
Resolution	3.94 LP/mm

10.8.PANORAMIC SENSOR CHARACTERISTICS(PAN)

Pixel size	48 x 48 μm
Sensitive area dimensions	6 x 146 mm
Resolution	5.2 LP/mm
Primary screen	> 0.5 mm Pb

10.9.TELERADIOGRAPHIC (CEPH) SENSOR CHARACTERISTICS (CEPH)


Pixel size	48 x 48 μm
Sensitive area dimensions	6 x 220 mm
Resolution	5.2 LP/mm
Primary screen	> 0.5 mm Pb

10.10. LASER CHARACTERISTICS

Optical power	Class 1 according to IEC 60825-1: 2003
Wavelength	635 nm
Output power	1.5 mW
Activation time	Continuous wave; Timing limited to 30"

10.11. DIMENSIONAL CHARACTERISTICS

Weight (basic machine)	180 Kg
Weight (cephalometric unit)	25 Kg
Maximum plan dimensions (basic machine)	1519 x 1312 mm
Maximum plan dimensions (with cephalometric unit installed)	1519 x 1829 mm
Height	Min 1590 mm Max 2380 mm

10.12. ENVIRONMENT CHARACTERISTICS

Operating conditions	Temperature + 10 - +35 °C
	Relative humidity 10 – 90%
	Pressure 710 – 1060 hPa
	Altitude a.s.l. < = 3000 m
	Pollution degree: 2
	CTI class: IIIb
Transport and storage conditions	Temperature -10 - +70 °C
	Relative humidity 10 – 90%
	Pressure 710 – 1060 hPa

10.13. PC REQUIREMENTS

Only for 3D devices.

CB3D acquisitions

Requirements for the PC workstation dedicated to CB3D primary reconstruction directly connected to the X-ray acquisition device.

COMPONENT	CONFIGURATION 1	CONFIGURATION 2	CONFIGURATION 3	CONFIGURATION 4
PROCESSOR	1 CPU INTEL I7 3770 - 3,4 GHZ – 4 CORE	1 CPU INTEL XEON E5-2630 - 2,3GHZ - 6 CORE	1 CPU INTEL XEON E3-1270 V3 – 3.50 GHZ – 4 CORE	CPU XEON E5-1620 - 3.60 GHZ - 4 CORE
HARD DISK SPACE ⁽¹⁾	1 X 500 GB SATA HARD DISK - 10K RPM	2 300 GB SATA HARD DISKS - 10K RPM – CONFIGURED IN RAID 0	1 X 500 GB SATA HARD DISK - 10KRPM	
RAM	8 GB DDR3 1600 MHZ, OR HIGHER	4 GB DDR3	8 GB DDR3 1600 MHZ, OR HIGHER	
POWER SUPPLY ⁽²⁾	>= 550 W	>= 800 W	>= 400 W	>= 600 W
VIDEO CARD ⁽³⁾	ATI RADEON HD 7770 OC GHZ EDIT. VAPORX 1GB GDDR5		ATI RADEON R7 250 1GB GDDR5	ATI RADEON HD 7770 OC GHZ EDIT. VAPORX 1GB GDDR5
OPERATING SYSTEM	WINDOWS 7 PROFESSIONAL 32 64 BIT SP1			
NETWORK CARD	2 ETHERNET INTERFACES (OF WHICH AT LEAST ONE INTEL PRO-1000)			
OTHER DEVICES	MINIMUM DISPLAY RESOLUTION 1280X1024, 64 MILLION COLOURS (OR HIGHER) DVD BURNER			

Notes:

- (1) Using a 7200rpm hard disk involves a considerable loss of performance.
- (2) To use the video cards tested according to specifications, you need to have at least power supply that meets the power consumption requirements of the chosen video card.
- (3) In addition to the video card indicated, see the complete list of tested video cards.

Validation tests were conducted with the following video cards:

VIDEO CARD	POWER SUPPLY
ATI RADEON HD 4850 / 5770 – 1GB	>= 450 W
ATI RADEON HD 4870 / 4890 / 5850 / 5870 – 1GB	>= 500 W
SAPPHIRE RADEON HD 6750 / 6770 – VAPORX – 1GB – RAM GDDR5	>= 450 W
SAPPHIRE RADEON HD 6850 / 6870 / 6950 / 7870 – VAPORX – 1GB – RAM GDDR5	>= 500 W
SAPPHIRE RADEON HD 7770 – GHZ ED. - OC - VAPORX – 1GB – RAM GDDR5	>= 500 W
SAPPHIRE RADEON HD 6970 – VAPORX – 2GB – RAM GDDR5	>= 550 W
SAPPHIRE RADEON HD 7850 - 1GB / 2GB - RAM GDDR5	>= 500 W
SAPPHIRE RADEON HD 7870 - 2GB - RAM GDDR5	>= 500 W
SAPPHIRE RADEON HD 7750 - OC - 1GB – RAM GDDR5	>= 400 W
AMD FIREPRO W7000 – 4GB GDDR5	>= 400 W
AMD/SAPPHIRE RADEON R7 250 1GB – RAM GDDR5	>= 400 W
AMD/SAPPHIRE RADEON R7 260X - RAM 1GB DDR5	>= 500 W
AMD/SAPPHIRE RADEON R9 270 - RAM 2GB DDR5	>= 500 W

2D acquisitions


Requirements for the PC workstation directly connected to the X-ray acquisition device if it is dedicated to two-dimensional examinations only. The minimum recommended requirements are indicated below. Anything less than indicated may result in unsatisfactory performance or impossibility to acquire X-ray examinations from this workstation.

Processor	Intel Core2-Duo / AMD Athlon X2 (or higher)
Operating system	(No Home Ed.): Windows XP 32bit only \ Windows Vista Pro. 32bit only \ Windows 7 Pro. 32 64 bit \ Windows 8 Professional 32 64 bit
Hard disk space	100 GB - 7200 rpm (min)
RAM memory	4 GB Vista-Windows7 2GB WindowsXP
Video card	3D VideoCard / 1 GB RAM
Network card	1 Ethernet 100Mbit interface, HUB-switch recommended
Other devices	Display resolution 1280x1024, 64 million colours (or higher); CD-ROM (recommended)

Additional workstations

For the minimum requirements of the workstations not directly connected to the X-ray acquisition device and where NNT Station or NNT Viewer will be installed, consult the NNT manual.

10.14. IDENTIFICATION LABEL POSITION

- 1 PAN sensor label
- 2 Warnings and laser hazard label
- 3 CB3D sensor label
- 4 Generator label
- 5 CEPH arm label
- 6 CEPH sensor identification label
- 7 WARNING and DHHS label
- 8 Main label
- 9 Certifications label

 Plate images are purely illustrative; refer to the plate placed on the device.

11. ERROR MESSAGES

Code	Message	Solution
USER COMMUNICATIONS		
	The shielded door connected to the X-ray device is open	Close the door or any other device to make X-ray emission safe for the operator.
WARNING		
0.1	X-ray Button released during exposure	Hold the X-ray emission button until the end of the procedure.
0.2	X-ray Button not released at the end of examination	<ol style="list-style-type: none"> 1) X-ray button was held for too long after the end of exposure (>15s). 2) A reset command was sent while the xray button was being pressed. 3) The device exited the stand-by status while the xray button was being pressed.
0.3	Emergency Button pressed	Release the emergency stop button.
0.4	POWER Circuit Board unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.5	KEYBOARD Circuit Board unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.6	CONVERTER Circuit Board unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.8	Circuit Board of CEPH-Detector Lift (CSENS) unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.9	Circuit Board of Headrest Lift (MSENS) unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.10	Circuit Board of Detectors Carousel (GSENS) unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.11	Circuit Board of COLUMN-INVERTER unable to communicate with MAIN Board	Turn off the unit, check CAN bus connections, wait 30sec and turn it back on. If the problem persists contact tech support.
0.13	2D Detector failed to turn ON	Check the removable sensor's contacts. Unplug the sensor and plug it back.
0.14	X-Ray Start-Command aborted	<ol style="list-style-type: none"> 1) Check Converter board diagnostic LEDs (>15s). 2) Verify CAN connector (KCB1) on converter board.
0.16	Relocate 2D Detector to PAN position	Relocate 2D Detector to PAN.

Code	Message	Solution
0.17	Relocate 2D Detector to CEPH position	Relocate 2D Detector to CEPH position
0.18	Check correct position of Craniostat	Check correct position of Craniostat
0.19	2D Detector is not recognized (run "Sensor registration")	Please serialize the sensor (from unit calibration software on PC)
0.20	2D Detector not calibrated	Calibrate the sensor (from unit calibration software on PC)
0.21	Serial Number Mismatch: not the same 2D Detector used for calibration	Please serialize the sensor (from unit calibration software on PC)
0.22	Detected Open Door (interlock is active)	Check interlock switch (door not closed?)
0.23	Image not downloaded from equipment, transfer to PC or USB before proceeding	New exposure requested while previous has not been downloaded yet
0.24	Please wait for confirmation before press	The xray button was pressed before starting the cbct exam procedure. For CBCT exams, launch "Patient scan" from within NNT and wait for an explicit request before pressing the xray button.
0.25	CEPH enabled	Only at first startup after CEPH arm installation, call tech support if persistent.
0.26	CMOS Time and Date Not Set. CMOS battery low?	Set time and date from PC/Check if clock battery on mainboard is dead
0.27	Device not configured	Set configuration from PC
0.28	Device not calibrated	Perform unit calibration
0.29	Please, let the patient release the bite, then press reset to move the chinrest to the new position	Chinrest is going to perform a long movement, so the patient must release the bite and let the chinrest free to move
0.30	Buttom pressed on startup	Detected buttons pressed during startup. Verify that all the keyboard buttons are working properly.
0.31	USB Key not found	
0.32	USB error on writing	USB key not compatible / Corrupted / Full
0.33	Unable to move detector's carousel automatically. Please, move it manually	

12. USER'S LICENCE CONTRACT

IMPORTANT: READ CAREFULLY

12.1. GENERAL CONDITIONS OF THE IMAGE SOFTWARE LICENCE

This licence exclusively applies to the software, understood as specific drivers and libraries for connection to and control of the digital X-ray system, and to the image display and storage software altogether identified as "NNT" and "NNT viewer" (hereinafter referred to as "Software") designed by Cefla S.C. - Imola (Italy), (hereinafter referred to as "Author") and handed to the customer (hereinafter referred to as "User"). These conditions are understood fully recognised and accepted at the time of program installation.

Proper installation and subsequent use of the program therefore implies unconditional acceptance of all the conditions.

12.1.1. USER'S LICENCE

Use of the software, which is and remains the property of the author to every effect, is granted to the user under a non-exclusive, free-of-charge user's licence for an indeterminate period of time.

The user's licence does not include the right to obtain the software in its source format or to make use of the related design or logic system documentation. The user therefore acknowledges that the author is and remains the sole, exclusive and legitimate holder of all rights of property, copyright and economic use of the software as well as related documentation. The user is aware and expressly acknowledges that it is impossible, in general, to develop software products to be entirely error-free as regards all possible applications and uses and that upgraded versions of the Software may introduce new errors which were not present in the previous version. The User therefore expressly accepts that the Software may contain errors. To this end, the User agrees, in his own interest, to make back-up copies of the files (data and images) regularly and diligently in order to be prepared for any possible malfunctioning of the Software product or of the systems the software is installed in.

Acceptance that the Software may contain errors is a condition that must be met in order for the Author to grant the User the Licence to use the Software Product.

This licence grants the Software end-user the sole right to use it in the ways described below and according to the additional conditions and warnings given in the "on-line" documentation.

12.1.2. COPYRIGHT

This Software is protected by all applicable laws, decrees, regulations and any other national provision and, in particular, by laws regarding rights of authorship and intellectual property rights as well as European Union directives and applicable international conventions.

The Software product, every single copy of it and every part or section of it (including – as mere examples and not comprehensively – every image, icon, photograph, video, text), any text or image and the electronic and "on-line" documentation regarding the Software product are the exclusive property of the Author or its suppliers and cannot be reproduced by the User in any form or by any means whatsoever. All rights regarding use of trademarks, product names and any other distinctive sign are reserved by the Author or respective legitimate owners.

In particular, the user agrees not to copy or reproduce, entirely or in part, the software he has been granted the right to use in printed form or in any form which can be read by a processor. This also applies to any other documentation regarding the software unless for safety reasons or to ensure proper functioning and then only with the express authorisation of the author.

12.1.3. USE OF THE SOFTWARE PRODUCT AND EXPRESS TERMINATION CLAUSE

- a. Software.** The User can install the Software and can use it under the conditions and within the limits established in this licence.

The User is not authorised to distribute the Software, whether free-of-charge or upon payment, and cannot provide services or develop products or software applications based on the Image Software or which make use of it in any way whatsoever.

It is expressly forbidden to make any changes to the software or incorporate it partially or fully in other programs. In the same way, no trademarks/brand names/numbers or any other item included in the software or its storage mediums may be removed, altered or copied.

- b. Restrictions regarding decoding, decompilation and disassembly.** The User cannot decode, decompile, disassemble, modify or translate the Software product except when expressly laid down by the law.

- c. **Prohibition of transfer / sub-licence / leasing.** The User cannot transfer, grant under a sub-licence, whether free-of-charge or upon payment, or let out under a leasing or any other arrangement, this Software to third parties.
Unless otherwise agreed, the user cannot allow third parties the use of the software even for promotional, demonstration or educational purposes.

If the User does not adhere to the conditions of use of the Software and, therefore, in the event of failure to observe, or violation of, the terms laid out in this paragraph at points a, b and c the licence shall be withdrawn. If this should happen, the Author has the right to request that the User destroys all copies of the Software and all component parts in its possession.

Under these circumstances the author shall have the right to claim for damages.

12.1.4. GUARANTEE AND EXCLUSIONS FROM THE GUARANTEE

The Software is supplied in its original state and the Author does not offer and shall not acknowledge any guarantee for original defects or defects which might subsequently arise and shall not make any guarantee of quality or trouble-free functioning as regards the Software. Furthermore the Author shall not offer or acknowledge any guarantee regarding compliance of the Software with that described in the electronic documentation made available "on-line" or in any other manner, except for the guarantee for the actual storage medium (USB PEN DRIVE) if this should be damaged or unusable.

Any other guarantee is excluded even if the Software is installed in or is part of other software applications developed by third parties. In addition, in relation to these applications, the Author expressly states that no checks of any kind have or will be made as regards the functioning of this product.

12.1.5. LIMITS OF RESPONSIBILITY

Under no circumstances shall the Author or its suppliers be held responsible for direct or indirect damage (including damage resulting from lost earnings, lack of earnings or savings, interruption of business, loss of information or data and other economic loss) caused to the User or third parties by the use or improper use of the Software product even if the author has been informed of the possibility of such damage. This limit of responsibility is not only applicable when the Software is used improperly and not as instructed in the Author's indications but also when the Software is used properly and in compliance with the Author's directions.

12.1.6. APPLICABLE LAWS, JURISDICTION AND COURT OF JURISDICTION

The conditions of this licence are subject to Italian law. Notwithstanding other standards bringing together individual laws or International Conventions, any disputes which might arise between the parts regarding the interpretation or execution of this agreement shall come under the exclusive jurisdiction of Italian law and the exclusive court of jurisdiction shall be the Law Court of Bologna.

12.1.7. WRITTEN FORM

Any modifications should be made in written form. Failure to observe the above shall result in the agreement being made null and void. The parties agree not to recognise the validity in law of any verbal agreements made by any person previously, at the same time or subsequently to the written contract and state that at no time from this moment on shall any attempt be made to validate any verbal agreement which might alter that laid down in the present general conditions. The author reserves the right to modify the present general conditions at any time and unilaterally and to communicate any such changes by the most suitable means.

12.1.8. TRANSLATION

As regards the translation of this licence into other languages it remains understood that, should any inconsistencies arise regarding interpretation, the Italian version shall prevail.

12.1.9. SOFTWARE COMPLIANCE WITH THE REGULATIONS IN FORCE

If the Author declares compliance of the Software with specific laws or regulatory standards, such compliance shall be deemed valid on the date the Author places the product on the market.

It being evident that legislative or regulatory modifications are always possible as are changes in the interpretation of regulations, and also considering that any modifications and changes would make the Software in possession of the User no longer compliant to such regulations, the User is obliged to from time to time ascertain whether the use he intends to make of the Software is legitimate (or in any case still legitimate) and abstain from using the Software in case of doubt about compliance with the law for the use he intends to make of the Software, promptly notifying the Author of the circumstances.

13. INSPECTION AND MAINTENANCE
13.1.USER INSPECTION

These instructions describe the maintenance procedures for the extraoral X-ray system.
These instructions apply to all the revisions of said equipment, as well as all the accessories that may have been provided, therefore the description of some parts may not correspond to your equipment.
Inspection and preventive maintenance must be performed at scheduled intervals to protect the health and safety of patients, users and other persons in accordance with national regulations regarding the use, maintenance of dental x-ray devices that are in force in the country where the device is installed.
In order to ensure the operational safety and functional reliability of your product, the system owner should check the equipment at regular intervals (at least once a year) or commission an authorized technician to do so.
If one or more checks to be performed are not satisfactory, please contact your dealer tech support.

			Answer questions with yes (✓) or no (–)				
			Inspection DATE				
Step	Description	Reference in User Manual	__/__/20__	__/__/20__	__/__/20__	__/__/20__	__/__/20__
1	Check that all labels located <ul style="list-style-type: none">- on the detector/detectors,- on the 3D Panel (if present)- at the base of the column,- on the X-ray tube housing, are intact, well attached and legible.	Section Identification label position					
2	Check there are no external damages to the equipment which could reduce protection against radiation.	Section Description of operations					
3	Check the remote emission button cable is not broken or scratched.	Section X-ray emission remote control					
4	Check the power switch verifying that the switch is working properly and the touch panel display illuminates when the switch is in the on position.	Section Starting the system					
5	Verify that exposure is immediately interrupted when X-rays button is released	Section X-ray emission remote control					
6	Check the touch console functionality: functions must respond to interrogation.	Section Touch screen display					
7	Check proper functioning of X-ray exposure LED and exposure buzzer	Section X-ray emission remote control					
8	Check proper functioning of X-ray penetrating laser traces	Section Laser traces					
9	Check if detector fits easily in the guides Pan and CEPH sides – only for unit with Cephalometric arm and one removable detector.	Section Sensor positioning					
10	Perform a motors movement test running the Dummy panoramic procedure, making sure the movement on the 3 axis X, Y,R (rotation and translation simultaneously) is smoothly and without noise. See figure 1, next page.	To execute a dummy run, exam simulation without emitting radiation, select any panoramic exam and reset the unit. Once the reset position is reached, hold down button “Reset” throughout the entire exam procedure, otherwise movement will stop.					
11	In case a Cephalometric arm is installed, perform a motors movement test of the cephalometric detector using the Dummy Ceph exam procedure, making sure the movement on the axe H (translation) is smoothly and without noise. See figure 1, next page.	To execute a dummy run, exam simulation without emitting radiation, select any cephalometric exam and reset the unit. Once the reset position is reached, hold down button “Reset” throughout the entire exam procedure, otherwise movement will stop.					
12	Perform column movements, up and down on the Z axis and check proper functioning. See figure 2, next page.	Section Control Panel (console onboard the machine)					
13	Check the emergency button functionality. Emergency button is used to stop X-ray device operation, it's located under the patient support arm in proximity of the telescopic column Emergency button).	User Manual, section Emergency button					
14	Check the X-Ray Generator functionality performing a complete trial exam exposure. Select any panoramic exam and reset the unit. Hold down the emission button throughout the entire exam procedure. Absence of error messages assures a proper generator functionality.	Section Performing a 2-D X-Ray examination					
If Quality Phantoms are not available at installation site, then contact your dealer tech support to have the Quality Assurance procedure performed at your premises							
15	Perform an exam on the 3D Quality Phantom and evaluate the quality of the outcome volumetric study.	Acquisition Operations with Annex					

Operator Name					
Signature					

The undersigned confirms that the equipment was checked for the above criteria and that, in case of any malfunction, an authorized technician of the local dealer was informed.

 All inspection and maintenance work performed by the system owner and/or service engineer must be recorded in this document and kept near the unit!

13.2.TECHNICAL MAINTENANCE

These instructions describe the maintenance procedures for the extraoral X-ray system. These instructions apply to all the revisions of said equipment.
In order to ensure the operational safety and functional reliability of the equipment installed, at least once a year an authorized service technician must perform a full inspection of the device.
When taking measurements that require a multimeter, always use a calibrated digital multimeter.
All the following tests will be carried out. Customer should be notified prior to replacing parts.

			Answer questions with yes (√) or no (–)		
			Inspection DATE		
Step	Description	Reference in User Manual	__/__/20__	__/__/20__	__/__/20__
1	"Check that all labels located - on the detector/detectors, - on the 3D Panel (if present) - at the base of the column, - on the X-ray tube housing, are intact, well attached and legible."	User Manual, section Identification label position			
2	Check there are no external damages to the equipment which could reduce protection against radiation.	User Manual, section Description of operations			
3	Open the generator covers and check there are no oil leakage on the tube-head; replace it if there are signs of leakages.	Setting Tube head assembly replacement			
4	Remove the dust accumulated inside the generator casing using a vacuum cleaner and reattach the covers.	Setting Tube head assembly replacement			
5	Check the remote emission button cable is not broken or scratched	User Manual, section X-ray emission remote control			
6	Power off the unit, disconnect it from the mains power supply and check the condition of the main power supply cable. Replace it in case of damage. Connect it back making sure the safety ground is securely connected.	Technical manual, section Electrical connections on the PFC board			
7	Check the power switch verifying that the ON/OFF button is working properly and the touch panel display illuminates when the switch is in the on position.	User Manual, section Starting the system			
8	Check the touch console functionality: functions must respond to interrogation.	User Manual, Touch screen display			
9	Check proper functioning of X-ray exposure LED and exposure buzzer	User Manual, X-ray emission remote control			
10	Check proper functioning of X-ray centring laser traces	User Manual, section Laser traces			
11	Check if detector fits easily in the guides Pan and CEPH sides – only for unit with Cephalometric arm and one removable detector.	User Manual, section Sensor positioning			
12	Perform a motors movement test running the Dummy panoramic procedure, making sure the movement on the 3 axis X, Y,R (rotation and translation simultaneously) is smoothly and without noise. See figure 1, next page.	To execute a dummy run, exam simulation without emitting radiation, select any panoramic exam and reset the unit. Once the reset position is reached, hold down button “Reset” throughout the entire exam procedure, otherwise movement will stop.			
13	In case a Cephalometric arm is installed, perform a motors movement test of the cephalometric detector using the Dummy Ceph exam procedure, making sure the movement on the axe H (translation) is smoothly and without noise. See figure 1, next page.	To execute a dummy run, exam simulation without emitting radiation, select any cephalometric exam and reset the unit. Once the reset position is reached, hold down button “Reset” throughout the entire exam procedure, otherwise movement will stop.			
14	Perform column movements, up and down on the Z axis and check proper functioning. See figure 2, next page.	User Manual, section Control Panel (console onboard the machine)			
15	Check the emergency button functionality. Emergency button is used to stop X-ray device operation, it's located under the patient support arm in proximity of the telescopic column	User Manual, section Emergency button			
16	Verify that exposure is immediately interrupted when X-rays button is released	User Manual, section X-ray emission remote control			
17	Check the X-Ray Generator functionality performing a complete trial exam exposure. Select any panoramic exam and reset the unit. Hold down the emission button throughout the entire exam procedure. Absence of error messages assures a proper generator functionality.	User Manual, section Performing a 2-D X-Ray examination			
18	Perform a complete 2D calibration of the unit : column calibration, PAN x-ray alignment, PAN detector calibration, PAN mechanical alignment, Laser test. In case of Cephalometric option perform also the followings calibration steps: Ceph x-ray alignment, Ceph detector calibration, Ceph mechanical alignment, Nasion calibration, Ear guide rings alignment.	Technical manual, section 2D Calibration			
19	At the end make a calibration backup	Technical manual, section Calibration backup			
20	Performa complete 3D calibration of the unit: Beam Limiter Test, Daily Check, Cylindrical test phantom acquisition (“Mattarellum”).	Technical manual, section 3D Calibration			
21	At the end make a 3D calibration backup	Technical manual, section 3D Calibration backup			
22	Perform an exam on the 3D Quality Phantom and evaluate the quality of the outcome volumetric study.	Acquisition Operations with Annex			

Operator Name			
Signature			

The undersigned confirms that he/she has checked the unit for the above criteria and that he/she has handed over the unit in fully functional condition.


 All inspection and maintenance work performed by the system owner and/or service engineer must be recorded in this document and kept near the unit!!

Figure 1



Figure 2

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