RESTORATION



#COLTENE

WHAT DO YOU CALL A CONVICTION THAT EVERYTHING CAN ALWAYS BE A LITTLE BETTER, EASIER AND MORE RELIABLE?

WE CALL IT Upgrade Dentistry



TABLE OF CONTENTS

1 ETCH & BOND

Etchant Gel S

One Coat	10
ONE COAT 7 UNIVERSAL	11
One Coat 7.0 Activator	11
One Coat Bond	12
One Coat Self-Etching Bond	12
ParaBond®	13
A.R.T. Bond	13
2 CONVENTIONAL COMPOSITES	
MIRIS*2	16
MIRIS ² Kits	16
MIRIS ² Refills	17
MIRIS ² Accessories	17
BRILLIANT EverGlow [®]	18
BRILLIANT EverGlow Kits	19
BRILLIANT EverGlow Refills	20
BRILLIANT EverGlow Accessories	20
BRILLIANT EverGlow Flow Kits	21
BRILLIANT EverGlow Flow Refills	21
BRILLIANT EverGlow Flow Accessories	22
SYNERGY® D6	22
SYNERGY D6 Intro Kits	23
SYNERGY D6 Refills	23
SYNERGY D6 Accessories	24
SYNERGY D6 Flow Kits	24
SYNERGY D6 Flow Refills	25
SYNERGY D6 Flow Accessories	25
SYNERGY® Nano Formula	25
Dispenser for composite tips	26
ROEKO matrix bands	26



3 BULK-FILL COMPOSITE

FIII-Upl Accessories 31 4 CAD/CAM SOLUTIONS BRILLIANT Crios 34 BRILLIANT Crios Block Refilis 35 BRILLIANT Crios Block Refilis 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM 36 BRILLIANT COMPONEER RIT Syringes 36 BRILLIANT COMPONEER RIT Tips 36 BRILLIANT COMPONEER Refilis 35 6 LUTING MATERIALS 35 DuoCem* 42 SoloCem* & DuoCem* Accessories 42 7 CORE BUILD-UP 46 ParaCore* 5 ml 46 ParaCore* 5 ml Refilis 47 ParaCore* 25 ml Refilis 48	Fill-Up!		30
### CAD/CAM SOLUTIONS #### BRILLIANT Crios ##### BRILLIANT Crios Block Refilis ##################################	Fill-Up! f	Refills	31
BRILLIANT Crios 34 BRILLIANT Crios Block Refills 35 BRILLIANT Crios Disc Refills 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 38 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER RETTips 35 BRILLIANT COMPONEER Refills 35 BRILLIANT COMPONEER Refills 35 BRILLIANT COMPONEER Refills 35 BOLOCem® 42 SoloCem® & DuoCem® Accessories 42 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	Fill-Up!	Accessories	31
BRILLIANT Crios 34 BRILLIANT Crios Block Refills 35 BRILLIANT Crios Disc Refills 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 38 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER RETTips 35 BRILLIANT COMPONEER Refills 35 BRILLIANT COMPONEER Refills 35 BRILLIANT COMPONEER Refills 35 BOLOCem® 42 SoloCem® & DuoCem® Accessories 42 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	4		
BRILLIANT Crios Block Refills 35 BRILLIANT Crios Disc Refills 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 38 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER Sets 35 BRILLIANT COMPONEER Refills 35 6 LUTING MATERIALS 42 SoloCem® 43 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP 46 ParaCore® 5 ml 46 ParaCore® 7 ml Refills 47 ParaCore® 25 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	4	CAD/CAM SOLUTIONS	
BRILLIANT Crios Disc Refilis 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 38 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER Fets 35 BRILLIANT COMPONEER Refilis 35 6 LUTING MATERIALS 42 SoloCem* 43 SoloCem* & DuoCem* Accessories 43 7 CORE BUILD-UP 46 ParaCore* 5 ml 46 ParaCore* 5 ml Refilis 47 ParaCore* 25 ml 48 ParaCore* 25 ml Refilis 48 ParaCore* 25 ml Refilis 48	BRILLIA	ANT Crios	34
BRILLIANT Crios Disc Refilis 35 BRILLIANT Crios Block 35 5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 36 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER Refils 35 BRILLIANT COMPONEER Refils 35 BRILLIANT COMPONEER Refils 35 6 LUTING MATERIALS 42 SoloCem® 42 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refils 47 ParaCore® 25 ml 48 ParaCore® 25 ml 48 ParaCore® 25 ml Refills 48			35
5 COMPOSITE VENEERING SYSTEM BRILLIANT COMPONEER KIT Syringes 38 BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER Sets 35 BRILLIANT COMPONEER Refilis 35 6 LUTING MATERIALS DuoCem* 42 SoloCem* 43 TOCORE BUILD-UP 46 ParaCore* 5 ml 46 ParaCore* 5 ml Refilis 47 ParaCore* 25 ml 48 ParaCore* 25 ml Refilis 48 ParaCore* 25 ml Refilis 48	BRILLIAN	NT Crios Disc Refills	35
BRILLIANT COMPONEER KIT Syringes BRILLIANT COMPONEER KIT Tips BRILLIANT COMPONEER Sets BRILLIANT COMPONEER Refills 33 6 LUTING MATERIALS DuoCem* 42 SoloCem* 42 SoloCem* & DuoCem* Accessories 43 7 CORE BUILD-UP ParaCore* 5 ml ParaCore* 5 ml ParaCore* 25 ml ParaCore* 25 ml ParaCore* 25 ml Refills 48 ParaCore* 25 ml Refills 48	BRILLIAN	NT Crios Block	35
BRILLIANT COMPONEER KIT Tips 38 BRILLIANT COMPONEER Sets 39 BRILLIANT COMPONEER Refilis 39 6 LUTING MATERIALS 42 SoloCem® 43 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refilis 47 ParaCore® 25 ml 48 ParaCore® 25 ml Refilis 48 ParaCore® 25 ml Refilis 48	5	COMPOSITE VENEERING SYSTEM	
BRILLIANT COMPONEER Sets 35 BRILLIANT COMPONEER Refills 35 6 LUTING MATERIALS 42 DuoCem® 43 SoloCem® 43 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP ParaCore® 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	BRILLIAN	NT COMPONEER KIT Syringes	38
BRILLIANT COMPONEER Sets 36 BRILLIANT COMPONEER Refilis 35 6 LUTING MATERIALS 42 DuoCem° 42 SoloCem° 43 SoloCem° & DuoCem° Accessories 43 7 CORE BUILD-UP 46 ParaCore° 5 ml 46 ParaCore° 5 ml Refills 47 ParaCore° 25 ml Refills 48 ParaCore° 25 ml Refills 48 ParaCore° 25 ml Refills 48			38
BRILLIANT COMPONEER Refills 6 LUTING MATERIALS DuoCem® 42 SoloCem® 43 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	BRILLIAN	NT COMPONEER Sets	39
DuoCem® 42 SoloCem® 43 SoloCem® & DuoCem® Accessories 43 7 CORE BUILD-UP ParaCore® 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48			39
SoloCem® & DuoCem® Accessories 7 CORE BUILD-UP ParaCore® ParaCore® 5 ml ParaCore® 5 ml Refills ParaCore® 25 ml Refills ParaCore® 25 ml Refills 48	6	LUTING MATERIALS	
SoloCem® & DuoCem® Accessories 7 CORE BUILD-UP ParaCore® 46 ParaCore® 5 ml ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48	DuoCen	m [®]	42
7 CORE BUILD-UP ParaCore® 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	SoloCen	n®	43
ParaCore® 46 ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml 48 ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills 48	SoloCen	n® & DuoCem® Accessories	43
ParaCore® 5 ml 46 ParaCore® 5 ml Refills 47 ParaCore® 25 ml 48 ParaCore® 25 ml Refills 48	7	CORE BUILD-UP	
ParaCore® 5 ml Refills ParaCore® 25 ml ParaCore® 25 ml Refills 48 ParaCore® 25 ml Refills	ParaCo	re®	46
ParaCore® 25 ml ParaCore® 25 ml Refills 48	ParaCore	e ^o 5 ml	46
ParaCore® 25 ml Refills 48	ParaCore	e® 5 ml Refills	47
	ParaCore	e° 25 ml	48
Accessory Core 49	ParaCore	e® 25 ml Refills	48
	Accesso	ory Core	49

TABLE OF CONTENTS

8 TEMPORARY MATERIALS

Coltosol® F	52
DuoTEMP®	52
Cool Temp® NATURAL	53
Cool Temp® NATURAL Refills	53
Cool Temp® NATURAL Accessories	53
TempoSIL® 2	54
9 CURING LIGHTS	
S.P.E.C.® 3	58
S.P.E.C.® 3 Accessories	58
Coltolux® LED	59
Coltolux® LED Accessories	59
10 PARAPULPAL PINS	
MAX® Pin	62
MAX® Pin Accessories	62
TMS	62
TMS LINK PLUS®	63
TMS LINK® Series	64
TMS®	65
KODEX®-Drills	66
Accessories Parapulpal Pins	66







SPOTLIGHT



ONE COAT 7 UNIVERSAL

11

ONE COAT 7 UNIVERSAL is a universal, light-cured, one-component bonding agent able to be used with self-etch, selective etch or total etch technique for adhesive restorations.

Etchant Gel S



35 % phosphoric acid etching gel with precise placement, usable for enamel and dentin.

ltem	
Kit	
$2 \times 2.5 \text{ ml}$	Etchant Gel S
$1 \times 20 \text{ pcs}$	Application needles
	Kit 2 × 2.5 ml

REF	Accessories
7876	Application Needles Black, 20 pcs

One Coat

Dentists from all over the world rely for years on the proven One Coat adhesives of COLTENE. The light-curing bonding systems provide adhesion to various surfaces and distinguish with outstanding adhesion values and excellent marginal integrity.









ONE COAT 7 UNIVERSAL



60019538



60019540

ONE COAT 7 UNIVERSAL is a universal, light-curing, one-component bonding agent for adhesive restoration technique, applicable in self-etch, selective-etch or total-etch technique.

Optionally, ONE COAT 7 UNIVERSAL can also be used with One Coat 7.0 Activator for chemically curing materials.

- Universal bonding to dentin, enamel, ceramics, zirconia, metals and composite
- Excellent and reliable bond strength no matter which etching technique used
- Great marginal sealing for long term stability
- User-friendly handling
- Antibacterial (Proliferation test for antimicrobial activity with S.epidermidis by QualityLabs, Germany)

REF	Accessories
7634	Microbrushes Black, 50 pcs
7876	Application Needles Black, 20 pcs

One Coat 7.0 Activator



7054

One Coat 7.0 Activator is an activator to provide the use with chemically and dual cured materials. One Coat 7.0 Activator is used exclusively in combination with One Coat 7.0 or ONE COAT 7 UNIVERSAL.

REF	ltem
7054	Activator Kit
	1 × 3 ml One Coat 7.0 Activator 1 × 50 pcs Microbrush, black

One Coat Bond



60011757



7879

Total-Etch light curing dentin/enamel bonding system for universal indications, prevents postoperative sensitivity.

- Gel-like consistency does not slump or drip for safe application
- Solvent free no evaporation, no drying out
- Super concentrated in average 100 drops per syringe

REF	Item
60011757	Intro Kit 1 × 1.2 ml One Coat Bond
	1 × 2.5 ml Etchant Gel S 1 × 50 pcs Microbrush, black 1 × 20 pcs Application needles
7873	Refill, 2 × 1.2 ml
7879	Refill, 1 × 1.2 ml

REF	Accessories
7634	Microbrushes Black, 50 pcs
7876	Application Needles Black, 20 pcs

One Coat Self-Etching Bond



7630

Formulated to provide an easy technique, One Coat Self-Etching Bond eliminates the need for an extra etching step. Primer 1 solubilizes the smear layer, without totally removing it as total etch bonds do. Bond 2 penetrates through the smear layer into the tubules to create the so called resin tags, which are responsible for its outstanding retention

- Consistent and homogenous bond layer which prevents postoperative sensitivity
- High Shear Bond Strength
- Excellent Marginal Integrity

REF	ltem	
7630	Intro Kit	
	$1 \times 5 \text{ ml}$	One Coat Self-Etching Bond Primer
	$1 \times 5 \text{ ml}$	One Coat Self-Etching Bond
	$1 \times 50 \text{ pcs}$	Microbrush, black
	$1 \times 50 \text{ pcs}$	Microbrush, yellow
	1 ×	Mixing well
7635	Refill Primer, 1 × 5 ml	
7636	Refill Bon	nd, 1 × 5 ml

REF	Accessories
7634	Microbrushes Black, 50 pcs
7633	Microbrushes Yellow, 50 pcs



ParaBond®



7486

Chemically cured, self-conditioning adhesive system for enamel and dentine which consists of a non-rinse conditioner and a chemically curing adhesive (adhesives A and B).

- Seals the tooth surface and minimizes the risk of postoperative sensitivities
- Ideal for situations where light might not penetrate, such as for post cementation, metal crowns, deep cavities or opaque ceramic
- High shear bond strengths on dentine and enamel
- Universally applicable with chemically-, dual- and light-curing materials

REF	ltem	
7486	Adhesive	e Kit
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A
	1 × 3 ml	ParaBond Adhesive B
7493	Adhesive	e Refill A&B
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B
7494	Non-Rin	se Conditioner Refill, 1 $ imes$ 3 ml

A.R.T. Bond



7362

Classical light curing dentine/enamel bonding system for universal indications, high bond strength, withstands the test of time, scientifically recommended.

REF	ltem		
7362	Kit		
	$1 \times 5 \text{ ml}$	A.R.T BOND Primer A	
	$1 \times 5 \text{ ml}$	A.R.T BOND Primer B	
	$1 \times 5 \text{ ml}$	A.R.T BOND	
	2 × 50 pcs	Brushes	
	1×	Brushholder black	
	1×	Mixing well	
	1 ×	Cover mixing well	
	1 ×	Brushholder white	
7363	Refill Bor	nd	
	1 × 5 ml	A.R.T BOND	
	40 ×	Mixing well	
7364	Refill Prir	ner A&B	
	$1 \times 5 \text{ ml}$	A.R.T BOND Primer A	
	$1 \times 5 \text{ ml}$	A.R.T BOND Primer B	
	40 ×	Mixing well	

SPOTLIGHT



BRILLIANT EverGlow[®]

SYNERGY® D6

18

22

BRILLIANT EverGlow is the latest composite line of COLTENE. Developed to provide an ideal mixture of high, long-lasting gloss and shade aesthetics, handling convenience and durability, it is a true all-round material. BRILLIANT EverGlow complies with the VITA based duo shade system and is availabe in 3 opaque, 7 universal and 2 translucent shades as well as in a flowable consistency.

SYNERGY D6 is a nano-hybrid composite for everyday dentistry. The simple selection of shades is a special feature. SYNERGY D6 complies with the duo shade system. The product line is available in two consistencies.

MIRIS^{®2}

Artificial Reality

MIRIS2 offers a unique shade system on the basis of the natural layering concept, which works with dentin and enamel masses whose optical properties are comparable with those of the original tissue. For this purpose, the shades and translucency of natural teeth in adolescent, adult and elderly patients were measured and emulated. For achieving an aesthetically optimal definition of shade, the MIRIS2 Shade Guide is separated into composite-based dentin cores and enamel shells that can be freely interposed on one another.

In addition to dentine and enamel shades, special effects shades complete the system to allow customised restorations, taking individual aspects into account, and thus achieve unsurpassed aesthetics - even in highly complex cases.

MIRIS² Kits

ltem	Content																	
			DENTIN								ENAMEL					EFFECT		
		S0	S 1	S2	S3	S4	S5	S6	S7	WB	WR	NR	NT	IR	W/O	Blue	Gold	White
8431	Syringe Set 13 Composite Syringes, 4 g each (52 g) 4 Effect Shade Syringes, 2.3 g each (9.2 g) 8 Effect Shade Application Needles 1 Shade Guide 2 Mixing Pads	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
8465	Tips Set 78 Composite Tips, 0.25 g each (19.5 g) 4 Effect Shade Syringes, 2.3 g each (9.2 g) 8 Effect Shade Application Needles 1 Shade Guide 2 Mixing Pads	6	6	6	6	6	6	6	6	6	6	6	6	6	1	1	1	1







8431 8465 8411



MIRIS² Refills

REF		Content																	
	DENTIN											ENAMEL				EFF	ECT		
			SO	S1	52	53	S4	S5	S6	S7	WB	WR	NR	NT	IR	W/O	Blue	Gold	White
Syringe	1 pcs	4 g	8409	8432	8433	8434	8435	8436	8437	8438	8440	8439	8441	8442	8443	=	-	=	-
Syringe	1 pcs	2.3 g	-	-	-	-	-	-	-	-	-	-	-	-	-	8446	8403	8447	8445
Tip	10 pcs	0.25 g	8406	8472	8473	8474	8475	8476	8477	8478	8480	8479	8481	8482	8483	-	-	-	-
Tip	20 pcs	0.25 g	8408	8486	8487	8488	8489	8490	8491	8492	8494	8493	8495	8496	8497	-	-	-	-



MIRIS² Accessories



REF	Item
8411	Shade Guide, 1 pcs
60022485	Composite Tip Dispenser, 1 pcs

BRILLIANT EverGlow®

Enduring gloss - made brilliant

With BRILLIANT EverGlow highly aesthetic and long-lasting restorations can be achieved effortlessly.

- Exceptional polishability and long-lasting brilliance
- Aesthetic single-shade restorations
- Versatile shade system with three translucency levels
- Excellent sculptability and smooth consistency
- Good wettability with minimal stickiness to the instrument





BRILLIANT EverGlow Kits

The universal submicron hybrid composite BRILLIANT EverGlow excels through convenient sculptability and form stability, excellent polishability and a long-lasting high lustre. Its sophisticated shade system offers three translucency levels and enhanced Duo Shades that allow for a flexible use for single-shade and multi-shade layered restorations.

ltem	Content											
		OPA	\QUE			ļ	JNIVERSA	L			TRANS	LUCENT
		OA1	OA3	BL	A1/B1	A2/B2	A3.5/B3	A3/D3	A4/C4	C2/C3	BLT	TS
60019720	Trial Kit Syringes 4 Syringes, 3 g each (12 g) 1 Shade Guide 1 Mixing Pad	-	-	-	-	1	1	1	-	-	-	1
60019721	Trial Kit Tips 32 Tips, 0.2 g each (6.4 g) 1 Shade Guide 1 Mixing Pad	-	-	-	-	8	8	8	-	-	-	8
60019722	System Kit Syringes 9 Syringes, 3 g each (27 g) 1 Shade Guide 2 Mixing Pads	-	-	1	1	1	1	1	1	1	1	1
60019723	System Kit Tips 60 Tips, 0.2 g each (12 g) 1 Shade Guide 2 Mixing Pads	-	-	5	8	8	8	8	5	5	5	8
60019726	Regular Kit Syringes 9 Syringes, 3 g each (27 g) 1 Shade Guide 2 Mixing Pads	1	1	-	1	1	1	2	-	1	-	1
60019727	Regular Kit Tips 60 Tips, 0.2 g each (12 g) 1 Shade Guide 2 Mixing Pads	5	5	_	5	10	10	12	-	5	8	-
60020090	Restoration Highlight Kit Syringes 3 BRILLIANT EverGlow Syringes, 3 g each (9 g) 1 ONE COAT 7 UNIVERSAL (5 ml) + 50 Brushes 2 DIATECH ShapeGuard Polishers	-	-	-	1	1	-	1	-	-	_	-
60020091	Restoration Highlight Kit Tips 30 BRILLIANT EverGlow Tips, 0.2 g each (6 g) 1 ONE COAT 7 UNIVERSAL (5 ml) + 50 Brushes 2 DIATECH ShapeGuard Polishers	-	-	-	10	10	-	10	_	-	_	-







60019722 6002090 60019727

BRILLIANT EverGlow Refills

REF			Conten	t										
				OPAQUE					UNIVERSAL				TRANS	LUCENT
		OA1	OA3	OBL	BL	A1/B1	A2/B2	A3/D3	A3.5/B3	C2/C3	A4/C4	BLT	TS	
Syringe	1 pcs	3 g	60019744	60019745	60019743	60019697	60019698	60019699	60019700	60019701	60019702	60019703	60019705	60019704
Tip	20 pcs	0.2 g	60019747	60019748	60019746	60019706	60019707	60019708	60019709	60019710	60019716	60019717	60019719	60019718





60019707



60019708

BRILLIANT EverGlow Accessories



REF	Item
60021919	Shade Guide, 1 pcs
60022485	Composite Tip Dispenser, 1 pcs



BRILLIANT EverGlow Flow Kits

BRILLIANT EverGlow Flow is based on the same filler technology as BRILLIANT EverGlow and is therefore the perfect complement. Thanks to its well-adjusted thixotropy, BRILLIANT EverGlow Flow has a low viscousity yet offers high stability. A new, extra thin application needle extends the range of accessories and allows for highest application precision.

Item	Content											
		UNIVERSAL										
		A1/B1	A2/B2	A3.5/B3	A3/D3							
60019760	Intro Kit 4 Syringes, 2 g each (8 g) 40 Flow Application Needles (20 x Standard, 20 x Fine)	1	1	1	1							



60019760

BRILLIANT EverGlow Flow Refills

Item			REF										
			OPAQUE	OPAQUE UNIVERSAL									
		OA2	A1/B1	A2/B2	A3.5/B3	A3/D3	A4/C4	BL	TS				
Syringe	1 pcs	2 g	60019751	60019753	60019754	60019756	60019755	60019757	60019752	60019758			



BRILLIANT EverGlow Flow Accessories



REF	Item
60022137	Flow Application Needles Fine, ID 0.4mm, 20 pcs
8469	Flow Application Needles Standard, ID 0.8 mm, 20 pcs

60022137

SYNERGY® D6

Beauty made simple

SYNERGY D6 is a nano-hybrid composite for the everyday use. Its consistency is firmer compared to BRILLIANT EverGlow.

- Simple and economic shade range
- Excellent blend-in effect
- Firm consistency with good sculptability
- Quick polishability
- Proven quality







SYNERGY D6 Intro Kits

Content	REF								
			ENAMEL						
		A1/B1	A2/B2	A3.5/B3	A3/D3	C2/C3	WB	wo	Univ.
8362	Intro Kit Syringes 8 Syringes, 4 g each (32 g) 1 Shade Guide 2 Mixing Pads	1	1	1	1	1	1	1	1
8363	Intro Kit Tips 80 Tips, 0.25g each (20 g) 1 Shade Guide 2 Mixing Pads	10	10	10	10	10	10	10	10





8363

SYNERGY D6 Refills

REF			Content									
			DENTIN ENAMEL									
			WB	BL/O	A1/B1	A2/B2	A3/D3	A3.5/B3	C2/C3	A4/C4	Univ.	wo
Syringe	1 pcs	4 g	8365	60013814	8366	8367	8368	8369	8370	8384	8371	8372
Tip	10 pcs	0.25 g	8375	60013813	8376	8377	8378	8379	8380	8383	8381	8382





SYNERGY D6 Accessories



REF	Item
8410	Shade Guide, 1 pcs
60022485	Composite Tip Dispenser, 1 pcs

SYNERGY D6 Flow Kits

The formula of Synergy D6 Flow is based on the technology of SYNERGY D6. The material is stable, yet flows gently under pressure to ensure easy handling and precision application.

REF	ltem						
				DENTIN			ENAMEL
		A1/B1	A2/B2	A3.5/B3	A3/D3	White Bleach	Universal
7650	Intro Kit 4 Syringes, 2.3 g each (9.2 g) 32 Flow Application Needles (Standard, Ø 0.8 mm)	-	1	-	1	1	1
7660	Dentin Intro Kit 4 Syringes, 2.3 g each (9.2 g) 32 Flow Application Needles (Standard, Ø 0.8 mm)	1	1	1	1	-	-







7650 7650 7650



SYNERGY D6 Flow Refills

REF		Content	Content									
				DEI	ITIN			ENAMEL				
		A1/B1	A2/B2	A3/D3	A3.5/B3	A4/C4	WB	Univ.				
Syringe	1 pcs 2.3 g	7651	7652	7653	7654	7655	7656	7657				



7651

SYNERGY D6 Flow Accessories



8469

REF Item

8469 Flow Application Needles Standard, ID 0.8 mm, 20 pcs

SYNERGY® Nano Formula

REF		Content									
		A1/D2	A2/B2	A3/D3	A3.5/B3	A4/M5	C2/C3				
Tip	20 pcs 0.25 g	7831	7832	7833	7834	7835	7836				



7831

Dispenser for composite tips



REF	Item
60022485	Composite Tip Dispenser, 1 pcs

ROEKO matrix bands




590305



REF	ltem
600009	Strip, Transparent Polyester Matrix Band, Length 15m, Width 9mm, 1 pcs
590305	Steel Matrix Band, Stainless Steel, Thickness 0.03mm, Length 3m, Width 5mm, 1 pcs
590306	Steel Matrix Band, Stainless Steel, Thickness 0.03mm, Length 3m, Width 6mm, 1 pcs
590307	Steel Matrix Band, Stainless Steel, Thickness 0.03mm, Length 3m, Width 7mm, 1 pcs
590505	Steel Matrix Band, Stainless Steel, Thickness 0.05mm, Length 3m, Width 5mm, 1 pcs
590506	Steel Matrix Band, Stainless Steel, Thickness 0.05mm, Length 3m, Width 6mm, 1 pcs
590507	Steel Matrix Band, Stainless Steel, Thickness 0.05mm, Length 3m, Width 7mm, 1 pcs
591505	Steel Matrix Band, soft, Stainless Steel, Thickness 0.05mm, Length 3m, Width 5mm, 1 pcs
591506	Steel Matrix Band, soft, Stainless Steel, Thickness 0.05mm, Length 3m, Width 6mm, 1 pcs
591507	Steel Matrix Band, soft, Stainless Steel, Thickness 0.05mm, Length 3m, Width 7mm, 1 pcs



NOTES







Fill-Up! 30

Fill-Up! is a dual curing, medium viscous bulk composite of the newest generation applied in one single step.

Fill-Up!



60019342



60019355

Welcome on the safe side!

Fill-Up! is a dual curing, medium viscous bulk composite for posterior Class I and II restorations. Owing to its dual curing property, fillings with Fill-Up! can be administered in arbitrary filling depth in a single step! An additional covering layer can optionally be applied to enhance aesthetics but is not necessary.

Conventional light curing bulk filling materials are limited in their curing depth, putting into question whether the restoration fully cures. On the contrary, the dual curing property of Fill-Up! guarantees a thorough curing down to the bottom of every cavity. Moreover, chemical curing minimises shrinkage stress dramatically, preventing microfractures and postoperative sensitivities. As supported by study, the bonding partners ParaBond and ONE COAT 7 UNIVER-SAL (together with the activator) cater for a perfect marginal seal comparable to the one of conventional composite - before and after thermocycling and chewing abrasion – ensuring a secure long-term restorative solution.

Fill-Up! is your perfect choice for all Class I and II fillings, cavity lining and core build-ups – every time when a fast restoration is needed or for cavities difficult to reach with the curing light.

- Unlimited depth polymerisation due to dual curing system
- Guaranteed single-layer technique, top layer optional but not necessary
- Minimal shrinkage stress with chemical curing
- Optimal marginal sealing with the matched ParaBond or ONE COAT 7 UNIVERSAL bonding systems
- Convenient automix application

REF	Item
60019342	Intro Kit with ParaBond 2 x 4.5 g Fill-Up! Universal 1 x 3 ml ParaBond Adhesive A 1 x 3 ml ParaBond Adhesive B 6 x Mixing Tip Short Super Fine 8 x Mixing Tip Short Fine 1 x 2.5 ml Etchant Gel S
60019973	Intro Kit with ONE COAT 7 UNIVERSAL 2 × 4.5 g Fill-Up! Universal 8 × Mixing Tip Short Fine 6 × Mixing Tip Short Super Fine 1 × 5 ml ONE COAT 7 UNIVERSAL 1 × 3 ml One Coat 7.0 Activator
60019355	Economy Kit 5 x 4.5 g Fill-Up! Universal 20 x Mixing Tip Short Fine 15 x Mixing Tip Short Super Fine



Fill-Up! Refills



REF	Item
60019341	Refill Syringe, 1 × 4.5 g

Fill-Up! Accessories



REF	Item
6747	Mixing Tip Short Fine, 40 pcs
6759	Mixing Tip Short Super Fine, 40 pcs





SPOTLIGHT



BRILLIANT Crios

34

Composite has long been the gold standard in direct filling therapy. With BRILLIANT Crios, the well-proven dental material entered the indirect CAD/CAM restoration technology. Compared to ceramics, CAD/CAM composites offer considerable advantages.

BRILLIANT Crios

BRILLIANT Crios is a reinforced composite. The material combines all the advantages of an innovative submicron hybrid composite material with those of a CAD/CAM fabrication process for reliable, aesthetic and fast restorations without a separate firing process. Consistently controlled fabrication with constant thermal curing, gives the reinforced composite excellent mechanical properties. Two sizes and two translucencies with 13 shades offer a broad range for a daily use in the dental practice.

- Tooth-like modulus of elasticity for a shock-absorbing effect and pleasant bite feeling
- High resistance and low abrasion on antagonist
- Exceptional grinding precision for more freedom of preparation
- Effortless polishing for fast gloss
- Blends in extremely well to give natural aesthetics





BRILLIANT Crios Block



REF	ltem	
60020097	Intro Kit	CEREC 14
	1 ×	BRILLIANT Crios A2 LT 14
	1 ×	BRILLIANT Crios A3 LT 14
	1 ×	BRILLIANT Crios A3.5 LT 14
	1 ×	BRILLIANT Crios A2 HT 14
	1 ×	BRILLIANT Crios A3 HT 14
	$5 \times 0.1 \text{ ml}$	ONE COAT 7 UNIVERSAL
60020095	Intro Kit	CEREC 12
	1 ×	BRILLIANT Crios A2 LT 12
	1 ×	BRILLIANT Crios A3 LT 12
	1 ×	BRILLIANT Crios A3.5 LT 12
	1 ×	BRILLIANT Crios A2 HT 12
	1 ×	BRILLIANT Crios A3 HT 12
	5 × 0.1 ml	ONE COAT 7 UNIVERSAL



60020095

BRILLIANT Crios Block Refills

Item		REF	REF											
	Low Translucent High Translucent													
Size	Content	BL	A1	A2	А3	A3.5	B1	B2	В3	C2	A1	A2	А3	B1
12	5 pcs	60019998	60019999	60020000	60020001	60020002	60020003	60020004	60020005	60020006	60020007	60020008	60020009	60020010
14	5 pcs	60019985	60019986	60019987	60019988	60019989	60019990	60019991	60019992	60019993	60019994	60019995	60019996	60019997

BRILLIANT Crios Disc Refills

Item			REF						
			Low Translucent				High Translucent		
Height	Ø mm	Content	A1	A2	А3	A3.5	A1	A2	А3
14 mm	98.5	1 pcs	60022899	60022900	60022901	60022902	60022907	60022908	60022909
18 mm	98.5	1 pcs	60022886	60022887	60022888	60022889	60022894	60022895	60022896





60022908





BRILLIANT COMPONEER®

38

The new BRILLIANT COMPONEER veneering system is comprised of prefabricated composite veneers for aesthetic anterior restorations in a single visit. The BRILLIANT COMPONEER shells are based on the innovative BRILLIANT EverGlow sub-micron filler technology, which results in exceptional polishability and long-lasting gloss retention. The composite veneers are extremely thin and allow a minimal-invasive treatment to open up new economic possibilities for both dentists and patients.

BRILLIANT COMPONEER KIT Syringes

Direct Composite Veneering System - A smile made Brilliant!

The new BRILLIANT COMPONEER veneering system is comprised of prefabricated composite veneers for anterior restoration. Due to the smooth consistency of BRILLIANT EverGlow and excellent shade matching with BRILLIANT COMPONEER, the system is intuitive and easy to use. As only one visit is necessary, the innovative Composite Veneering System BRILLIANT COMPONEER opens up new economic possibilities for both dentists and patients.

REF		Item		
KLI		item		
6002321	4	Intro Kit Syringes	Accessor	ies
			25×	COMPONEER Brushes
			1 ×	COMPONEER Contour Guide - Lower Medium
	107 W		1 ×	Placer
1	RILLIANT COMPONELS		1×	COMPONEER Contour Guide - Lower Small
/ Br		A COPLEME	1 × 2.5 ml	Etchant Gel S
\		and the same of th	1×	Modelling Instrument MBC
\			2×	Placer Adapter
1			10×	Placer Stamp
1	01/016		1×	Silicone polisher 2303RA
		F	1×	Silicone polisher 2311RA
			1×	Silicone polisher 2403RA
			1×	Silicone polisher 2411RA
			1 ×	COMPONEER Contour Guide - Upper Large
			1 ×	COMPONEER Contour Guide - Upper Medium
			1×	COMPONEER Contour Guide - Upper Small
			1 × 3 g	BRILLIANT EverGlow A1/B1
1×6 pcs	Veneers Up	pper Jaw, Bleach Translucent, M, 11-13/21-23	1 × 3 g	BRILLIANT EverGlow A2/B2
1×6 pcs	Veneers Up	pper Jaw, Bleach Translucent, L, 11-13/21-23	1 × 3 g	BRILLIANT EverGlow Translucent (Trans)
1×6 pcs	Veneers Up	pper Jaw, Translucent, M, 11-13/21-23	1×3g	BRILLIANT EverGlow Bleach Translucent (BL Trans)
1×6 pcs	Veneers Up	pper Jaw, Translucent, L, 11-13/21-23	$1 \times 5 \text{ ml}$	ONE COAT 7 UNIVERSAL
1×2 pcs	Veneers Up	pper Jaw Premolars, Bleach Translucent, L, 14/24	1 × 20 pcs	Application needles
1×2 pcs	Veneers Up	pper Jaw Premolars, Translucent, L, 14/24	1 ×	BRILLIANT COMPONEER Shade Guide

BRILLIANT COMPONEER KIT Tips

REF	Item	
60023215	Intro Kit Tips	Accessories
	IMIT CONTROLL	25 × COMPONEER Brushes 1 × COMPONEER Contour Guide - Lower Medium 1 × Placer 1 × COMPONEER Contour Guide - Lower Small 1 × 2.5 ml Etchant Gel S 1 × Modelling Instrument MBC 2 × Placer Adapter 10 × Placer Stamp 1 × Silicone polisher 2303RA 1 × Silicone polisher 2311RA 1 × Silicone polisher 2403RA 1 × Silicone polisher 2403RA 1 × COMPONEER Contour Guide - Upper Large 1 × COMPONEER Contour Guide - Upper Medium 1 × COMPONEER Contour Guide - Upper Small
1 × 6 pcs 1 × 6 pcs 1 × 6 pcs 1 × 6 pcs 1 × 2 pcs 1 × 2 pcs	Veneers Upper Jaw, Bleach Translucent, M, 11-13/21-23 Veneers Upper Jaw, Bleach Translucent, L, 11-13/21-23 Veneers Upper Jaw, Translucent, M, 11-13/21-23 Veneers Upper Jaw, Translucent, L, 11-13/21-23 Veneers Upper Jaw Premolars, Bleach Translucent, L, 14/24 Veneers Upper Jaw Premolars, Translucent, L, 14/24	8 × 0.2 g BRILLIANT EverGlow A1/B1 8 × 0.2 g BRILLIANT EverGlow A2/B2 8 × 0.2 g BRILLIANT EverGlow Translucent (Trans) 8 × 0.2 g BRILLIANT EverGlow Bleach Translucent (BL Trans) 1 × 5 ml ONE COAT 7 UNIVERSAL 1 × 20 pcs Application needles 1 × BRILLIANT COMPONEER Shade Guide



BRILLIANT COMPONEER Sets

Content				REF						
				Bleach Tr	anslucent			Trans	lucent	
			S	М	L	XL	S	М	L	XL
BRILLIANT COMPONEER Set Upper Jaw	6 pcs	11, 12, 13, 21, 22, 23	60023198	60023200	60023201	60023203	60023204	60023206	60023207	60023209
BRILLIANT COMPONEER Set Upper Jaw Premolars	2 pcs	14, 24	60023194	-	60023195	-	60023196	-	60023197	-
BRILLIANT COMPONEER Set Lower Jaw	6 pcs	31, 32, 33, 41, 42, 43	60023210	60023211	-	-	60023212	60023213	-	-







BRILLIANT COMPONEER Refills

Content			REF							
				Bleach Tr	anslucent			Trans	lucent	
			S	М	L	XL	S	М	L	XL
Upper Jaw	11	1 pcs	60023122	60023128	60023134	60023140	60023158	60023164	60023170	60023176
	12	1 pcs	60023123	60023129	60023135	60023141	60023159	60023165	60023171	60023177
	13	1 pcs	60023124	60023130	60023136	60023142	60023160	60023166	60023172	60023178
	21	1 pcs	60023125	60023131	60023137	60023143	60023161	60023167	60023173	60023179
	22	1 pcs	60023126	60023132	60023138	60023144	60023162	60023168	60023174	60023180
	23	1 pcs	60023127	60023133	60023139	60023145	60023163	60023169	60023175	60023181
Lower Jaw	31	1 pcs	60023146	60023152			60023182	60023188		
	32	1 pcs	60023147	60023153			60023183	60023189		
	33	1 pcs	60023148	60023154			60023184	60023190		
	41	1 pcs	60023149	60023155	-	-	60023185	60023191	-	_
	42	1 pcs	60023150	60023156			60023186	60023192		
	43	1 pcs	60023151	60023157			60023187	60023193		





SPOTLIGHT



SoloCem[®] 43 DuoCem[®] 42

SoloCem is a self-adhesive resin cement that ensures reliable, tight restoration, in concert with a particularly low level of shrinkage. This dual-curing cement achieves outstanding adhesion values that guarantee durable stability on the enamel and dentine without bonding.

DuoCem is a dual curing resin cement in 5ml-automix-syringes, which is available in Dentin, Translucent and White Opaque. Thanks to the excellent technical values and the easy handling DuoCem is perfect for daily work in a dental practice.

DuoCem®



60014381



60014382

Dual curing resin cement.

Coltène cements have been proving their worth for many years, convincing dentists the world over. DuoCem thus builds on the strengths of these recognised brands, combining them with the benefits of an automix syringe. This combination results in an all-purpose resin cement that offers excellent technical characteristics and can be applied with ease. Thus making DuoCem the ideal cement for routine bonding tasks in the dental practice

- Strong adhesion and properties similar to composites
- All-purpose material for a variety of indications
- Convenient flow properties and ideal processing times
- Easy removal of excess material

REF	Item
60014381	Intro Kit 1 × 5 ml DuoCem Translucent 1 × 5 ml DuoCem Dentin 1 × 10 pcs Mixing Tip Short Fine 1 × 10 pcs Mixing Tip Short Super Fine
60020092	Trans Refill 1 × 5 ml DuoCem Translucent 1 × 5 pcs Mixing Tip Short Fine 1 × 5 pcs Mixing Tip Short Super Fine
60014382	Dentin Refill 1 × 5 ml DuoCem Dentin 1 × 5 pcs Mixing Tip Short Fine 1 × 5 pcs Mixing Tip Short Super Fine
60014383	White Opaque Refill 1 × 5 ml DuoCem White 1 × 5 pcs Mixing Tip Short Fine 1 × 5 pcs Mixing Tip Short Super Fine



SoloCem®



60013986



60014062

Self-adhesive dual curing resin cement.

Secure cementing does not depend on the number of working steps. Self-adhesive SoloCem achieves reliable values - with no additional bonding. This is the result of monomers contained and the composite-like formulation. The convenient automix syringe and the ready-to-use mixing tips offer a simple and time-saving application of SoloCem.

- A true self-adhesive cement: Etching, bonding and cementation in one step
- Long-lasting cementation without requiring a separate adhesive
- Fast and easy handling
- Provides high bond strength on a variety of materials
- Unique combination of MDP & 4-META

REF	ltem		
60013986	Intro Kit		
	$1 \times 5 \text{ ml}$	SoloCem Translucent	
	$1 \times 5 \text{ ml}$	SoloCem Dentin	
	$1 \times 10 \text{ pcs}$	Mixing Tip Short Fine	
	1 × 10 pcs	Mixing Tip Short Super Fine	
60019840	Trans Re	fill	
	$1 \times 5 \text{ ml}$	SoloCem Translucent	
	1×5 pcs	Mixing Tip Short Fine	
	1 × 5 pcs	Mixing Tip Short Super Fine	
60014062	Dentin R	efill	
	$1 \times 5 \text{ ml}$	SoloCem Dentin	
	1×5 pcs	Mixing Tip Short Fine	
	1 × 5 pcs	Mixing Tip Short Super Fine	
60014061	White O	paque Refill	
	$1 \times 5 \text{ ml}$	SoloCem White Opaque	
	1×5 pcs	Mixing Tip Short Fine	
	1×5 pcs	Mixing Tip Short Super Fine	

SoloCem® & DuoCem® Accessories



6759

Accessories

REF	Item
6759	Mixing Tip Short Super Fine, 40 pcs
6747	Mixing Tip Short Fine, 40 pcs









ParaCore®

46

ParaCore simplifies the post & core restorative technique with its ability to be used as a 3-in-1 material for post cementation, core build-ups and crown & bridge cementation. Using one material for cementation and core build-ups provides an optimal «monoblock bond interface» between dentin, post and crown.

ParaCore®

3 Indications - 1 Material

ParaCore simplifies the post & core restorative technique with its ability to be used as a 3-in-1 material for post cementation, core build-ups and crown & bridge cementation. Using one material for cementation and core build-ups provides an optimal «monoblock bond interface» between dentin, post and crown.

- One Material for 3 Key Indications (post cementation, core-build ups, crown & bridge cementation)
- Optimal monoblock bond interface between the dentin-post-crown
- Superior bond strength and durability
- Stable consistency







ParaCore® 5 ml

Dual curing core material and resin cement in 5 ml automix-syringe.



5885



60013753

REF	ltem	
5885	Intro Kit	
	$1 \times 5 \text{ ml}$	ParaCore Dentin
	$1 \times 5 \text{ ml}$	ParaCore White
	$1 \times 5 \text{ ml}$	ParaCore Dentin SLOW
	$1 \times 5 \text{ ml}$	ParaCore Trans SLOW
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B
	1 × 20 pcs	Mixing Tip Short Super Fine
	$1 \times 50 \text{ pcs}$	Brushes
	1 ×	Mixing well
	1 ×	Brushholder white
60013753	Slow Intr	o Kit
	$2 \times 5 \text{ ml}$	ParaCore Dentin SLOW
	$2 \times 5 \text{ ml}$	ParaCore Trans SLOW
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B
	1 × 20 pcs	Mixing Tip Short Super Fine
	1 × 50 pcs	Brushes
	1 ×	Mixing well
	1 ×	Brushholder white



ParaCore® 5 ml Refills



60011391

REF	ltem	
5886	White Re	fill ParaCore White
	1 × 20 pcs	Mixing Tip Short Super Fine
5887	Dentin Re	efill
	2 × 5 ml	ParaCore Dentin
	1 × 20 pcs	Mixing Tip Short Super Fine
60011391	Slow Der	ntin Refill
	$2 \times 5 \text{ ml}$	ParaCore Dentin SLOW
	1 × 20 pcs	Mixing Tip Short Super Fine
60011392	Slow Trar	ns Refill
	$2 \times 5 \text{ ml}$	ParaCore Trans SLOW
	$1 \times 20 \text{ pcs}$	Mixing Tip Short Super Fine

ParaCore® 25 ml



5853




5854

Dual curing core build-up material in 25 ml cartridge.

REF	ltem					
1500	Dentin Ir	ntro Kit with Dispenser				
	$1 \times 25 \text{ ml}$	ParaCore Dentin				
	1 ×	COLTENE Dispenser 25ml				
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B				
	$1 \times 20 \text{ pcs}$	Mixing Tip Yellow				
	1 × 20 pcs	Oral Tip Yellow				
	1 × 50 pcs	Brushes				
	1 ×	Mixing well				
	1 ×	Brushholder white				
1501	Dentin Ir	ntro Kit				
	$1 \times 25 \text{ ml}$	ParaCore Dentin				
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B				
	$1 \times 20 \text{ pcs}$	Mixing Tip Yellow				
	$1 \times 20 \text{ pcs}$	Oral Tip Yellow				
	$1 \times 50 \text{ pcs}$	Brushes				
	1 ×	Mixing well				
	1 ×	Brushholder white				
1502	White Intro Kit with Dispenser					
	$1 \times 25 \text{ ml}$	ParaCore White				
	1 ×	COLTENE Dispenser 25ml				
	$1 \times 3 \text{ ml}$	ParaBond Non-Rinse Conditioner				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive A				
	$1 \times 3 \text{ ml}$	ParaBond Adhesive B				
	1 × 20 pcs	Mixing Tip Yellow				
	1 × 20 pcs	Oral Tip Yellow				
	1×50 pcs	Brushes				
	1 ×	Mixing well				
	1 ×	Brushholder white				
1503	White Int	tro Kit				
	1 × 25 ml	ParaCore White				
	1 × 3 ml	ParaBond Non-Rinse Conditioner				
	1 × 3 ml	ParaBond Adhesive A				
	1 × 3 ml	ParaBond Adhesive B				
	1 × 20 pcs	Mixing Tip Yellow				
	1 × 20 pcs	Oral Tip Yellow				
	1 × 50 pcs	Brushes				
	1×	Mixing well				
	1×	Brushholder white				

ParaCore® 25 ml Refills

REF	ltem	
5853	Dentin Re	fill
	$1 \times 25 \text{ ml}$	ParaCore Dentin
	$1 \times 20 \text{ pcs}$	Mixing Tip Yellow
		Oral Tip Yellow
5854	White Ref	
	$1 \times 25 \text{ ml}$	ParaCore White
	1 × 20 pcs	Mixing Tip Yellow
	$1 \times 20 \text{ pcs}$	Oral Tip Yellow



Accessory Core

REF	Item
4470	Dispenser Cartridge, 1 pcs
6550	Mixing Tip Yellow, 40 pcs
6555	Oral Tip Yellow, 100 pcs
6747	Mixing Tip Short Fine, 40 pcs
6759	Mixing Tip Short Super Fine, 40 pcs



SPOTLIGHT



Cool Temp® NATURAL

53

Eugenol-free material for temporary fillings. For over 40 years, Coltosol F proofs its reliability on a daily basis by delivering satisfying results all over the world.

Cool Temp NATURAL is a temporary crowns and bridge material based on bis-acryl composite.

Coltosol® F



5906



5930

Coltosol F is a zinc oxide / zinc sulphate-based eugenol-free filling material for temporary fillings and temporary root canal seals. Available in universal white.

- High marginal seal due to slight expansion
- Contains zinc oxide and fluoride
- Eugenol-free and radiopaque
- Self-curing under the action of saliva
- No mixing, easy to remove

REF	Item
5906	Single Pack Jar, 1×38 g
5911	Eco Pack Jar, 3 × 38 g
5930	Syringe, 5 × 8 g
50003797	Spindle for Syringe, 1 pcs

DuoTEMP®



5831

Dual-curing temporary filling material for temporary fillings, temporary endodontic sealing and temporary fillings for inlay- and onlay preparations. The material is based on zinc oxide / zinc sulphate and is available in universal white.

- Dual-curing (light and saliva)
- High marginal seal due to slight expansion
- Contains zinc oxide and fluoride
- Eugenol-free and radiopaque
- Good retention, easy removal



7090

REF	Item
5831	Single Pack Syringe, 1 \times 5 g
7090	Eco Pack Syringe, 5 × 5 g



Cool Temp® NATURAL

Temporary crown and bridge material based on bis-acryl composite, available in automix cartridges and automix syringes.

- Low polymerisation temperature
- Three fluorescent shades
- Low shrinkage
- Excellent polishability
- Easy and fast processing time

REF	ltem	
5800	Starter Ki	t
	1 × 85 g	Cool Temp NATURAL A2
	1 × 85 g	Cool Temp NATURAL A3.5
	1 × 2.3 g	SYNERGY D6 Flow Dentin A2/B2
	1 × 2.3 g	SYNERGY D6 Flow Dentin A3.5/B3
	2 × 10 pcs	Mixing Tip Blue pointed
	4×	Flow Application Needles 0.8 mm

Cool Temp® NATURAL Refills



5805



5824

REF	Item
5815	A1 Refill Cartridge 1 × 85 g Cool Temp NATURAL A1 1 × 10 pcs Mixing Tip Blue pointed
5805	A2 Refill Cartridge 1 × 85 g Cool Temp NATURAL A2 1 × 10 pcs Mixing Tip Blue pointed
5810	A3.5 Refill Cartridge 1 × 85 g Cool Temp NATURAL A3.5 1 × 10 pcs Mixing Tip Blue pointed
5820	A1 Refill Syringe 2 × 8.5 g Cool Temp NATURAL A1 8 × Mixing Tip Brown pointed
5824	A2 Refill Syringe 2 × 8.5 g Cool Temp NATURAL A2 8 × Mixing Tip Brown pointed
5825	A3.5 Refill Syringe 2 × 8.5 g Cool Temp NATURAL A3.5 8 × Mixing Tip Brown pointed

Cool Temp® NATURAL Accessories

Accessories

REF	Item
6748	Mixing Tip Brown pointed, 40 pcs
7774	Mixing Tip Blue pointed, 40 pcs
6900	Dispenser 50/75 ml, 1 pcs

53

TempoSIL® 2



Secure adhesion to natural tooth structure, core build up material and metal alloys as well

Secure adhesion to natural tooth structure, core build up material and metal alloys as well
as zirconia

TempoSIL 2 is a temporary, radiopaque A-Silicone cement for cementing temporary restorations or temporary cementing final restorations. Available in two shades: White and Dentin.

- Superior marginal integrity
- Fast setting time
- Eugenol-free
- Simple application as well as easy removal



3686


6730

REF	Item			
3686	Trial Pack 1 × 5 ml TempoSIL 2 Dentin 1 × 5 ml TempoSIL 2 White 1 × 15 pcs Mixing Tip Brown Long Fine			
6727	Dentin Intro Kit 2 × 5 ml TempoSIL 2 Dentin 1 × 15 pcs Mixing Tip Brown Long Fine			
6729	White Intro Kit 2 × 5 ml TempoSIL 2 White 1 × 15 pcs Mixing Tip Brown Long Fine			
6728	Dentin Refill, 4 × 5 ml			
6730	White Refill, 4 × 5 ml			

REF	Accessories
6749	Mixing Tip Brown Long Fine, 40 pcs



NOTES





The S.P.E.C. 3 curing light defines performance without compromise. Chair time can be effectively reduced due to rapid and deep polymerization on various light cured materials.

The Coltolux LED is a simple yet reliable LED curing light. The cordless and lightweight pen-style design offers good handling properties for the daily use in the practice.

S.P.E.C.® 3



60013942

S.P.E.C. 3 meets the highest demands of performance in the dental practice. The stable, light-weight aluminium handpiece and the 360° rotating light guide provide convenient handling. Three different curing modes provide performance without compromise, enabling variable use to meet the individual requirements of dentists. A high performance, interchangeable lithium polymer battery provides over 300 10s cures between recharges.

- Cures 2 mm deep in one second
- Ergonomic lightweight design
- Three distinct power modes

REF	ltem				
60013942	Kit Europe				
	1×	Handpiece			
	1×	Charger Base			
	1×	Power Supply Plug Type C			
	1×	Clear Light Guide 8 mm Turbo-Tip			
	1 × 2 pcs	Light Shield			
	1 × 100 pcs	Barrier Sleeve			
	1 ×	Eye Shield			
60013943	Kit United	d Kingdom			
	1×	Handpiece			
	1×	Charger Base			
	1×	Power Supply Plug Type G			
	1×	Clear Light Guide 8 mm Turbo-Tip			
	1 × 2 pcs	Light Shield			
	1 × 100 pcs	Barrier Sleeve			
	1 ×	Eye Shield			

S.P.E.C.® 3 Accessories



60019327



60013949

REF	Item
60013950	Clear Light Guide, Turbo-Tip 8 mm, 1 pcs
60013951	Clear Light Guide, 11 mm, 1 pcs
60019326	Black Light Guide, Turbo-Tip 8 mm, 1 pcs
60019327	Black Light Guide, 11 mm, 1 pcs
60013949	Barrier Sleeves, 100 pcs
60013948	Light Shield 8 x 11 mm, 20 pcs
60014360	Light Shield 11 x 11 mm, 20 pcs
60014444	Eye Shield, 1 pcs
60013952	Battery Module, 1 pcs
60013953	Charger Base, 1 pcs
60013955	Power Supply, 1 pcs



Coltolux® LED



sion and eliminates the need for active cooling.

- Dome shaped lens concentrates light output on the restoration
- Audible signal every ten seconds and automatic power saving sleep mode after 20 seconds of operation

The Coltolux LED combines a slim curing probe with high output in a cordless lightweight pen-style design. The Coltolux LED's heat sink technology reduces the amount of heat emis-

- Hundreds of cures with a fully charged battery
- LED light output matches the absorption spectrum of Camphorquinone (CQ) Photoinitiators

REF	ltem		
C7970	Kit 230/240V		
	1 ×	Handpiece	
	1 ×	Charger Base	
	1 ×	Power Supply	
	1×3 pcs	Light Shield	
	1×5 pcs	Light Lense	
	1 × 100 pcs	Barrier Sleeve	

C7970

Coltolux® LED Accessories

Accessories

REF	Item
C7972	Light Lenses, Press-on, 25 pcs
C7990	Light Shield, Snap-on, 25 pcs
C7957	Barrier Sleeves, 1000 pcs
C7970RK2	Charger Base 230/240V, 1 pcs

SPOTLIGHT



MAX Pins have been specifically designed to significantly reduce the stresses on the tooth associated with pin insertion and function. Yet, at the same time, retention of the restoration is optimized. MAX Pin is a complete redesign of pin technology – for safer, predictable restorative dentistry.

TMS pins are easy to use, self-threading retention pins. The well thought out construction of TMS, TMS LINK and TMS LINK PLUS is the result of decades of experience and an uncompromising effort to meet restorative needs.

MAX® Pin

Self-threading parapulpal retention pins with shoulder stop and rounded head design for good adaptation. Corrosion-free, titanium alloy. Self-alignment during placement.

Material	Color code	Ø mm	Pin Design	Sales Unit	Self-shearing	REF	Content
Titanium Alloy	Blue	0.425	Single	Complete-Kit	Yes	M61	25 Pins, 2 MAX Drills, 1 Hand Driver
				Bulk Kit	Yes	M63	60 Pins
	Purple	0.525	Single	Complete-Kit	Yes	M71	25 Pins, 2 MAX Drills, 1 Hand Driver
				Bulk Kit	Yes	M73	60 Pins



MAX® Pin Accessories

REF	
M91	Drills Blue, Depth-Limiting, 5 pcs
M92	Drills Purple, Depth-Limiting, 5 pcs

TMS

For permanent retention in the most demanding indications

TMS pins are easy to use, self-threading retention pins. The well thought out construction of TMS, TMS LINK and TMS LINK PLUS is the result of decades of experience and an uncompromising effort to meet restorative needs.



TMS LINK PLUS®

TMS LINK PLUS pins offer smooth self-alignment, tapered pin end, shoulder stop and automatic shear-points. They are available as single pins or as convenient 2-in-1 pin system for a time saving placement. The pins can be mechanically installed.

Material	Color code	Ømm	Pin Design	Sales Unit	Self-shearing	REF	Content
Titanium Alloy	Minim Silver	0.525	Single	Complete-Kit	Yes	EL82124	24 Pins, 1 KODEX Drill K92, 1 Hand Driver
				Bulk Kit	Yes	EL82260	60 Pins
			2-in-1	Complete-Kit	Yes	EL84124	24 Pins, 2 KODEX Drills K92, 1 Hand Driver
				Bulk Kit	Yes	EL84260	60 Pins
Gold-Plated Stainless Steel	Minim Silver	0.525	Single	Complete-Kit	Yes	EL72124	24 Pins, 1 KODEX Drill K92, 1 Hand Driver
				Bulk Kit	Yes	EL72260	60 Pins
			2-in-1	Complete-Kit	Yes	EL73124	24 Pins, 2 KODEX Drills K92, 1 Hand Driver
				Bulk Kit	Yes	EL73260	60 Pins
	Regular Gold	0.675	Single	Complete-Kit	Yes	EL74124	24 Pins, 1 KODEX Drill K96, 1 Hand Driver
				Bulk Kit	Yes	EL74260	60 Pins
			2-in-1	Complete-Kit	Yes	EL75124	24 Pins, 2 KODEX Drills K96, 1 Hand Driver
				Bulk Kit	Yes	EL75260	60 Pins



EL72124

TMS LINK® Series

The 2nd generation TMS pins can be mechanically installed and offer a ball joint link to the shank for a smooth self-alignment as well as a tapered pin end for an easier insertion. The TMS LINK pin is available as single pin or as convenient 2-in-1 system.

Material	Color code	Ø mm	Pin Design	Sales Unit	Self-shearing	REF	Content
Titanium Alloy	Minikin Red	0.425	Single	Complete-Kit	Yes	L811	20 Pins, 1 KODEX Drill K91, 1 Hand Driver
				Bulk Kit	Yes	L812	50 Pins
Gold-Plated	Minuta Pink	0.350	Single	Complete-Kit	Yes	L501	20 Pins, 1 KODEX Drill K89, 1 Hand Driver
Stainless Steel				Bulk Kit	Yes	L502	50 Pins
	Minikin Red	0.425	Single	Complete-Kit	Yes	L511	20 Pins, 1 KODEX Drill K91, 1 Hand Driver
				Bulk Kit	Yes	L512	50 Pins
	Minim Silver	0.525	Single	Complete-Kit	Yes	L521	20 Pins, 1 KODEX Drill K92, 1 Hand Driver
				Bulk Kit	Yes	L522	50 Pins
			2-in-1	Complete-Kit	Yes	L531	20 Pins, 2 KODEX Drills K92, 1 Hand Driver
				Bulk Kit	Yes	L532	50 Pins
	Regular Gold	0.675	Single	Complete-Kit	Yes	L541	20 Pins, 1 KODEX Drill K96, 1 Hand Driver
				Bulk Kit	Yes	L542	50 Pins
			2-in-1	Complete-Kit	Yes	L551	20 Pins, 2 KODEX Drills K96, 1 Hand Driver
				Bulk Kit	Yes	L552	50 Pins



L501



TMS®

The 1st generation of TMS was conceived for manual insertion and is available as single pin or as convenient 2-in-1 system.

Material	Color code	Ø mm	Pin Design	Sales Unit	Self-shearing	REF	Content
Gold-Plated Stainless Steel	Minuta Pink	0.350	Single	Complete-Kit	Yes	T111	20 Pins, 1 KODEX Drill K89, 1 Hand Wrench
				Refill	Yes	T112	50 Pins
	Minikin Red	0.425	Single	Complete-Kit	Yes	T107	20 Pins, 1 KODEX Drill K91, 1 Hand Wrench
				Refill	Yes	T108	50 Pins
				Bulk Kit	Yes	T109	100 Pins
	Minim Silver	0.525	Single	Complete-Kit	No	T04	20 Pins, 1 KODEX Drill K93, 1 Hand Wrench
				Refill	No	T05	50 Pins
				Bulk Kit	No	T06	100 Pins
			2-in-1	Complete-Kit	Yes	T10	20 Pins, 1 KODEX Drill K92, 1 Hand Wrench
				Refill	Yes	T11	50 Pins
				Bulk Kit	Yes	T12	100 Pins
	Regular	0.675	Single	Complete-Kit	Yes	T13	20 Pins, 1 KODEX Drill K97, 1 Hand Wrench
	Gold				No	T01	20 Pins, 1 KODEX Drill K97, 1 Hand Wrench
				Refill	No	T02	50 Pins
					Yes	T14	50 Pins
				Bulk Kit	No	T03	100 Pins
					Yes	T15	100 Pins
			2-in-1	Complete-Kit	Yes	T07	20 Pins, 1 KODEX Drill K96, 1 Hand Wrench
				Refill	Yes	T08	50 Pins
				Bulk Kit	Yes	T09	100 Pins





T02

KODEX®-Drills

KODEX-Drills complete the TMS, TMS LINK and TMS LINK PLUS retention pin systems. A colour coding helps to select the appropriate size, thereby ensuring an optimal preparation of the canal.

Color code	Ø mm	Drill length	Depth-limiting	Shank length	REF	Content
Minuta Pink	0.350	1.3 mm	Yes	Long	K89	6 pcs
Minikin Red	0.425	1.5 mm	Yes	Short	K90	6 pcs
		1.5 mm	Yes	Long	K91	6 pcs
Minim Silver	0.525	2.0 mm	Yes	Regular	K92	6 pcs
		4.0 mm	No	Regular	K93	6 pcs
Regular Gold	0.675	2.0 mm	Yes	Regular	K96	6 pcs
		5.0 mm	No	Regular	K97	6 pcs



Accessories Parapulpal Pins

Item
Hand Wrench, Minuta Pink, 1 pcs
Hand Wrench, Regular Gold & Minikin Red, 1 pcs
Hand Wrench, Minim Silver, 1 pcs
Bending Tool, Minim Silver, 1 pcs
Bending Tool, Regular Gold, 1 pcs
Universal Hand Driver, 1 pcs



LOCATIONS

BR Vigodent SA Indústria e Comércio

Rua Pesqueira 26 Bonsucesso Rio de Janeiro CEP 21041-150 / Brazil T +5521 3865 5600 F +5521 2280 8411 sac@vigodent.com.br

CN

Coltène/Whaledent China Trading CO. Ltd. Rm. 1103/1104, Bldg. B2, Wu Dong Plaza #9 Chegongzhuang Street, West District Beijing 100044 / P.R. China T +86 10 6212 3251 F +86 10 6212 3252 info.cn@coltene.com

ES Coltène Iberia S.L.U.

c/Doctor Fleming, 3 - Bajo A 28220 Majadahonda, Madrid / Spain PT T +34 91 63 86 492 F +34 91 63 86 279 info.ib@coltene.com

IN Coltène/Whaledent Private Ltd.

106/107, Hallmark, Vasant Oscar L.B.S Marg, Mulund West Mumbai 400080, Maharashtra / India T +91 22 25923626 F +91 22 25623165 info.in@coltene.com

JΡ

Coltène Japan LLC 2-2-28 2nd Floor, Nishikicho 2-2-2-8 210 Floor, NISHINCHO 190-0022 Tachikawa Tokyo / Japan T +81 42 595 6945 F +81 42 595 6946 info.jp@coltene.com

US Coltène/Whaledent Inc. 235 Ascot Parkway

Cuyahoga Falls, Ohio 44223 / USA T+1 330 916 8800 F +1 330 916 7077 info.us@coltene.com

CH Coltène/Whaledent AG

Feldwiesenstrasse 20 9450 Altstätten / Switzerland T +41 71 757 53 00 F +41 71 757 53 01 info.ch@coltene.com

Coltène/Whaledent GmbH & Co. KG Raiffeisenstraße 30 DE

89129 Langenau / Germany T +49 7345 805 0 F +49 7345 805 201 info.de@coltene.com

FR Coltène/Whaledent S.a.r.l.


Synergie Park 3 avenue Louis Néel 59260 Lezennes / France T +33 320 64 64 00 F +33 320 64 64 03 info.fr@coltene.com

IT Coltène Italia S.r.l

c/o Polo Tecnologico della Brianza Via Alessandro Volta, 94 - Cond. Bianchi 20832 Desio (MB) / Italy T +39 0362 62 67 32 F +39 0362 63 25 80 info.it@coltene.com

UK Coltène/Whaledent Ltd.

The President Suite Kendal House, Burgess Hill Victoria Way West Sussex, RH15 9NF / UK T +44 1444 235486 F +44 1444 870640 info.uk@coltene.com



© COLTENE - www.coltene.com

