

Now with the new generation of smart endodontic devices

Catalogue  
**2015**





**Access Preparation, Pulp Extraction**

**4**

**NiTi Systems and Motors**

**6**

**Conventional Preparation**

**14**

**Length Determination**

**20**

**Ultrasound and Irrigation**

**21**

**Obturation**

**23**

**Posts**

**33**

**Boxes and Accessories**

**35**

**Your VDW Contacts**

**37**

### VDW®STERILE Gates Enlargers

#### For enlargement of the coronal third

- Stainless steel
- Recommended speed 800-1,200 rpm
- Non-cutting tip
- Blister pack of 6 pieces



Size	19 mm ●	19 mm ▲
1 / 50	0380 019 001	0180 019 001
2 / 70	0380 019 002	0180 019 002
3 / 90	0380 019 003	0180 019 003
4 / 110	0380 019 004	0180 019 004
5 / 130	0380 019 005	0180 019 005
6 / 150	0380 019 006	0180 019 006
1-6 / 50-150	0380 019 300	0180 019 300



### Peeso Enlargers

#### For enlargement of the coronal third

- Stainless steel
- Recommended speed 800-1,200 rpm
- Non-cutting tip
- Blister pack of 6 pieces



Size	19 mm
1 / 70	0183 019 001
2 / 90	0183 019 002
3 / 110	0183 019 003
4 / 130	0183 019 004
5 / 150	0183 019 005
6 / 170	0183 019 006
1-6 / 70-170	0183 019 300



### Beutelrock Reamers

#### For enlarging coronal and straight root canal sections

- Stainless steel, cylindrical
- Recommended speed 450-800 rpm
- Blister pack of 6 pieces



Size		18 mm
1 / 35	○	0944 018 001
2 / 45	●	0944 018 002
3 / 60	●	0944 018 003
4 / 75	●	0944 018 004
5 / 90	●	0944 018 005
6 / 105	●	0944 018 006
1-6 / 35-105		0944 018 300



## Beutelrock Enlargers

### For flaring of root canal orifices

- Stainless steel
- Recommended speed 800-1,200 rpm
- Blister pack of 6 pieces



Size		11 mm
1 / 90	○	0947 011 001
2 / 100	●	0947 011 002
3 / 120	●	0947 011 003
4 / 140	●	0947 011 004
5 / 160	●	0947 011 005
6 / 180	●	0947 011 006
1-6 / 90-180		0947 011 300



## VDW®STERILE Barbed Broaches

### For easy and safe pulp tissue removal

Flexible barbs are arranged spirally on the cone shaped core. The size corresponds to 1/2 of the core diameter.

- Stainless steel, CC-Cord® Handle
- 150 pieces, in blister packs



Size		21 mm
XXXXF / 0	●	0333 021 410
XXXF / 1	○	0333 021 420
XXF / 2	●	0333 021 430
XF / 3	●	0333 021 440
F / 4	●	0333 021 450
M / 5	●	0333 021 460
C / 6	●	0333 021 470
XXF-F / 2-4		0333 021 480



## RECIPROC® Instruments



### System designed specifically for root canal preparation in reciprocation

- Root canal preparation with only one instrument
- Retreatment of gutta-percha filling material and carrier based obturators with R25
- Single-use convenience; sterile instruments, individually packed
- Three instrument sizes for narrow, medium and wide canals

STERILE

Blister pack of 6 instruments		21 mm	25 mm	31 mm	
R25	●	0212 021 025	0212 025 025	0212 031 025	
R40	●	0212 021 040	0212 025 040	0212 031 040	
R50	●	0212 021 050	0212 025 050	0212 031 050	
3 x R40, 3 x R50 each		0212 021 233	0212 025 233	0212 031 233	
Blister pack of 4 instruments		21 mm	25 mm	31 mm	
R25	●	0012 021 025	0012 025 025	0012 031 025	
2 x R25, 1 x R40, 1 x R50 each		0012 021 200	0012 025 200	0012 031 200	

### RECIPROC® System Kit 24

- 2 blister packs of 6 instruments R25, sterile
- 1 blister pack of 6 instruments R40, sterile
- 1 blister pack of 6 instruments R50, sterile
- Training Model
- RECIPROC® Paper Points, assorted, sizes R25, R40, R50, sterile
- RECIPROC® Gutta-Percha, assorted, sizes R25, R40, R50
- User card, brochure, directions for use

#### REF

1213 025 000  
25 mm




# Mtwo® Instruments



## The efficient NiTi system

- One sequence for all types of root canal
- Highly flexible due to small instrument core and s-shaped cross-section



Mtwo® is licensed by sweden & martina

## Basic Sequence Blister pack of 6 instruments



Size/Taper	Mtwo® working part 16 mm		Mtwo® working part 21 mm			
	21 mm	25 mm	25 mm	31 mm		
10/.04	●	0234 021 010	0234 025 010	1234 025 010	1234 031 010	
15/.05	○	0235 021 015	0235 025 015	1235 025 015	1235 031 015	
20/.06	●	0236 021 020	0236 025 020	1236 025 020	1236 031 020	
25/.06	●	0236 021 025	0236 025 025	1236 025 025	1236 031 025	

## Shaping of Larger Canal Anatomies Blister pack of 6 instruments



Size/Taper	Mtwo® working part 16 mm		Mtwo® working part 21 mm			
	21 mm	25 mm	31 mm			
30/.05	●	0235 021 030	0235 025 030	1235 031 030		
35/.04	●	0234 021 035	0234 025 035	1234 031 035		
40/.04	●	0234 021 040	0234 025 040	1234 031 040		
45/.04	○	0234 021 045	0234 025 045	1234 031 045		
50/.04	●	0234 021 050	0234 025 050	1234 031 050		
60/.04	●	0234 021 060	0234 025 060	1234 031 060		
Size/Taper	Mtwo® working part 16 mm		Mtwo® working part 21 mm			
	21 mm	25 mm	31 mm			
25/.07	●	0237 021 025	0237 025 025	1237 031 025		
30/.06	●	0236 021 030	0236 025 030	1236 031 030		
35/.06	●	0236 021 035	0236 025 035	1236 031 035		
40/.06	●	0236 021 040	0236 025 040	1236 031 040		

## Assortments



Blister pack of 4 instruments	Mtwo® working part 16 mm		Mtwo® working part 21 mm	
	21 mm	25 mm	25 mm	31 mm
10/.04, 15/.05, 20/.06, 25/.06	0007 021 701	0007 025 711	0007 025 701	0007 031 701
30/.05, 35/.04, 40/.04, 25/.07	0007 021 702	0007 025 702	-	0007 031 702
Blister pack of 6 instruments	21 mm	25 mm	25 mm	31 mm
2 x 10/.04, 15/.05, 20/.06 each	0230 021 456	0230 025 456	-	-
2 x 25/.06, 30/.05, 35/.04, 40/.04, 25/.07	0230 021 457	0230 025 457	-	-
2 x 10/.04, 15/.05, 20/.06 each	-	-	1230 025 456	
2 x 25/.06, wp 21 mm, 30/.05, 35/.04, 40/.04, 25/.07, wp 16 mm	-	1230 025 457		-
2 x 30/.06, 35/.06, 40/.06 each	-	0236 025 203	-	-
2 x 45/.04, 50/.04, 60/.04 each	-	0234 025 246	-	-

## Mtwo® Retreatment Files

For efficient removal of filling material

- Cutting tip
- Working part 16 mm

Blister pack of 6 instruments



Size/Taper		21 mm	
R15/.05	○	0232 021 015	
R25/.05	●	0232 021 025	
3 x R15/.05, R25/.05 each		0232 021 202	

### Mtwo® Basic Kit

REF  
0238 025 000

Kit with instruments in 25 mm working length:

- 1 blister pack of Mtwo® 2 x 10/.04, 15/.05, 20/.06 each, 21 mm working part, sterile
- 1 blister pack of Mtwo® 2 x 25/.06, 21 mm working part, 1 x 30/.05, 35/.04, 40/.04, 25/.07 each, 16 mm working part, sterile
- 3 plastic training blocks
- CD-ROM and user card



### Mtwo® SystemBox

0239 000 000

- With module for easy organisation and storage of Mtwo® instruments
- Instruments not included, autoclavable



### Mtwo® Modul

0408

- Instrument module, without box
- Instruments not included




### Mtwo® Control Stickers

0491 000 000

For practical monitoring of instrument use

- 16 Control Stickers for the basic sequence
- 12 Control Stickers for additional sizes
- self-adhesive, autoclavable





## FlexMaster® Instruments

# FlexMaster®

### The established crown-down NiTi system

- Three sequences to choose from for every type of canal
- Safe use due to optimal stability of file even with low tapered files (2 %)
- Special K-type cutting edges for efficient removal of debris and dentine chips

### FlexMaster® Taper .02 Blister pack of 6 instruments



Size/Taper		21 mm	25 mm	
15/.02	○	0341 021 015	0341 025 015	
20/.02	●	0341 021 020	0341 025 020	
25/.02	●	0341 021 025	0341 025 025	
30/.02	●	0341 021 030	0341 025 030	
35/.02	●	0341 021 035	0341 025 035	
40/.02	●	0341 021 040	0341 025 040	
45/.02	○	0341 021 045	0341 025 045	
50/.02	●	-	0341 025 050	
60/.02	●	-	0341 025 060	
70/.02	●	-	0341 025 070	

### FlexMaster® Taper .04 Blister pack of 6 instruments



Size/Taper		21 mm	25 mm	31 mm	
15/.04	○	0342 021 015	0342 025 015	-	
20/.04	●	0342 021 020	0342 025 020	0342 031 020	
25/.04	●	0342 021 025	0342 025 025	0342 031 025	
30/.04	●	0342 021 030	0342 025 030	0342 031 030	
35/.04	●	0342 021 035	0342 025 035	0342 031 035	
40/.04	●	0342 021 040	0342 025 040	0342 031 040	

### FlexMaster® Taper .06 Blister pack of 6 instruments



Size/Taper		21 mm	25 mm	
15/.06	○	0343 021 015	0343 025 015	
20/.06	●	0343 021 020	0343 025 020	
25/.06	●	0343 021 025	0343 025 025	
30/.06	●	0343 021 030	0343 025 030	
35/.06	●	0343 021 035	0343 025 035	
40/.06	●	0343 021 040	0343 025 040	

### FlexMaster® Assortments Blister pack of 6 instruments



Size/Taper	21 mm	25 mm	31 mm
20/.02, 25/.02, 30/.02, 35/.02, 40/.02, 45/.02	-	0341 025 231	0341 031 231
20/.04, 25/.04, 30/.04, 20/.06, 25/.06, 30/.06	-	0344 025 232	-
2x 50/.02, 60/.02, 70/.02 each	0341 021 241	-	0341 031 241

REF

## FlexMaster® IntroFile

0357 019 000

- 11 % taper
- 19 mm working length, 9 mm active working part
- ISO size 22 at the tip



- FlexMaster® Basic Kit** \_\_\_\_\_ REF 0346 025 500
- 10 FlexMaster® instruments
    - Taper .02: sizes 20, 25, 30, 35
    - Taper .04: sizes 20, 25, 30
    - Taper .06: sizes 20, 25, 30
  - 1 IntroFile for enlargement of the canal orifice
  - 1 C-PILOT® File ISO 10
  - Product overview
  - Instruments are non-sterile




- FlexMaster® SystemBox** \_\_\_\_\_ REF 0340
- The perforated base allows efficient sterilisation with autoclave paper (REF 0191)
  - Without instruments



- FlexMaster® AccessoryBox** \_\_\_\_\_ REF 0345
- The perforated base allows efficient sterilisation with autoclave paper (REF 0191)
  - Without instruments



- FlexMaster® CombiBox** \_\_\_\_\_ REF 0445
- The perforated base allows efficient sterilisation with autoclave paper (REF 0192)
  - Without instruments



- FlexMaster® Control Stickers** \_\_\_\_\_ REF 0489
- for SystemBox
- For practical monitoring of instrument use
- Dispenser with adhesive labels, sterilisable in autoclave
  - 250 stickers



## Accessories

- FileCare® EDTA** \_\_\_\_\_ REF 1010 000 002
- Chelating agent
- Active substances: 15 % EDTA and 10 % urea peroxide gel
  - Removes smear layer
  - Improved tubular penetration of endodontic medication and filling material
  - Easier penetration of the instrument, as FileCare®EDTA decalcifies dentine by binding the calcium in an electronegative complex
  - Pulp tissue and dentinal debris are actively flushed out of the canal
- 2 x 3 ml  
REF 1010 000 002
- 5 x 3 ml  
REF 1010 000 005



## VDW.CONNECT Drive®

## VDW.CONNECT Drive®

**Cordless endo motor with iPad application**

- iPad application for rotary and reciprocating file systems
- Ready for use with RECIPROC®. No App or iPad required! Unique!
- Powerful drive up to 5 Ncm
- Ergonomic weight balance
- Contra-angle with mini head rotates 360°
- Powerful lithium-ion battery

iPad mini not included in delivery scope

	REF
<b>■ VDW.CONNECT Drive®</b> • Incl. 6:1 contra-angle and handpiece holder	1183 000 000 (for EU countries)  1184 000 000 (for non-EU countries)
<b>■ VDW.CONNECT Accessories for iPad mini</b> • Stand with case for iPad mini • 5 sleeves for iPad mini	1185 000 550
<b>■ VDW.CONNECT Drive® Set</b> • VDW.CONNECT Drive® • VDW.CONNECT Accessories for iPad mini	1183 000 500
<b>■ VDW.CONNECT Drive® Set Medium</b> • VDW.CONNECT Drive® Set + RECIPROC® System Kit 24* *Contents of the RECIPROC® System Kit 24, see page 6	1181 000 000



## Accessories for VDW.CONNECT Drive®

	REF
■ 6:1 contra-angle	1185 000 504
■ Handpiece sleeves, 100 pieces	1185 000 503
■ Sleeves for iPad mini, 100 pieces	1185 000 507
■ Charger incl. adaptors	1185 000 501
■ Contra-Angle Adaptor for Assistina	1185 000 511 (corresponds with W&H REF 06037100)

## VDW.CONNECT Application

# VDW.CONNECT App

### Compatible with iPad mini or iPad

- For wireless communication with VDW.CONNECT Drive® via Bluetooth
- Supports continuous rotary and reciprocating preparation
- Various pre-programmed file systems; settings can be changed
- Intuitive navigation

The VDW.CONNECT application is compatible with both iPad mini and iPad, equipped with Bluetooth 4.0 low energy. Your iOS version must not be lower than 8.0.

## ■ Home

- Info bar: \_\_\_\_\_
  - Direct link to the Endo Community VDWuser.net
  - Set up of individual user profiles
  - Quick access to user manual and settings
- Control panel: \_\_\_\_\_
  - Icons to select pre-set file systems
  - Icons to set up individual file sequences and instruments
  - Shortcuts to set favourites
- Status bar: \_\_\_\_\_
  - Permanent control over actual status of connectivity and battery charging capacity



## ■ Preparation

- Files with ISO colour coding \_\_\_\_\_
- Keypad \_\_\_\_\_



**Free download**  
in the App Store incl. free updates  
to extended functions in the future

Apple, the Apple logo and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

The current version of the VDW.CONNECT App may vary slightly from the screenshots shown on this page.

**VDW.GOLD®RECIPROC®****VDW.GOLD®RECIPROC®****Endo motor with integrated apex locator in reciprocating and continuous rotary mode**

- Extensive range of reciprocating and continuous rotary file systems
- Perfect interplay in the RECIPROC® mode
- Unique RECIPROC REVERSE comfort function
- Torque control with Automatic-Stop-Reverse

■ **VDW.GOLD®RECIPROC®** \_\_\_\_\_ REF  
1173 000 000

- Incl. 6:1 contra-angle

■ **VDW.GOLD®RECIPROC® + RECIPROC® System Kit** \_\_\_\_\_ REF  
1173 025 611

- 1 blister pack of 6 instruments R25, sterile
- 1 blister pack of 6 instruments assorted, 3 x R40 and R50 each, sterile
- RECIPROC® Paper Points and Gutta-Percha, assorted

**Accessories for VDW.GOLD®RECIPROC® and VDW.GOLD®**

REF

■ Contra-Angle Sleeve, 5 pieces _____	1107 000 512
■ Lip Clip, 5 pieces _____	1107 000 514
■ File Clip, 2 pieces _____	1107 000 515
■ Lip Clip Cable with Ferrite Ring _____	1177 000 516
■ File Clip Cable _____	1107 000 519
■ Charger _____	1107 000 517
■ Foot Switch _____	1107 000 510

**VDW.SILVER®RECIPROC®****VDW.SILVER®RECIPROC®****Ideal for reciprocating and classical preparation technique**

- Well-structured file library
- Perfect interplay in the RECIPROC® mode
- Unique RECIPROC REVERSE comfort function
- Torque control with Automatic-Stop-Reverse

■ **VDW.SILVER®RECIPROC®** \_\_\_\_\_ REF  
1163 000 000

- Incl. 6:1 contra-angle

■ **VDW.SILVER®RECIPROC® + RECIPROC® System Kit** \_\_\_\_\_ REF  
1163 025 611

- 1 blister pack of 6 instruments R25, sterile
- 1 blister pack of 6 instruments assorted, 3 x R40 and R50 each, sterile
- RECIPROC® Paper Points and Gutta-Percha, assorted



■ **6:1 Contra-Angle** \_\_\_\_\_ REF  
1079 000 000

- For VDW.GOLD®RECIPROC®, VDW.GOLD®, VDW.SILVER®RECIPROC®, VDW.SILVER®





### VDW®STERILE C-PILOT® Files

#### For severely curved or calcified root canals

- Stainless steel, CC+ handle
- High bending strength due to a special thermal hardening process
- Enables creation of a smooth glide path prior to every rotary preparation
- The inactive pilot tip conducts the instrument safely along the canal
- Blister pack of 6 pieces, sterile



Size		19 mm	21 mm	25 mm	
06	●	0368 019 006	0368 021 006	0368 025 006	
08	●	0368 019 008	0368 021 008	0368 025 008	
10	●	0368 019 010	0368 021 010	0368 025 010	
06-10		0368 019 210	0368 021 210	0368 025 210	
12.5	●	0368 019 012	0368 021 012	0368 025 012	
15	○	0368 019 015	0368 021 015	0368 025 015	

### VDW®STERILE K-Reamers

#### For safe root canal preparation with a reaming motion

- Stainless steel
- ISO 06-40 square cross-section for better fracture resistance
- ISO 45-140 triangular cross-section for better cutting efficiency
- Special steel alloy with uniform structure for maximum resistance to fracture and maximum flexibility
- The non-cutting tip guides the instrument safely along the canal
- Blister pack of 6 pieces



Size		21 mm	25 mm	31 mm	
06	●	0353 021 006	0353 025 006	0353 031 006	
08	●	0353 021 008	0353 025 008	0353 031 008	
10	●	0353 021 010	0353 025 010	0353 031 010	
06-10		0353 021 210	0353 025 210	0353 031 210	
15	○	0353 021 015	0353 025 015	0353 031 015	
20	●	0353 021 020	0353 025 020	0353 031 020	
25	●	0353 021 025	0353 025 025	0353 031 025	
30	●	0353 021 030	0353 025 030	0353 031 030	
35	●	0353 021 035	0353 025 035	0353 031 035	
40	●	0353 021 040	0353 025 040	0353 031 040	
15-40		0353 021 230	0353 025 230	0353 031 230	
45	○	0353 021 045	0353 025 045	0353 031 045	
50	●	0353 021 050	0353 025 050	0353 031 050	
55	●	0353 021 055	0353 025 055	0353 031 055	
60	●	0353 021 060	0353 025 060	0353 031 060	
70	●	0353 021 070	0353 025 070	0353 031 070	
80	●	0353 021 080	0353 025 080	0353 031 080	
45-80		0353 021 240	0353 025 240	0353 031 240	
90	○	0353 021 090	0353 025 090	0353 031 090	
100	●	0353 021 100	0353 025 100	0353 031 100	
110	●	0353 021 110	0353 025 110	0353 031 110	
90-110		0353 021 250	-	-	
120	●	-	0353 025 120	0353 031 120	
130	●	-	0353 025 130	0353 031 130	
140	●	-	0353 025 140	0353 031 140	
90-140		-	0353 025 260	0353 031 260	

**K-Reamers**



Single sizes	21 mm	25 mm	28 mm	31 mm
06-140	-	0053 025 XXX	-	0053 031 XXX
06-110	0053 021 XXX	-	0053 028 XXX	-
<b>Assortments</b>				
06-10	0053 021 210	0053 025 210	0053 028 210	0053 031 210
15-40	0053 021 230	0053 025 230	0053 028 230	0053 031 230
45-80	0053 021 240	0053 025 240	0053 028 240	0053 031 240
90-110	0053 021 250	0053 025 250	0053 028 250	0053 031 250
90-140	-	0053 025 260	-	0053 031 260
120-140	-	0053 025 270	-	0053 031 270

**VDW®STERILE K-Files**

For all preparation techniques with filing motion

- Stainless steel
- ISO 06-25 square cross-section for better fracture resistance
- ISO 30-140 triangular cross-section for better cutting efficiency
- Special steel alloy with uniform structure for maximum resistance to fracture and maximum flexibility
- The inactive tip conducts the instrument safely along the canal
- Blister pack of 6 pieces



Size		21 mm	25 mm	31 mm	
06	●	0363 021 006	0363 025 006	0363 031 006	
08	●	0363 021 008	0363 025 008	0363 031 008	
10	●	0363 021 010	0363 025 010	0363 031 010	
06-10		0363 021 210	0363 025 210	0363 031 210	
12.5	●	0363 021 012	0363 025 012	-	
15	○	0363 021 015	0363 025 015	0363 031 015	
17.5	●	0363 021 017	0363 025 017	-	
20	●	0363 021 020	0363 025 020	0363 031 020	
25	●	0363 021 025	0363 025 025	0363 031 025	
30	●	0363 021 030	0363 025 030	0363 031 030	
35	●	0363 021 035	0363 025 035	0363 031 035	
40	●	0363 021 040	0363 025 040	0363 031 040	
15-40		0363 021 230	0363 025 230	0363 031 230	
45	○	0363 021 045	0363 025 045	0363 031 045	
50	●	0363 021 050	0363 025 050	0363 031 050	
55	●	0363 021 055	0363 025 055	0363 031 055	
60	●	0363 021 060	0363 025 060	0363 031 060	
70	●	0363 021 070	0363 025 070	0363 031 070	
80	●	0363 021 080	0363 025 080	0363 031 080	
45-80		0363 021 240	0363 025 240	0363 031 240	
90	○	0363 021 090	0363 025 090	0363 031 090	
100	●	0363 021 100	0363 025 100	0363 031 100	
110	●	0363 021 110	0363 025 110	0363 031 110	
90-110		0363 021 250	-	-	
120	●	-	0363 025 120	0363 031 120	
130	●	-	0363 025 130	0363 031 130	
140	●	-	0363 025 140	0363 031 140	
90-140		-	0363 025 260	0363 031 260	

K-Files



Single sizes	21 mm	25 mm	28 mm	31 mm
06-140	-	0063 025 XXX	-	0063 031 XXX
06-110	0063 021 XXX	-	0063 028 XXX	-
<b>Assortments</b>				
06-10	0063 021 210	0063 025 210	0063 028 210	0063 031 210
15-40	0063 021 230	0063 025 230	0063 028 230	0063 031 230
45-80	0063 021 240	0063 025 240	0063 028 240	0063 031 240
90-110	0063 021 250	0063 025 250	0063 028 250	0063 031 250
90-140	-	0063 025 260	-	0063 031 260
120-140	-	0063 025 270	-	0063 031 270

VDW®STERILE Hedstroem Files

Highly efficient enlargement of root canals in filing motion

- Excellent cutting ability
- Special steel alloy with uniform structure for maximum resistance to fracture and maximum flexibility
- The inactive tip conducts the instrument safely along the canal
- Blister pack of 6 pieces



Size		21 mm	25 mm	31 mm
08	●	0373 021 008	0373 025 008	0373 031 008
10	●	0373 021 010	0373 025 010	0373 031 010
08-10		0373 021 220	0373 025 220	0373 031 220
15	○	0373 021 015	0373 025 015	0373 031 015
20	●	0373 021 020	0373 025 020	0373 031 020
25	●	0373 021 025	0373 025 025	0373 031 025
30	●	0373 021 030	0373 025 030	0373 031 030
35	●	0373 021 035	0373 025 035	0373 031 035
40	●	0373 021 040	0373 025 040	0373 031 040
15-40		0373 021 230	0373 025 230	0373 031 230
45	○	0373 021 045	0373 025 045	0373 031 045
50	●	0373 021 050	0373 025 050	0373 031 050
55	●	0373 021 055	0373 025 055	0373 031 055
60	●	0373 021 060	0373 025 060	0373 031 060
70	●	0373 021 070	0373 025 070	0373 031 070
80	●	0373 021 080	0373 025 080	0373 031 080
45-80		0373 021 240	0373 025 240	0373 031 240
90	○	0373 021 090	0373 025 090	0373 031 090
100	●	0373 021 100	0373 025 100	0373 031 100
110	●	0373 021 110	0373 025 110	0373 031 110
90-110		0373 021 250	-	-
120	●	-	0373 025 120	0373 031 120
130	●	-	0373 025 130	0373 031 130
140	●	-	0373 025 140	0373 031 140
90-140		-	0373 025 260	0373 031 260





### Hedstroem Files



Single sizes	21 mm	25 mm	28 mm	31 mm
08-140	-	0073 025 XXX	-	0073 031 XXX
08-110	0073 021 XXX	-	0073 028 XXX	-
<b>Assortments</b>				
08-10	0073 021 220	0073 025 220	0073 028 220	0073 031 220
15-40	0073 021 230	0073 025 230	0073 028 230	0073 031 230
45-80	0073 021 240	0073 025 240	0073 028 240	0073 031 240
90-110	0073 021 250	0073 025 250	0073 028 250	0073 031 250
90-140	-	0073 025 260	-	0073 031 260
120-140	-	0073 025 270	-	0073 031 270

### VDW®STERILE Flexicut® Files

For preparation of severely curved and narrow root canals with a filing motion

- Stainless steel
- Special steel alloy with uniform structure for maximum resistance to fracture and maximum flexibility
- Triangular cross-section for highest cutting efficiency
- Cutting angle effective with filing motion
- More cutting edges than K-Reamers
- The inactive tip conducts the instrument safely along the canal
- Blister pack of 6 pieces



Size		21 mm	25 mm	
15	○	0364 021 015	0364 025 015	
20	●	0364 021 020	0364 025 020	
25	●	0364 021 025	0364 025 025	
30	●	0364 021 030	0364 025 030	
35	●	0364 021 035	0364 025 035	
40	●	0364 021 040	0364 025 040	
15-40		0364 021 230	0364 025 230	

### Flexicut® Files



Single sizes	21 mm	25 mm
15-40	0064 021 XXX	0064 025 XXX
<b>Assortments</b>		
15-40	0064 021 230	0064 025 230

## Sterile NiTi Hand Instruments


### VDW®STERILE NiTi K-Files

For the enlargement of severely curved canals

- Super flexible files, follow any canal curvature
- Blister pack of 6 pieces



Size		21 mm	25 mm
15	○	0362 021 015	0362 025 015
20	●	0362 021 020	0362 025 020
25	●	0362 021 025	0362 025 025
30	●	0362 021 030	0362 025 030
35	●	0362 021 035	0362 025 035
40	●	0362 021 040	0362 025 040
15-40		0362 021 230	0362 025 230
45	○	0362 021 045	0362 025 045
50	●	0362 021 050	0362 025 050
55	●	0362 021 055	0362 025 055
60	●	0362 021 060	0362 025 060



### NiTi K-Files



Single sizes	21 mm	25 mm
15-60	0263 021 XXX	0263 025 XXX
Assortments		
15-40	0263 021 230	0263 025 230


## Non-sterile Stainless Steel Instruments

### MC K-Files

- Localisation and exposure of root canal orifices
- Ergonomic handle
- Removing gutta-percha, sealer and cement residue out of isthmi
- Set of 6 instruments
- Optimal working angle in molar region



Size		21 mm
08	●	0060 004 008
10	●	0060 004 010
15	○	0060 004 015
08-15		0060 004 221




### MC Hedstroem Files

- Removing/scraping off gutta-percha residue from canal walls
- Ergonomic handle
- Optimal working angle in molar region
- Set of 6 instruments



Size		21 mm
08	●	0070 004 008
10	●	0070 004 010
15	○	0070 004 015
08-15		0070 004 221



## K-Reamers

### For Contra-Angle

- Max. 800 rpm
- Blister pack of 6 instruments



Size		21 mm	25 mm
08	●	0058 021 008	0058 025 008
10	●	0058 021 010	0058 025 010
08-10		0058 021 220	0058 025 220
15	○	0058 021 015	0058 025 015
20	●	0058 021 020	0058 025 020
25	●	0058 021 025	0058 025 025
30	●	0058 021 030	0058 025 030
35	●	0058 021 035	0058 025 035
40	●	0058 021 040	0058 025 040
15-40		0058 021 230	0058 025 230
45	○	0058 021 045	0058 025 045
50	●	0058 021 050	0058 025 050
55	●	0058 021 055	0058 025 055
60	●	0058 021 060	0058 025 060
70	●	0058 021 070	0058 025 070
80	●	0058 021 080	0058 025 080
45-80		0058 021 240	0058 025 240



## Hedstroem Files

### For Contra-Angle

- Quarter turn only
- 450-800 rpm
- Blister pack of 6 instruments



Size		21 mm	25 mm
10	●	0078 021 010	0078 025 010
15	○	0078 021 015	0078 025 015
20	●	0078 021 020	0078 025 020
25	●	0078 021 025	0078 025 025
30	●	0078 021 030	0078 025 030
35	●	0078 021 035	0078 025 035
40	●	0078 021 040	0078 025 040
15-40		0078 021 230	0078 025 230
45	○	-	0078 025 045
50	●	-	0078 025 050
55	●	-	0078 025 055
60	●	-	0078 025 060
70	●	-	0078 025 070
80	●	-	0078 025 080
45-80		-	0078 025 240



## Silicone Stoppers

### For safe control of instrument working length

- Single-use recommended
- X-ray visible
- Dispenser of 200 pieces

for length		
21 mm	○	0165
25 mm	●	0166
28 mm	●	0167
31 mm	●	0168



## RAYPEX®6

# RAYPEX® 6

High-tech apex locator with “Touch-and-Zoom” for precise length determination

	REF
<b>RAYPEX®6</b> <ul style="list-style-type: none"> <li>• Unique color touch screen</li> <li>• “Plug-and-Start” and “Plug-and-Check” features</li> </ul>	1113 000 000



## Accessories

RAYPEX®6	REF
■ Measuring Cable	1119 000 506
■ Charger	1119 000 501
■ Battery	1119 000 507

## Raypex®5

■ Measuring Cable	0141 000 506
■ Charger	0141 000 501
■ Battery	0141 000 507

## Raypex®4

■ Measuring Cable	0141 000 516
■ Charger	0141 000 517
■ Battery	0141 000 513

## Accessories for all RAYPEX® apex locators Raypex®4, Raypex®5 and RAYPEX®6

■ Lip Clip*, 5 pieces	0141 000 514
■ File Clip*, 2 pieces	0141 000 515
■ Touch Probe*, 2 pieces	0141 000 518

\*autoclavable

## VDW.ULTRA®

# VDW.ULTRA®

### High performance endo ultrasonic device

- Endo focus with clearly marked endo power ranges in large scale
- Correct power setting within the power range allows optimal balance between efficient and gentle operation
- Patented piezo technology



**VDW.ULTRA®** \_\_\_\_\_ REF 1400 000 501

- Incl. ultrasonic handpiece
- Incl. 1 ultrasonic endo tip

**Ultrasonic Handpiece** \_\_\_\_\_ REF 1401 000 000

## VDW Endo Ultrasonic tips and files

### Practical series of endo ultrasonic tips and files



#### IRRI files

##### IRRI S (smooth wire)

- Passive ultrasonic irrigation (Pull)
- Each in blister pack of 4 files

##### IRRI K (K-type)

- Exposure of delicate anatomical structures
- Each in blister pack of 4 files

	Length/ISO			
IRRI S	21/25	1460 021 025		IRRI S
IRRI S	25/25	1460 025 025		
IRRI K	21/15	1463 021 015		IRRI K
IRRI K	21/25	1463 021 025		
IRRI K	25/15	1463 025 015		
IRRI K	25/25	1463 025 025		


#### CAVI tips (diamond-coated)

##### CAVI 1-BD (ball)

- Access refinement
- Each in blister pack of 1 tip

##### CAVI 2-D (swan neck, optimised visibility)

##### CAVI 3-D (standard)

CAVI 1-BD	1410 000 001		CAVI 1
CAVI 2-D	1410 000 002		CAVI 2
CAVI 3-D	1410 000 003		CAVI 3

#### REDO tips

##### REDO 1 (short)

- Retirements
- Each in blister pack of 1 tip

##### REDO 2 (medium, special Ti-Nb alloy)

##### REDO 5 (long)

REDO 1	1411 000 001		REDO 1
REDO 2	1411 000 002		REDO 2
REDO 5	1411 000 005		REDO 5

#### MAXI tip

##### MAXI MPR

- Dislodging metal posts at high power level
- In blister pack of 1 tip

MAXI MPR	1412 000 000		MAXI
----------	--------------	--	------

## EDDY™

**The innovative endo irrigation tip – highly efficient, safe and simple to use**

- Sonic powered polyamide tip
- Powerful activation at max. 6000 Hz via Airscaler
- Forging technique (no ledges/steps)
- Highly efficient, even in curved canals and complex canal anatomies
- Oscillating, 3D movement enables acoustic streaming and cavitation

### ■ EDDY™

- 5 blister packs of 2 Irrigation Tips, sterile

REF

1441 000 000

### ■ EDDY™ Starter Kit

- 1 W&H Airscaler Handpiece, Type Proxeo ZA-55-LM with light and power setting
- 20 blister packs of 2 EDDY™ Irrigation Tips, sterile

1443 000 000



## GUTTAFUSION® for RECIPROC®

# GUTTAFUSION®

Suitable for the RECIPROC® System. The use of Size Verifiers is not necessary.

### GUTTAFUSION® for RECIPROC® Obturators

Obturators made entirely of gutta-percha

- Good radiopacity
- Colour coded according to the RECIPROC® system
- Blister pack of 6 pieces

Blister pack of 6 obturators		1 blister pack	5 blister packs	
R25	●	1531 000 025	1532 000 025	
R40	●	1531 000 040	1532 000 040	
R50	●	1531 000 050	1532 000 050	

## GUTTAFUSION®


Universal system for 3-dimensional obturation after preparation with rotary NiTi systems, for example Mtwo®.

GUTTAFUSION® Size Verifiers are necessary with these NiTi systems in order to ensure correct selection of the appropriate obturator.

### GUTTAFUSION® Obturators

Obturators made entirely of gutta-percha

- Good radiopacity
- ISO colour coded
- Blister pack of 6 pieces


Blister pack of 6 obturators		1 blister pack	5 blister packs	
20	●	1541 000 020	1542 000 020	
25	●	1541 000 025	1542 000 025	
30	●	1541 000 030	1542 000 030	
35	●	1541 000 035	1542 000 035	
40	●	1541 000 040	1542 000 040	
45	○	1541 000 045	-	
50	●	1541 000 050	-	
55	●	1541 000 055	-	

### GUTTAFUSION® Size Verifier

Instruments to determine the correct obturator size after preparation with rotary NiTi systems

- Nickel-Titanium
- Radiopaque
- ISO colour coded
- Blister pack of 6 pieces



Blister pack of 6 size verifiers		1 blister pack	
20	●	1508 000 020	
25	●	1508 000 025	
30	●	1508 000 030	
35	●	1508 000 035	
40	●	1508 000 040	
45	○	1508 000 045	
50	●	1508 000 050	
55	●	1508 000 055	
25-35		1508 000 305	

■ **GUTTAFUSION® for RECIPROC® Basic Kit** \_\_\_\_\_ REF 1530 000 000

- 3 blister packs of 6 GUTTAFUSION® for RECIPROC® Obturators R25
- 1 blister pack of 6 GUTTAFUSION® for RECIPROC® Obturators R40
- 1 blister pack of 6 GUTTAFUSION® for RECIPROC® Obturators R50
- Obturation Practice Block suitable for GUTTAFUSION® for RECIPROC® Obturators R25
- 1 x RECIPROC® Paper Points assorted, sizes R25, R40, R50, sterile
- User card, directions for use



■ **GUTTAFUSION® for RECIPROC® Starter Kit** \_\_\_\_\_ 1535 000 000  
230 V

- GUTTAFUSION® Oven
- GUTTAFUSION® for RECIPROC® Basic Kit



■ **GUTTAFUSION® Basic Kit for Mtwo®** \_\_\_\_\_ 1540 000 000

- 2 blister packs of 6 GUTTAFUSION® Obturators size 25
- 1 blister pack of 6 GUTTAFUSION® Obturators size 30
- 1 blister pack of 6 GUTTAFUSION® Obturators size 35
- 1 blister pack of 6 GUTTAFUSION® Size Verifiers assorted sizes 25, 30, 35
- Obturation Practice Block suitable for GUTTAFUSION® Obturators size 25
- 1 x Mtwo® Paper Points assorted, sizes 25/.06, 30/.05, 35/.04, 40/.04, 25/.07, sterile
- User card, directions for use



■ **GUTTAFUSION® Starter Kit for Mtwo®** \_\_\_\_\_ 1545 000 000  
230 V

- GUTTAFUSION® Oven
- GUTTAFUSION® Basic Kit for Mtwo®



■ **GUTTAFUSION® Oven** \_\_\_\_\_ 1500 000 000  
230 V

- Oven for simultaneous heating of up to two obturators
- Visual and acoustic warning signals
  - Simple cleaning thanks to special mode



■ **GuttaMaster® PREPI® Burs** \_\_\_\_\_ 0515 000 000

- Burs for quick and easy cutting of the handle of GuttaMaster® obturators
- Available in different diameters: 1 x 0.7 mm / 0.9 mm each, 2 x 1.2 mm / 1.5 mm each; assortment 6 pack
  - Use for high speed turbine contra-angle 150,000 up to 200,000 rpm





## BeeFill®2in1

### Obturation device for Downpack and Backfill technique

- Excellent visibility of the root canal entrance thanks to slim handpieces and extra long canulas
- Light-weight handpieces (20.1 g and 127.4 g)
- Steady application thanks to precise and micromotor controlled release of gutta-percha
- Safe use due to visual and audible signals

# BeeFill®2in1

REF



### ■ BeeFill®2in1

0541 000 000

### ■ BeeFill®2in1 Obturation Kit

0545 000 000

- 1 x BeeFill®2in1 device
- 3 x BeeFill® heat pluggers: standard 40/.03, soft 50/.05, soft 60/.06
- 1 x BeeFill® thermal response tip
- 30 BeeFill® cartridges, 10 pieces each size:  
20G = Ø 0.80 mm | 23G = Ø 0.60 mm | 25G = Ø 0.45 mm
- 3 x Machtou hand pluggers, size: 0 (NiTi), 1-2 and 3-4
- 1 x 2Seal easymiX® root canal sealer
- 2 x alpha Gutta-Percha, assorted ISO 15-40, taper .04 and .06



### BeeFill® Heat Pluggers

- For Downpack and searing off gutta-percha
- Continuous Wave or Schilder technique
- Fast heating and cooling down
- Autoclavable

Heat plugger standard 40/.03	0547 000 040
Heat plugger soft 50/.05	0546 000 050
Heat plugger soft 60/.06	0546 000 060



REF

### ■ BeeFill® Thermal Response Tip

Pulp sensitivity tester

0548 000 000



### Machtou Hand Pluggers

**For warm, vertical compaction**

- Ergonomic silicone handle
- Autoclavable
- Stainless steel pluggers with length markings
- Flexible NiTi hand plugger (size 0) for narrow and curved canals

Size/ISO	
0/40 (NiTi)	1063 000 040
1-2/50+60	1063 001 002
3-4/80+100	1063 003 004



### BeeFill® Cartridges

**Injection of warm gutta-percha**

- Flexible and extra long canula – also for curved root canals
- Matching cartridge for each canal anatomy: 25G = Ø 0.45 mm | 23G = Ø 0.60 mm | 20G = Ø 0.80 mm
- One cartridge fills 4-6 root canals
- 10 cartridges per blister pack

Size	
20G / 0.80 mm	0532 053 080
23G / 0.60 mm	0532 053 060
25G / 0.45 mm	0532 053 045



### Paper Points Taper .02

- Sizes 15-80 slide cassette of 200 pieces
- Sizes 90-140 slide cassette of 90 pieces



Size/Taper		29 mm
15/.02	○	0550 029 015
20/.02	●	0550 029 020
25/.02	●	0550 029 025
30/.02	●	0550 029 030
35/.02	●	0550 029 035
40/.02	●	0550 029 040
15-40/.02		0550 029 230
45/.02	○	0550 029 045
50/.02	●	0550 029 050
55/.02	●	0550 029 055
60/.02	●	0550 029 060
70/.02	●	0550 029 070
80/.02	●	0550 029 080
45-80/.02		0550 029 240
90-140/.02		0550 029 260



### Paper Points Taper .04

- Sizes 15-40, slide cassette of 120 pieces



Size/Taper		29 mm
15/.04	○	0558 029 015
20/.04	●	0558 029 020
25/.04	●	0558 029 025
30/.04	●	0558 029 030
35/.04	●	0558 029 035
40/.04	●	0558 029 040
15-40/.04		0558 029 230



### Paper Points Taper .06

- Sizes 15-40, slide cassette of 120 pieces



Size/Taper		29 mm
15/.06	○	0559 029 015
20/.06	●	0559 029 020
25/.06	●	0559 029 025
30/.06	●	0559 029 030
35/.06	●	0559 029 035
40/.06	●	0559 029 040
15-40/.06		0559 029 230



### Paper Points sterile Taper .02

- Sterile packed in blister cells to support infection control
- Sizes 15-80, 180 pieces: 5 paper points per blister cell
- Sizes 90-140, 144 pieces: 4 paper points per blister cell
- Assortments: 6 blister cells of each size



Size/Taper		29 mm
15/.02	○	0560 029 015
20/.02	●	0560 029 020
25/.02	●	0560 029 025
30/.02	●	0560 029 030
35/.02	●	0560 029 035
40/.02	●	0560 029 040
15-40/.02		0560 029 230
45/.02	○	0560 029 045
50/.02	●	0560 029 050
55/.02	●	0560 029 055
60/.02	●	0560 029 060
70/.02	●	0560 029 070
80/.02	●	0560 029 080
45-80/.02		0560 029 240
90-140/.02		0560 029 260



### Mtwo® Paper Points

Suitable for the Mtwo® NiTi system

- Fully saturated points retain good material consistency
- Sterilised points in blister cells to support infection control
- 144 pieces: 4 paper points per blister cell



Size/Taper		29 mm
25/.06	●	0225 629 025
30/.05	●	0225 529 030
35/.04	●	0225 429 035
40/.04	●	0225 429 040
45/.04	○	0225 429 045
50/.04	●	0225 429 050
60/.04	●	0225 429 060
25/.07	●	0225 729 025
30/.06	●	0225 629 030
35/.06	●	0225 629 035
40/.06	●	0225 629 040
48 x 25/.06; 24 x 30/.05, 35/.04, 40/.04, 25/.07 each		0225 029 457



## RECIPROC® Paper Points

### Suitable for the RECIPROC® NiTi system

- Fully saturated points retain good material consistency
- Sterilised points in blister cells to support infection control
- 144 pieces: 4 paper points per blister cell

STERILE

Size/Taper		29 mm
R25	●	0216 029 025
R40	●	0216 029 040
R50	●	0216 029 050
96 x R25, 24 x R40, 24 x R50		0216 029 237



REF

## 2Seal®

1020 000 000

Root canal sealer on epoxy-amine resin basis

- Excellent biocompatibility
- Very low shrinkage
- Good sealing properties
- For all obturation methods: thermostable up to 250°C
- For manual mixing: 2-tube pack with 4 ml paste A and B



## 2Seal easymiX®

# 2Seal easymiX

Root canal sealer on epoxy-amine resin basis in convenient double-chamber mixing syringe. Easy application with interchangeable, flexible intra oral tips.

- Always mixed in ideal ratio
- Can be applied directly into the root canal
- Precise application even in canals with difficult access



REF

## 2Seal easymiX® Starter Kit

1025 000 000

- 1 x 2Seal easymiX® Mixing Syringe, 15 g
- 20 x Mixing Tips with Intra Oral Tips
- 1 x Mixing Block
- 1 x Illustrated Manual
- 1 x Organiser



## 2Seal easymiX® Mixing Syringe

1025 000 001

- 2 x 2Seal easymiX® Mixing Syringe, 15 g



## 2Seal easymiX® Mixing Tips

1025 000 002

- 40 x Mixing Tips with Intra Oral Tips



## Standardised Gutta-Percha Cones

- Best plasticity for lateral compaction
- High radiopacity

### Gutta-Percha Taper .02

- Sizes 15-80 slide cassette of 120 pieces
- Sizes 90-140 slide cassette of 60 pieces

Size/Taper		28 mm
15/.02	○	0525 028 015
20/.02	●	0525 028 020
25/.02	●	0525 028 025
30/.02	●	0525 028 030
35/.02	●	0525 028 035
40/.02	●	0525 028 040
15-40/.02		0525 028 230
45/.02	○	0525 028 045
50/.02	●	0525 028 050
55/.02	●	0525 028 055
60/.02	●	0525 028 060
70/.02	●	0525 028 070
80/.02	●	0525 028 080
45-80/.02		0525 028 240
90-140/.02		0525 028 260



### Gutta-Percha Taper .04 and .06

For fast obturation of tapered canals

- Slide cassette of 60 pieces

Size/Taper		28 mm
15/.04	○	0528 028 015
20/.04	●	0528 028 020
25/.04	●	0528 028 025
30/.04	●	0528 028 030
35/.04	●	0528 028 035
40/.04	●	0528 028 040
15-40/.04		0528 028 230



Size/Taper		28 mm
15/.06	○	0529 028 015
20/.06	●	0529 028 020
25/.06	●	0529 028 025
30/.06	●	0529 028 030
35/.06	●	0529 028 035
40/.06	●	0529 028 040
15-40/.06		0529 028 230



### alpha Gutta-Percha Taper .04 and .06

#### For thermoplastic obturation

- Ideal for use with BeeFill®2in1, BeeFill®Pack etc.
- With low melting point
- Slide cassette of 60 pieces

Size/Taper		28 mm
15/.04	○	0521 028 015
20/.04	●	0521 028 020
25/.04	●	0521 028 025
30/.04	●	0521 028 030
35/.04	●	0521 028 035
40/.04	●	0521 028 040
15-40/.04		0521 028 230



Size/Taper		28 mm
15/.06	○	0522 028 015
20/.06	●	0522 028 020
25/.06	●	0522 028 025
30/.06	●	0522 028 030
35/.06	●	0522 028 035
40/.06	●	0522 028 040
15-40/.06		0522 028 230



### Mtwo® Gutta-Percha

#### Suitable for the Mtwo® NiTi system

- $\alpha$ -phase with low melting point for cold and warm obturation techniques
- Slide cassette of 60 pieces

Size/Taper		28 mm
25/.06	●	0220 628 025
30/.05	●	0220 528 030
35/.04	●	0220 428 035
40/.04	●	0220 428 040
45/.04	○	0220 428 045
50/.04	●	0220 428 050
60/.04	●	0220 428 060
25/.07	●	0220 728 025
30/.06	●	0220 628 030
35/.06	●	0220 628 035
40/.06	●	0220 628 040
20 x 25/.06, 10 x 30/.05, 35/.04, 40/.04, 25/.07 each		0220 028 457



### RECIPROC® Gutta-Percha

#### Suitable for the RECIPROC® NiTi system

- $\alpha$ -phase with low melting point for cold and warm obturation techniques
- Slide cassette of 60 pieces

Size		28 mm
R25	●	0214 028 025
R40	●	0214 028 040
R50	●	0214 028 050
40 x R25, 10 x R40, 10 x R50		0214 028 237



## Conventional Gutta-Percha

### Additional cones for the masterpoint technique

- Sizes XXF-L, slide cassette of 100 pieces

Size	20 mm
xxfine	0542 020 431
xfine	0542 020 441
fine	0542 020 451
medium	0542 020 461
large	0542 020 471
xf-l	0542 020 477



## Silver Points

- Approx. 99.85 % pure silver
- Max. 0.12-0.13 % nickel content
- Practical slide cassette of 30 pieces

Size	30 mm
10	0098 030 010
15	0098 030 015
20	0098 030 020



REF

### GuttaCut®

0143

Cordless, rechargeable handpiece with loading station and 4 interchangeable cutting tips for fast, precise and safe cutting of excess gutta-percha

- The loop on the cutting tip provides easy access to the cavity and excess gutta-percha can be cut off exactly at the canal orifice
- The handpiece is inserted cold and the cutting tip heats up instantly by pressing the button. When the button is released, the tip cools down within seconds.



### GuttaCut® Handpiece

0511

- A second handpiece for the charger



## GuttaCut® Cutting Tips

- Autoclavable
- Set of 4 pieces

Size	
1	0512 000 001
2	0512 000 002
3	0512 000 003
4	0512 000 004
1-4	0512 000 150



### VDW®STERILE Root Filler “L”

For filling root canals with paste, cement or sealer

- Blister of 4 pieces, sterile



Size	21 mm ●	25 mm ●	21 mm ▲	25 mm ▲
25	● 0393 021 025	● 0393 025 025	▲ 0093 021 025	▲ 0093 025 025
30	● 0393 021 030	● 0393 025 030	▲ 0093 021 030	▲ 0093 025 030
35	● 0393 021 035	● 0393 025 035	▲ 0093 021 035	▲ 0093 025 035
40	● 0393 021 040	● 0393 025 040	▲ 0093 021 040	▲ 0093 025 040
25-40	0393 021 290	0393 025 290	0093 021 290	0093 025 290



### VDW®STERILE NiTi Finger Spreaders

For lateral condensation in curved canals

- Blister of 6 pieces, sterile



Size	25 mm ●	25 mm ▲
15 ○	● 0392 025 015	▲ 0295 025 015
20 ●	● 0392 025 020	▲ 0295 025 020
25 ●	● 0392 025 025	▲ 0295 025 025
30 ●	● 0392 025 030	▲ 0295 025 030
35 ●	● 0392 025 035	▲ 0295 025 035
40 ●	● 0392 025 040	▲ 0295 025 040
15-40	0392 025 230	0295 025 230



### VDW®STERILE Finger Spreaders

For lateral condensation of gutta-percha cones (with pointed tip)

- Stainless steel
- Length of tapered section 16 mm
- Blister of 6 pieces, sterile



Size	25 mm ●	25 mm ▲
15 ○	● 0395 025 015	▲ 0095 025 015
20 ●	● 0395 025 020	▲ 0095 025 020
25 ●	● 0395 025 025	▲ 0095 025 025
30 ●	● 0395 025 030	▲ 0095 025 030
35 ●	● 0395 025 035	▲ 0095 025 035
40 ●	● 0395 025 040	▲ 0095 025 040
15-40	0395 025 230	0095 025 230



### VDW®STERILE Finger Pluggers

For vertical condensation of gutta-percha cones (with blunt tip)

- Stainless steel
- Blister of 6 pieces, sterile



Size	25 mm ●	25 mm ▲
15 ○	● 0399 025 015	▲ 0099 025 015
20 ●	● 0399 025 020	▲ 0099 025 020
25 ●	● 0399 025 025	▲ 0099 025 025
30 ●	● 0399 025 030	▲ 0099 025 030
35 ●	● 0399 025 035	▲ 0099 025 035
40 ●	● 0399 025 040	▲ 0099 025 040
15-40	0399 025 230	0099 025 230





## DT ILLUSION™XRO®SL

- Invisible in the tooth – but clearly visible during treatment
- Easier placement and retreatment thanks to colour pigments, enabling posts to be distinguished when required – even without magnification
- Easy post size identification at all times, without any risk of confusion



**DT ILLUSION™SL**  
XRO®

REF

### DT ILLUSION™XRO®SL Basic Kit \_\_\_\_\_ 1940 000 000

- 12 posts: 3 per size: #0 / #1 / #2 / #3
- 4 DT Finishing Drills: 1 per size: #0 (Universal Drill) / #1 / #2 / #3



### DT ILLUSION™XRO®SL Refill Blister Pack

- Individually packed for convenience
- 6 posts of one size in each blister pack

Size	
DT ILLUSION™XRO®SL #0	1941 000 000
DT ILLUSION™XRO®SL #1	1941 000 001
DT ILLUSION™XRO®SL #2	1941 000 002
DT ILLUSION™XRO®SL #3	1941 000 003




## DT LIGHT®SL

### Unique Safety Lock® Coating

Industrially preconditioned post with silane-silicate coating and protective polymer layer

- Reliable bonding with the composite cement
- Long-term bonding adhesion
- Compatible with all adhesive systems



**SAFETY LOCK®**  
TECHNOLOGY



REF

### DT LIGHT®SL Basic Kit \_\_\_\_\_ 1930 000 000

- 12 posts: 3 per size: #0 / #1 / #2 / #3
- 4 DT Finishing Drills: 1 per size: #0 (Universal Drill) / #1 / #2 / #3



### DT LIGHT®SL Refill Blister Pack

- Individually packed for convenience
- 6 posts of one size in each blister pack

Size	
DT LIGHT®SL #0	1931 000 000
DT LIGHT®SL #1	1931 000 001
DT LIGHT®SL #2	1931 000 002
DT LIGHT®SL #3	1931 000 003



## DT LIGHT®

### Clinically tried and tested quartz fibre posts

- Root-friendly: elasticity modulus similar to dentine
- Aesthetics: high translucency
- Minimally invasive: Double Taper design
- Clinical control: radiopaque
- Durable and resistant



### DT LIGHT® Basic Kit

- 12 posts: 3 per size: #0 / #1 / #2 / #3
- 4 DT Finishing Drills: 1 per size: #0 (Universal Drill) / #1 / #2 / #3

REF

0931 000 500



### DT LIGHT® Refill Blister Pack

- Individually packed for convenience
- 6 posts of one size in each blister pack

Size	
DT LIGHT® #0	0931 000 000
DT LIGHT® #1	0931 000 001
DT LIGHT® #2	0931 000 002
DT LIGHT® #3	0931 000 003



## DT Accessories

### DT Universal Drill (Finishing Drill #0)

Drill for removing the root canal filling

- Pack of 2 pieces

REF

0935



### DT Finishing Drill

Drill for post fitting in the root canal

- Pack of 2 pieces

DT Finishing Drill #1  
0935 001

DT Finishing Drill #2  
0935 002

DT Finishing Drill #3  
0935 003



### Composibrush

Endo brush for applying bonding material to the root canal

- Pack of 50 pieces

0936 000












### DT Post Removal Kit

For the fast, safe removal of DT Posts with no unnecessary dentine removal

- Cassette of 2 drills

0938 001



	REF
<p>■ <b>Endo Training Model Castillo</b> _____ 0247 000 000</p> <p>Training model for root canal preparation and electronic length determination</p> <ul style="list-style-type: none"> <li>• Plastic model of an upper jaw (11-17)</li> <li>• Position 16: Endo Training Tooth</li> <li>• Position 17: Plastic Training Block</li> <li>• Fixture for electronic length determination</li> <li>• Integrated Minifix Measuring Gauge</li> </ul>	
<p>■ <b>Plastic Training Blocks</b> _____ 0245</p> <ul style="list-style-type: none"> <li>• 3 pieces</li> </ul>	
<p>■ <b>Endo Training Tooth</b> _____ 0246 000 500</p> <p>Plastic tooth for endodontic training</p> <ul style="list-style-type: none"> <li>• Replica of a maxillary first molar with accessed pulp chamber</li> <li>• 4 root canals, 3 canals with stained pulp, MB2 pulp unstained</li> </ul>	
<p>■ <b>Minifix Measuring Gauge</b> _____ 0179</p> <p>For precise adjustment of the working length of endo instruments, gutta-percha and paper points</p>	
<p>■ <b>LavEndo®</b> _____ 0479</p> <p>Patented box for convenient wash-disinfection of root canal instruments</p> <ul style="list-style-type: none"> <li>• Without module, without instruments</li> </ul>	
<p>■ <b>LavEndo® with Mini Step Module</b> _____ 0477</p> <ul style="list-style-type: none"> <li>• For stainless steel instruments</li> <li>• Without instruments</li> </ul>	
<p>■ <b>LavEndo® with FlexMaster® Module</b> _____ 0478</p> <ul style="list-style-type: none"> <li>• For FlexMaster® instruments</li> <li>• Without instruments</li> </ul>	
<p>■ <b>Interim Stand</b> _____ 0495</p> <p>For quick chairside storing and cleaning of root canal instruments during treatment</p> <ul style="list-style-type: none"> <li>• Foam disc autoclavable for one single use</li> <li>• Interim Stand sterilisable without foam disc, resistant to temperatures up to 200°C/390°F</li> <li>• Without instruments</li> </ul>	
<p>■ <b>Refills for the Interim Stand</b> _____ 0496</p> <ul style="list-style-type: none"> <li>• Foam discs autoclavable for one single use</li> <li>• 55 foam discs</li> </ul>	

### MiniBox 2100

#### Box with lid

##### With Mini Step Module

- Holds 24 instruments
- Replaceable plates for marking each instrument row

or

##### With Mini Socket Module

- With 8 colour coded sockets for obturation points

MiniBox 2100, perforated	
Box with Mini Step Module	0471
Mini Step Module	0403
Box with Mini Socket Module	0472
Mini Socket Module	0404



### BasicBox 2100

#### Box with lid

##### With Step Module

- Holds 48 instruments
- Replaceable plates for marking each instrument row

or

##### With Socket Module

- With 18 colour coded sockets for obturation points

BasicBox 2100, perforated	
Box with Step Module	0441
Step Module	0401
Box with Socket Module	0442
Socket Module	0402



### SemiBox 2000

#### Box with lid

##### With 2 Step Modules

- Holds 48 instruments each
- Replaceable plates for marking each instrument row

or

##### With Step Module and Socket Module

- Step Module holds 48 instruments and has replaceable plates for marking each instrument row
- Socket module has 18 colour coded sockets for obturation points

SemiBox 2000, perforated	
Box with 2 Step Modules	0421
Box with Step Module and Socket Module	0422
Step Module	0401
Socket Module	0402



### Autoclave Paper

- Box of 250 sheets

For MiniBox 2100, FlexMaster® SystemBox, FlexMaster® AccessoryBox and Mtwo® SystemBox	0191
For BasicBox 2100 and FlexMaster® CombiBox	0192
For SemiBox 2000	0193



## Your International Contacts



**Leo Satanovskij**  
Head of  
International Sales 1  
Tel.: +49 89 - 62 73 41 64  
Central and Eastern  
Europe, CIS, Russia,  
Ukraine, Croatia,  
Rumania

leo.satanovskij@vdw-dental.com



**Alexandra Magg**  
Assistant  
International Sales 1  
Tel.: +49 89 - 62 73 41 62  
Central and Eastern  
Europe, CIS, Scandinavia,  
Baltic States, The Balkans,  
Middle East, Africa, UK,  
Ireland

alexandra.magg@vdw-dental.com



**Wolfgang Münch**  
Export Manager  
International Sales 1  
Tel.: +49 89 - 62 73 41 65  
UK, Ireland, Greece,  
Turkey, Middle East,  
Africa

wolfgang.muench@vdw-dental.com



**Hartmut Weiss**  
Export Manager  
International Sales 1  
Tel.: +49 89 - 62 73 41 35  
Scandinavia, Baltic  
States, The Balkans,  
Hungary, Slovakia,  
Italy, France

hartmut.weiss@vdw-dental.com



**Rolf zum Stickling**  
Head of  
International Sales 2  
Tel.: +49 89 - 62 73 41 63  
Australia, New Zealand,  
USA

rolf.zumstickling@vdw-dental.com



**Lydia Bonenberger**  
Assistant  
International Sales 2  
Tel.: +49 89 - 62 73 41 88  
Asia, Australia, Iberia,  
Latin America

lydia.bonenberger@vdw-dental.com



**Di Hu**  
Export Manager  
International Sales 2  
Tel.: +49 89 - 62 73 41 69  
Asia

di.hu@vdw-dental.com



**Rosalba Gutierrez**  
Territory Manager  
International Sales 2  
Tel.: +593 98 - 422 7247  
Latin America

rosalba.gutierrez@vdw-dental.com



**David Rios Zurdo**  
Regional Sales Manager  
International Sales 2  
Tel.: +34 634 296 686  
Iberia

david.rios@vdw-dental.com



**Sergio González Mije**  
Promoter  
International Sales 2  
Tel.: +34 607 673 979  
Iberia

sergio.gonzalez@vdw-dental.com



**Tino Joris**  
Territory Manager  
Sales Benelux  
Tel.: +31 165 - 30 86 87  
Netherlands

tino.joris@vdw-dental.com




**Sabien Baisier**  
Territory Manager  
Sales Benelux  
Tel.: +32 484 - 15 35 18  
Belgium & Luxembourg

sabien.baisier@vdw-dental.com

# Product Index

2Seal®	28	GUTTAFUSION® Oven	24
2Seal easymiX® Mixing Syringe	28	GUTTAFUSION® Size Verifiers	23
2Seal easymiX® Mixing Tips	28	GUTTAFUSION® Starter Kit for Mtwo®	24
2Seal easymiX® Starter Kit	28	GuttaMaster® PREPI® Burs	24
6:1 Contra-Angle	13	Hedstroem Files	16
alpha Gutta-Percha .04 und .06	30	Hedstroem Files (for contra-angle)	19
Autoclave Paper	36	Interim Stand	35
Barbed Broaches	5	K-Files	15
BasicBox 2100	36	K-Reamers	14
BeeFill®2in1	25	K-Reamers (for contra-angle)	19
BeeFill®2in1 Obturation Kit	25	LavEndo®	35
BeeFill® Cartridges	26	Machtou-Hand Pluggers	26
BeeFill® Heat Pluggers	25	MC Hedstroem Files	18
BeeFill® Thermal Response Tip	25	MC K-Files	18
Beutelrock Enlargers	5	MiniBox 2100	36
Beutelrock Reamers	4	Minifix Measuring Gauge	35
C-PILOT® Files	14	Mtwo® Basic Kit	8
Composibrush	34	Mtwo® Control Stickers	8
Conventional Gutta-Percha	31	Mtwo® Gutta-Percha	30
DT Finishing Drill	34	Mtwo® Instruments	7
DT ILLUSION™XRO®SL Basic Kit	33	Mtwo® Modul	8
DT ILLUSION™XRO®SL Refill Blister Pack	33	Mtwo® Paper Points	27
DT LIGHT® Basic Kit	34	Mtwo® Retreatment Files	8
DT LIGHT® Refill Blister Pack	34	Mtwo® SystemBox	8
DT LIGHT®SL Basic Kit	33	NiTi Finger Spreader	32
DT LIGHT®SL Refill Blister Pack	33	NiTi K-Files	18
DT Post Removal Kit	34	Paper Points	26
DT Universal Drill (Finishing Drill #0)	34	Peeso Enlargers	4
EDDY™	22	Plastic Training Blocks	35
EDDY™ Starter Kit	22	Raypex®4 Accessories	20
Endo Training Model Castillo	35	Raypex®5 Accessories	20
Endo Training Tooth	35	RAYPEX®6	20
FileCare®EDTA	10	RAYPEX®6 Accessories	20
Finger Pluggers	32	RECIPROC® Gutta-Percha	30
Finger Spreaders	32	RECIPROC® Instruments	6
Flexicut® Files	17	RECIPROC® Paper Points	28
FlexMaster® AccessoryBox	10	RECIPROC® System Kit 24	6
FlexMaster® Basic Kit	10	Root Filler “L”	32
FlexMaster® CombiBox	10	SemiBox 2000	36
FlexMaster® Control Stickers	10	Silicone Stoppers	19
FlexMaster® Instruments	9	Silver Points	31
FlexMaster® IntroFile	9	Ultrasonic handpiece	21
FlexMaster® SystemBox	10	VDW.CONNECT Accessories for iPad mini	11
Gates Enlargers	4	VDW.CONNECT App	12
Gutta-Percha Cones	29	VDW.CONNECT Drive®	11
GuttaCut®	31	VDW.CONNECT Drive® Accessories	11
GuttaCut® Cutting Tips	31	VDW.CONNECT Drive® Set	11
GuttaCut® Handpiece	31	VDW.CONNECT Drive® Set Medium	11
GUTTAFUSION® Basic Kit for Mtwo®	24	VDW.GOLD®RECIPROC®	13
GUTTAFUSION® for RECIPROC® Basic Kit	24	VDW.GOLD®RECIPROC® Accessories	13
GUTTAFUSION® for RECIPROC® Obturators	23	VDW.SILVER®RECIPROC®	13
GUTTAFUSION® for RECIPROC® Starter Kit	24	VDW.ULTRA®	21
GUTTAFUSION® Obturators	23	VDW Endo Ultrasonic tips and files	21





**VDW GmbH**

Bayerwaldstr. 15 • 81737 Munich • Germany  
Tel. +49 89 62734-0 • Fax +49 89 62734-304  
[www.vdw-dental.com](http://www.vdw-dental.com) • [info@vdw-dental.com](mailto:info@vdw-dental.com)