

PROSTHETICS



DENITAL CATALOG

WHAT DO YOU CALL A CONVICTION
THAT EVERYTHING CAN ALWAYS BE
A LITTLE BETTER, EASIER AND MORE
RELIABLE?

WE CALL IT
Upgrade Dentistry

With expertise and enthusiasm, COLTENE turned this claim into reality. We know exactly our end customers needs and want and work tirelessly to provide them with innovative and smart solutions.

This is not an empty promise, it is our mission. Upgrading Dentistry makes dental treatments more effective, safe and economical for dental professionals and patients.

Better Quality. Better Reliability. Better Practise.



www.coltene.com

TABLE OF CONTENTS

1 RETRACTION

ROEKO Comprecap / Comprecap anatomic	10
ROEKO Stay-put	10
Retraction Accessories	11
ROEKO Retracto	11

2 BITE REGISTRATION

JET BLUE BITE	14
JET BITE	14

3 IMPRESSION A-SILICONE

AFFINIS®	18
AFFINIS light body	19
AFFINIS fast light body	19
AFFINIS regular body	19
AFFINIS fast regular body	20
AFFINIS PRECIOUS	20
AFFINIS MonoBody	21
AFFINIS heavy body	21
AFFINIS fast heavy body	22
AFFINIS BLACK EDITION	22
AFFINIS System 360 putty	22
AFFINIS putty soft / fast	23
AFFINIS putty super soft	23
PRESIDENT The Original	24
PRESIDENT The Original Xtra light body	25
PRESIDENT The Original light body	25
PRESIDENT The Original regular body	25
PRESIDENT The Original MonoBody	25
PRESIDENT The Original heavy body	26
PRESIDENT The Original putties	26

4 IMPRESSION C-SILICONE

Speedex	30
Speedex light body	31
Speedex medium	31
Speedex putty	31
Speedex putty soft	32
Speedex Activator	32
Pressure Spot Indicator	32

5 ALGINATE

HydroPrint Premium	36
--------------------	----

6 PROSTHETICS ACCESSORIES

COLTENE Adhesive	40
Impression Trays	41
Impression Trays Sets	41
Impression Trays Quadrant & Anterior	42
Mixing Devices	42
Mixing Tips	43
Other Accessories	43

7 LABORATORY

GI-MASK	46
GI-Mask 3D	47
GI-MASK Automix New Formula	47
GI-MASK	48
Lab putty	49
Lab putty / Lab putty hard	50
ACCU-TRAC	51
ACCU-TRAC®	52
PINDEX	53
PINDEX®	54
HYGENIC Waxes	55
HYGENIC yellow bite wax	55
HYGENIC bite-wafers	56

TABLE OF CONTENTS

HYGENIC boxing wax	56
HYGENIC medium-soft no. 3 pink/red wax	56
HYGENIC u-shaped occlusal rim wax	56
HYGENIC orthodontic tray wax	57
HYGENIC wax ropes	57
HYGENIC wax sticks	57
HYGENIC extra tough wax	57
Laboratory Accessories	58
Laboratory Accessories	59

SPOTLIGHT



ROEKO Comprecap

10

ROEKO Comprecap compression caps stop bleeding naturally by compression and control the moisture around the tooth stump.

ROEKO Comprecap / Comprecap anatomic



530000



530000

ROEKO Comprecap compression caps stop bleeding naturally by compression and control the moisture. The use is simple and supports impression preparation regardless of the retraction technique. ROEKO Comprecap anatomic are formed compression caps with semicircled shaped spaces on two opposite sides and can therefore be easily placed on adjacent teeth.

REF	Item
530000	Set, 150 pcs
530008	Size 1, 120 pcs
530010	Size 3, 120 pcs
530014	Size 5, 60 pcs
531008	Size 1, 600 pcs
531010	Size 3, 600 pcs
531014	Size 5, 400 pcs
531201	Anatomic Size 1, 120 pcs
531203	Anatomic Size 3, 120 pcs
531205	Anatomic Size 5, 60 pcs
531299	Anatomic Set, 120 pcs

ROEKO Stay-put



521000

Stay-put combines the advantages of a braided retraction cord with the adaptability of a fine copper filament. The pliable core is so effective that the cord is not only easy to place in the sulcus but it stays there.

REF	Item
521000	Non-Impregnated, X-Fine, 1 pcs
521001	Non-Impregnated, Fine, 1 pcs
521002	Non-Impregnated, Medium, 1 pcs
521003	Non-Impregnated, Thick, 1 pcs

Retraction Accessories

Double ended, stainless steel instruments facilitate placing of retraction cord around the tooth.



520000

- Serrated or smooth edges
- Structured handle
- Ergonomically designed to work around the tooth
- The balanced notches hold the cord in place
- Sterilizable

REF	Item
520000	Retraction Cord Packer, Smooth, 1 pcs
520001	Retraction Cord Packer, Serrated, 1 pcs

ROEKO Retracto

Retracto cords are available as non impregnated. The mechanical use of the retraction cord causes compression. The compression contracts the treated tissue and causes haemostasis.

Item				REF	
				Braided	Twisted
ROEKO Retracto	Non Impregnated	x-fine	1 pcs	511000	523000
ROEKO Retracto	Non Impregnated	fine	1 pcs	511001	523001
ROEKO Retracto	Non Impregnated	medium	1 pcs	511002	523002
ROEKO Retracto	Non Impregnated	thick	1 pcs	511003	-



523002

The background is a solid orange color with several large, overlapping, curved shapes in a lighter shade of orange. These shapes are positioned in the upper left, middle, and lower right areas, creating a dynamic, abstract composition.

2

Bite Registration

SPOTLIGHT



JET BLUE BITE

14

Jet Blue Bite is a material for bite registration. It offers a fast intra oral setting-time and sets rock hard for greater accuracy but is still flexible enough for easy trimming or cutting.

JET BITE

14

Jet Bite is a material for bite registration. The direct placement and quick intra oral setting time as well as the ideal stability of the material make it easy to record the patients occlusion accurately.

JET BLUE BITE



6491



6495

Jet Blue Bite bite registration material flows with a smooth whipped-cream consistency exactly where you want it in a predictable quantity and ideal ribbon pattern. It offers a fast oral setting-time and sets rock hard for greater accuracy but is still flexible enough for easy trimming or cutting. The superior handling characteristics are the key for a perfect occlusal fit on the first try.

- Stays exactly where you want it
- No resistance to closure, no rebound!
- Rock hard for accuracy, yet flexible for easy trimming

REF	Item
6491	Fast Refill 1 × 50 ml Cartridge 1 × 6 pcs Universal Mixing Tip 1 × 6 pcs Spreader Tips
6493	Superfast Refill 1 × 50 ml Cartridge 1 × 6 pcs Universal Mixing Tip 1 × 6 pcs Spreader Tips
6495	Fast Refill 4 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip 1 × 12 pcs Spreader Tips
6497	Superfast Refill 4 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip 1 × 12 pcs Spreader Tips

JET BITE



6400

Jet Bite gives a great occlusal fit on the first try. The patient bites down naturally every time because of Jet Bite's fluffy, mousse-like consistency. There is no need for time-consuming adjustments. Jet Bite flows smoothly - sets rock hard.

- Smooth mousse-like consistency - No resistance to closure, no rebound
- High dimensional stability - no distortion
- Not brittle - easy to grind and cut

REF	Item
6400	Refill 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip 1 × 12 pcs Spreader Tips
6410	Refill, 20 × 50 ml

NOTES

A large grid of small dots for taking notes, covering the majority of the page.

The background is a solid orange color with several large, overlapping, semi-transparent curved shapes in a lighter shade of orange. These shapes are arranged in a way that suggests a stylized letter 'E' or a similar graphic element. The shapes are positioned in the upper and lower portions of the page, leaving a large central area of solid orange.

3

Impression A-Silicone

SPOTLIGHT



AFFINIS®

18

AFFINIS is COLTENEs premium A-Silicone for dental impressions. Apart from the ideal working and setting time, AFFINIS convinces dentists with a high dimensional stability and accuracy. The extraordinary surfactant system results in outstanding flow properties for greater detail reproduction.

PRESIDENT The Original

24

The introduction of PRESIDENT in 1975 – the first A-Silicone in the market – is considered a milestone in the development of high-precision dental impression materials. On the occasion of the 40th anniversary of PRESIDENT, the first A-Silicone line has been relaunched under the new brand name PRESIDENT The Original with improved characteristics.

AFFINIS®

Perfect Impressions

A perfect AFFINIS impression is the key for accurately fitting, laboratory fabricated restorations. It provides the ideal basis for successful teamwork between the dentist and the dental technician.

- Innovative surface affinity
- Optimised flow properties
- Excellent thixotropy and the self-contouring consistency
- Outstanding model fabrication
- Autoclavable



AFFINIS light body



6501

AFFINIS wash materials wet the tooth and gingiva extremely quick and effective. Its distinct surface affinity and flowability virtually eliminates the risk of air voids and distortions.

- Low contact angle
- Excellent flow properties in a moist environment

REF	Item
6501	Refill 2 × 50 ml Cartridge 12 pcs Universal Mixing Tip
6580	Refill, 20 × 50 ml

AFFINIS fast light body



6601

AFFINIS fast light body provides an extremely quick setting time. The mixed material shows adequate stability and does not drip from the tooth.

REF	Item
6601	Refill 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6680	Refill, 20 × 50 ml

AFFINIS regular body



6511

AFFINIS wash materials wet the tooth and gingiva extremely quickly and effectively. Its hydro-colloid-like property ensures that the wash material flows into the narrowest areas without applying pressure. Beside the low contact angle thanks to its structure-breaking property of the surfactant agent, AFFINIS regular body provides a medium wash consistency, an adequate stability and does not drip from the tooth.

REF	Item
6511	Refill 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6582	Refill, 20 × 50 ml

AFFINIS fast regular body



6611

Beside the low contact angle thanks to its structure-breaking property of the surfactant agent, AFFINIS fast regular body provides a medium wash consistency, an adequate stability and sets extremely fast.

REF	Item
6611	Refill 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6682	Refill, 20 × 50 ml

AFFINIS PRECIOUS



6775

AFFINIS PRECIOUS polyvinylsiloxane (PVS) wash materials have a contact angle of only 10°. Its superior flow properties allow the impression material to easily capture all critical surface details of the preparations in a wet environment.

The unique silver and gold pigmentation drastically reduces the light scattering effect, thus improving the visibility of the margins.

Thermo-reactive formula leaves you in control of the total setting time

AFFINIS PRECIOUS wash materials have a variable working time from 5 to 60 seconds. The variable working time is possible due to the thermo-reactive formula that has an accelerated setting after the tray has been seated and the material has reached its final film thickness. This allows the total impression setting time to be in direct relationship to the number of preparations included in the impression. The fewer preparations you have, the less total setting time you need.

REF	Item
6775	Light Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip
6776	Regular Body 2 × 50 ml Cartridge 1 × 12 pcs Universal Mixing Tip

AFFINIS MonoBody



6480

AFFINIS MonoBody is tissue friendly and neutral in taste without side effects or patient discomfort. After a short oral setting time of 2:30 min the AFFINIS impression can be easily removed without residue or distortion.

AFFINIS MonoBody can also be used as the tray material for the one-step putty/wash technique.

REF	Item
6480	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
6482	System 360 Refill, 2 × 380 ml
6547	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6587	Refill, 20 × 75 ml

AFFINIS heavy body



6590

Ideal tray material for the one-step putty/wash technique as well as reproduction of buccal and palatal areas.

- High consistency and optimal thixotropy
- Hydrocolloid-like properties

REF	Item
6485	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
6487	System 360 Refill, 2 × 380 ml
6520	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6585	Refill, 20 × 75 ml
6590	Refill 2 × 50 ml Cartridge 1 × 6 pcs Mixing Tips
6595	Refill, 20 × 50 ml

AFFINIS fast heavy body



6620

The high consistency, rapid build up of pressure and fast working time of AFFINIS fast heavy body produces strong compression, making this tray material ideal for the one-step putty/wash technique. Buccal and palatal areas can also be easily reproduced with AFFINIS fast heavy body.

REF	Item
6620	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
6650	Refill 2 × 50 ml Cartridge 1 × 6 pcs Mixing Tips
6652	Refill, 20 × 50 ml
6685	Refill, 20 × 75 ml

AFFINIS BLACK EDITION



60019554

The black tray material AFFINIS heavy body BLACK EDITION reaches with gold or silver coloured AFFINIS PRECIOUS excellent readability of the impression.

REF	Item
60019554	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
60019774	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
60019776	System 360 Refill, 2 × 380 ml

AFFINIS System 360 putty



6471

AFFINIS System 360 putty has a strong "firm body" which is immediately noticeable once dispensed from the mixing unit with no signs of stickiness. The dentist is comforted with a reliable, automatically-mixed putty material for successful precision impressions.

- Putty available in cartridge
- Prevents distortions in the impression caused by saliva and blood residue.

REF	Item
6471	Putty System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
6472	System 360 Refill, 2 × 380 ml

AFFINIS putty soft / fast



6530

AFFINIS putty soft is a standard putty in the AFFINIS putty range. The compatible consistencies of the base and catalyst ensure for quick and easy mixing.

- Easy mixing consistency
- Prevents distortions in the impression caused by saliva and blood residue.

REF	Item
6530	Putty Soft
	1 x 300 ml Base
	1 x 300 ml Catalyst
	2 x Dosage spoon
6531	Fast Putty Soft
	1 x 300 ml Base
	1 x 300 ml Catalyst
	2 x Dosage spoon

AFFINIS putty super soft



6535

AFFINIS putty super soft is smooth, quick and easy to mix to a homogeneous consistency and with a viscosity similar to heavy body.

- Easy mixing consistency
- Prevents distortions in the impression caused by saliva and blood residue.

REF	Item
6535	Putty Super Soft
	1 x 300 ml Base
	1 x 300 ml Catalyst
	2 x Dosage spoon

PRESIDENT The Original

A small step for dentists. A great step for dentistry.

The introduction of PRESIDENT in 1975 – the first A-Silicone in the market – is considered a milestone in the development of high-precision dental impression materials. Through continuous development, the PRESIDENT brand continues to impress with its remarkable quality.

The new, high pigmented colours of the PRESIDENT the Original line provide better contrast and show every detail, even before the impression is poured. All consistencies offer excellent flow properties in a moist environment.

- High tear strength and elasticity
- Optimised flowing behavior and thixotropy
- Easier-to-mix putties
- New and more contrasting colours
- New consistencies



PRESIDENT The Original Xtra light body



60019937

PRESIDENT the Original Xtra Light Body is a light-flowing material that can easily capture even the smallest marginal details.

REF	Item
60019937	Refill
	2 × 50 ml Cartridge
	1 × 12 pcs Universal Mixing Tip

PRESIDENT The Original light body



60019938

PRESIDENT the Original Light body is characterized by its excellent flow properties and viscosity. As a real all-rounder it is suitable for any impression technique.

REF	Item
60019938	Refill
	2 × 50 ml Cartridge
	1 × 12 pcs Universal Mixing Tip
60019949	Tube Refill
	1 × 53 ml Tube Base
	1 × 53 ml Tube Catalyst

PRESIDENT The Original regular body



60019939

PRESIDENT the Original Regular Body is a medium-flowing material which is very suitable for the double mixing technique due to its stability.

REF	Item
60019939	Refill
	2 × 50 ml Cartridge
	1 × 12 pcs Universal Mixing Tip
60019950	Tube Refill
	1 × 53 ml Tube Base
	1 × 53 ml Tube Catalyst

PRESIDENT The Original MonoBody



60019945

PRESIDENT the Original MonoBody is ideally suited for implant and/or single-phase impressions due to its adapted flow properties. It is available in 75ml cartridges for use with a dispenser and 5:1 cartridges for use in the Automix System.

REF	Item
60019945	System 360 Starter Kit
	1 × 380 ml Cartridge
	1 × 10 pcs Dynamic Mixing Tips
60019946	System 360 Refill, 2 × 380 ml

PRESIDENT The Original heavy body



60019940

PRESIDENT the Original Heavy Body is a viscous tray material available in 75ml cartridges for use with a dispenser and 5:1 cartridges for use in the Automix System.

REF	Item
60019940	Refill 2 × 50 ml Cartridge 1 × 6 pcs Mixing Tips
60019941	Refill 2 × 75 ml Cartridge 1 × 8 pcs Mixing Tips
60019947	System 360 Starter Kit 1 × 380 ml Cartridge 1 × 10 pcs Dynamic Mixing Tips
60019948	System 360 Refill, 2 × 380 ml

PRESIDENT The Original putties



60019942



60019943

PRESIDENT the Original Putties are kneadable, highly viscous tray materials. The remarkable recoverability of Putty Soft and Putty super soft results in highly precise impressions even in areas with strong undercuts.

REF	Item
60019942	Putty 1 × 300 ml Base 1 × 300 ml Catalyst 2 × Dosage spoon
60019943	Putty Soft 1 × 300 ml Base 1 × 300 ml Catalyst 2 × Dosage spoon
60022166	Fast Putty Soft 1 × 300 ml Base 1 × 300 ml Catalyst 2 × Dosage spoon
60019944	Putty Super Soft 1 × 300 ml Base 1 × 300 ml Catalyst 2 × Dosage spoon

NOTES

A large grid of small dots for taking notes, covering the majority of the page.

4

Impression C-Silicone

SPOTLIGHT



Speedex

30

Speedex is a C-Silicone with A-Silicone qualities. Speedex has an elastic recoverability and reproduction, which is otherwise only achieved by A-Silicones. As a result of the low shrinkage, impressions remain dimensionally stable and can be cast for up to 7 days without problems.

Speedex

Bites of passion

Speedex has elastic recoverability and reproduction, which is otherwise only achieved by A-Silicones. As a result of the low shrinkage, impressions remain dimensionally stable and can be cast for up to 7 days without problems. Furthermore this premium C-Silicone shows an outstanding wettability and excellent flow characteristics. Saliva and blood are displaced by Speedex medium and have no negative side effects.

The Speedex system gives the advantage to work with only one activator for 4 consistencies.



Speedex light body



4980

The Speedex system gives the advantage to work with only one activator for 4 consistencies. Moreover, Speedex offers physical properties that approximate to those of addition cured silicones - polyvinylsiloxanes:

- 7 days dimension stability
- Controllable working time
- Surface active
- Moisture control
- Optimized precision

REF	Item
4980	Light Body, 1 × 140 ml

Speedex medium



4986

Speedex medium has outstanding wettability and excellent flow characteristics. Saliva and blood are displaced by Speedex medium and have no negative side effects. With the variable setting time of Speedex medium the important mucosal movements, which are necessary for the full and partial denture technique, can be taken very precisely. A perfect impression is the key for an accurate fitting denture.

REF	Item
4986	Medium, 1 × 140 ml

Speedex putty



4970

Speedex putty is a Putty material with ideal consistency which perfectly fits to Speedex wash materials.

REF	Item
4970	Putty, 1 × 910 ml

Speedex putty soft



60023572

Speedex putty soft is smooth and easy to mix to a homogeneous consistency. It fits perfectly to Speedex wash materials.

REF	Item
60023572	Putty soft, 1 × 910 ml

Speedex Activator



4990

Universal Activator for Speedex putty, Speedex putty soft, Speedex light body and Speedex medium.

REF	Item
4990	Universal Activator, 1 × 60 ml

Pressure Spot Indicator



4440

The polysiloxane-based material for the location of high spots in dentures.

- Recognition of pressure points on new or existing
- partial and total dentures
- Determining the need for relining
- Testing of the exact fitting of inlays, crowns, bridges and cast brackets

REF	Item
4440	PSI (Base & Catalyst) 1 × 55 ml Tube Base 1 × 20 ml Tube Catalyst
4445	PSI Activator, 1 × 20 ml

NOTES

A large grid of small dots for taking notes, covering the majority of the page.

The background is a solid orange color with several large, overlapping, curved shapes in a lighter shade of orange. These shapes are positioned in the upper left, middle right, and lower left areas, creating a dynamic, abstract composition.

5

Alginate

SPOTLIGHT



HydroPrint Premium

36

HydroPrint Premium is an alginate based Impression material. It is easy to mix and dust free. Its peppermint flavour and low viscosity make the experience more pleasant for patients.

HydroPrint Premium



60023007

HydroPrint Premium is suitable for anatomic impressions. Anatomic impressions have different indications. They are needed for study-, analysis-, planning- and working models. Important properties of the alginates are the flowability, detail reproduction, dimensional stability, elastic recovery and disinfection.

The outstanding elastic properties of HydroPrint Premium are characterized by the high tear strength of the material.

REF	Item
60023007	Regular set, 1 × 454 g

REF	Accessories
H00564	Mixing Bowl Medium, 1 pcs
H00565	Mixing Bowl Large, 1 pcs
H00566	Mixing Bowl Extra Large, 1 pcs
H06012	Spatula, 3 pcs

NOTES

A large grid of small dots for taking notes, covering the majority of the page.

6

Prosthetics Accessories

SPOTLIGHT



Impression Trays 41

PRESIDENT Trays are single use plastic impression trays. A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray. The PRESIDENT Trays AC are also autoclavable.

Mixing Tips 43

Universal Mixing Tip - a universal mixer developed especially for all light and regular body materials as well as bite registration and gingival mask materials from COLTENE. The reduced extrusion force allows for a comfortable and accurate dispensing of material.

COLTENE Adhesive



4405

A thin layer of COLTENE Adhesive ensures adhesion of tray material to the tray sides.

COLTENE Adhesive

COLTENE Adhesive is a universal adhesive for silicone-based dental impression materials and serves as an adhesive bonding agent between the impression material and the impression tray. Separation of the impression material during removal from the mouth can therefore be prevented

COLTENE Adhesive AC

Coltène Adhesive is an autoclavable adhesive for all AFFINIS impression materials and serves as an adhesive bonding agent between this material and the PRESIDENT disposable impression trays AC. Separation of the impression material during removal from the mouth can therefore be prevented. The special spray nozzle provides application with a uniform, thin layer without becoming clogged. Coltène Adhesive AC was specially developed for the autoclavable system (AFFINIS impression materials, PRESIDENT disposable impression trays AC).

REF	Item
4405	Adhesive Autoclavable, 1 × 10 ml
4410	Adhesive, 1 × 10 ml

Impression Trays

A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray. The anatomic form relieves the fitting and ensures a perfect impression. Autoclavable and non autoclavable.

- Anatomically correct: A statistical database accumulated over 20 years has enabled Coltène / Whaledent to define the ideal form for each tray system
- Built-in retention: The President Tray System's design aggressively retains all types of impression material.
- Transparent: Easy and precise positioning of the tray.
- Modifiable form: The tray can individually be adjusted with heat.
- Single use: No risk of contamination, no time-consuming cleaning.

Item			REF					
			Upper Jaw	Upper Jaw	Upper Jaw	Lower Jaw	Lower Jaw	Lower Jaw
			S	M	L	S	M	L
Impression Trays	12 pcs	Autoclavable	C6954	C6952	C6950	C6955	C6953	C6951
		Not autoclavable	C6905	C6903	C6901	C6906	C6904	C6902



C6950



C6901

Impression Trays Sets



C6955

REF	Item
C6910	Assortment Impression Trays 1 × 12 pcs Full Upper Jaw Large 1 × 12 pcs Full Lower Jaw Large 1 × 12 pcs Full Upper Jaw Medium 1 × 12 pcs Full Lower Jaw Medium 1 × 12 pcs Full Upper Jaw Small 1 × 12 pcs Full Lower Jaw Small 1 × Quadrant UJ le. / LJ ri. 1 × Quadrant UJ ri. / LJ le. 1 × Anterior UJ / LJ
C6959	Assortment President Trays, autoclavable 1 × 12 pcs Full Upper Jaw Large, autoclavable 1 × 12 pcs Full Lower Jaw Large, autoclavable 1 × 12 pcs Full Upper Jaw Medium, autoclavable 1 × 12 pcs Full Lower Jaw Medium, autoclavable 1 × 12 pcs Full Upper Jaw Small, autoclavable 1 × 12 pcs Full Lower Jaw Small, autoclavable 1 × Quadrant UJ le. / LJ ri., autoclavable 1 × Quadrant UJ ri. / LJ le., autoclavable 1 × Anterior UJ / LJ, autoclavable

Impression Trays Quadrant & Anterior



C6956

REF	Item
C6907	Quadrant UJ le. / LJ ri., 12 pcs
C6908	Quadrant UJ ri. / LJ le., 12 pcs
C6909	Anterior UJ / LJ, 12 pcs
C6956	Quadrant UJ le. / LJ ri., autoclavable, 12 pcs
C6957	Quadrant UJ ri. / LJ le., autoclavable, 12 pcs
C6958	Anterior UJ / LJ, autoclavable, 12 pcs

Mixing Devices



6985

Dispensers and other mixing devices for cartridges

REF	Item
6985	Plunger 50/75 ml, 1 pcs
60011147	Sympress II Europe, 1 pcs
60011150	Sympress II UK, 1 pcs

REF	Accessories
6900	Dispenser 50/75 ml, 1 pcs

Mixing Tips



6552

One new UNIVERSAL MIXING TIP - Always the right choice

Mixing Tips and Oral Tips for different wash and tray materials which ensure a homogenous ratio between base and catalyst.

With up to 55 % less extrusion force during the application of wash material the new UNIVERSAL MIXING TIP ensures an easier work with dispensers.

- For all 50 ml cartridges from COLTENE
- For all COLTENE light and regular body, bite registration and gingival mask materials
- Simple and precise application due to reduced extrusion forces
- Practical, flexible and economical

REF	Item
6162	Mixing Tip Turquoise, 40 pcs
6163	System 360 Fixation Ring yellow, 1 pcs
6210	Oral Tip white, 100 pcs
6223	Spreader Tip, 100 pcs
6552	Mixing Tip turquoise, 40 pcs
6555	Oral Tip yellow, 100 pcs
60019663	Universal Mixing Tip Refill, 40 pcs
60019694	Universal Mixing Tip Refill, 100 pcs

Other Accessories



4422

Coltopress, the squeezing aid for tubes, empties tubes without waste.

COLTENE Mixing Spatula supports a consistent and homogeneous mixing.

The new COLTENE cutting knife is a universal version for both left- and right handed customers offering the following:

- New ergonomic design which perfectly fits to the size of the hand
- Special undercuts in the knife handle avoid slipping of fingers
- Sharp blade makes cutting easy and efficient
- New packaging and colour
- Autoclavable

REF	Item
4420	Cutting Knife, 1 pcs
4422	Coltopress, 1 pcs
60010943	SYMPRESS-CASE, 1 pcs
4412	Mixing Pad, 1 pcs
4414	Mixing Pad, 10 pcs

SPOTLIGHT



GI-MASK

46

The silicone-based gingival mask for precise reproduction of gingival contours around the margins maintains its accuracy for days and right up to the moment until the restoration is finished.

Lab putty

49

Lab-Putty is the outstanding polysiloxane laboratory modelling material with a 62 - 67 Shore A hardness. Ideal for model repairs, isolation, denture heat protection, keys, block outs and matrixes for the formation of general occlusion and functional guides.

GI-MASK

No compromise in function and aesthetics

The silicone-based gingival mask for precise reproduction of gingival contours around the margins maintains its accuracy for days and right up to the moment until the restoration is finished.

The GI-MASK family consists of the classical C-Silicone handmix, the A-Silicones GI-MASK Automix New Formula and GI-MASK 3D in cartridges.

GI-MASK

The advantage of GI-MASK is the exclusive accessory system which offers accurate dispensing, neat mixing and simple application. Its short, time-saving setting time makes GI-Mask a unique aid for dental laboratories.

GI-MASK Automix New Formula

GI-MASK Automix New Formula provides accurate reproduction of soft tissue structure. The A-Silicone based gingival mask provides all information for implant supra-structures and frameworks. The aesthetic coloring of the gingival mask supports the clinical situation. It facilitates the preparation of esthetically excellent and perfectly functioning crown and bridge work. In addition it helps to prepare ideal conditions for periodontal hygiene.

GI-MASK 3D

GI-MASK 3D has very similar properties to GI-MASK Automix New Formula. In addition, the gingival mask material can be scanned without prior surface treatment.



GI-Mask 3D



6002251

The scannable GI-MASK 3D allows digitising reproduced gingiva profiles in the preparation region without scan spray. The use of a gingival mask is particularly important in the fabrication of implant-supported supra-structures to achieve the best fit for dentures.

- Scannable without surface treatment
- Excellent flowability
- Extremely high elasticity and tear strength

REF	Item
6002250	STARTER KIT
	2 x 50 ml Cartridge
	1 x 50 ml Universal Separator
	1 x Trimmer for Silicones
	1 x 12 pcs Universal Mixing Tip
1 x 12 pcs Oral Tips	
6002251	Refill
	2 x 50 ml Cartridge
	1 x 12 pcs Universal Mixing Tip
	1 x 12 pcs Oral Tips

GI-MASK Automix New Formula



8060

A-Silicone GI-MASK Automix New Formula facilitates the preparation of esthetically excellent and perfectly functioning crown and bridge work. In addition it helps to prepare ideal conditions for periodontal hygiene.

It is particularly suitable for single crowns, bridges, implants, combination and telescopic work.

REF	Item
8060	Starter Kit
	2 x 50 ml Cartridge
	1 x 50 ml Universal Separator
	1 x Trimmer for Silicones
	1 x 12 pcs Universal Mixing Tip
	1 x 12 pcs Oral Tips
8065	Refill
	2 x 50 ml Cartridge
	1 x 12 pcs Universal Mixing Tip
	1 x 12 pcs Oral Tips

GI-MASK



8030



8036

The GI-MASK is the C-Silicone based gingival mask. It has been conceived for the quick and accurate reconstruction of the complete contours of the gingiva around the preparation on the model. It is particularly suitable for:

- single crowns
- bridges
- implants
- combination and telescopic work

REF	Item
8030	Starter Kit 1 × 150 ml Base 1 × 18 ml Activator 1 × 50 ml Universal Separator
8036	Activator, 1 × 18 ml
8037	Refill 1 × 150 ml Base 1 × 18 ml Activator

Lab putty

COLTENEs Lab-Putty was especially designed for the use in laboratory. It can be used for model repairs, isolation, denture heat protection, keys, block outs and matrixes for formation of general occlusion and functional patterns. Lab-Putty hard is characterized by its high heat-resistance (300°C) and its high end hardness (90 Shore A).



Lab putty / Lab putty hard



8805



8840

Lab-Putty is the outstanding polysiloxane laboratory modelling product material with 63 - 70 Shore A for model repairs, isolation and denture heat protection, keys and block outs and matrixes for formation of general occlusion and functional patterns.

A smooth acrylic surface is guaranteed with the use of Lab Putty hard for relining, keys and repair models. Separate isolation is eliminated. Its high heat-resistance (300° C) and perfect dimensional stability makes it suitable for soldering models (90 Shore A).

REF	Item
8805	Lab Putty 1 × 900 ml Base 1 × 40 ml Activator 1 × Dosage spoon
8815	Lab Putty 1 × 5000 ml Base 5 × Dosage spoon
8840	Activator, 1 × 40 ml
8880	Lab Putty Hard 1 × 2600 ml Base 5 × Dosage spoon
8882	Lab Putty Hard 1 × 900 ml Base 1 × 40 ml Activator 1 × Dosage spoon

ACCU-TRAC

The ACCU-TRAC system makes the fabrication of precision model and dies easier. It creates a precise base for your prosthodontic work. The model and base are poured as one, eliminating the additional work when using pins and sleeves. Time and material can be saved without sacrificing precision.



ACCU-TRAC®



ATK25



ATA20

The ACCU-TRAC® System makes the fabrication of precision models and dies easier: The model and base are poured as one, eliminating the additional work when using pins and sleeves. You save time and material, without sacrificing precision.

- One universal tray size, accomodates all maxillary and mandibular sizes, including quadrants
- Low base height, model fits any articulator
- One pour Technique saves time and materials
- Economical, all components reusable
- Numerical die location, immediate recognition of dies, easy die repositioning

REF	Item
AT6	Introductory Kit 1 × 6 pcs Precision Die Trays 1 × 6 pcs Pouring/Ejector Bases 1 × 6 pcs Articulation Spacers 1 × 25 pcs Articulation Keepers
AT72	Economy Pack, 72 pcs
ATR12	Acrylic frame Refill 1 × 12 pcs Precision Die Trays 1 × 50 pcs Articulation Keepers
ATR144	Economy Pack, 144 pcs
ATA20	Locking Arms Refill, 20 pcs
ATE12	Pouring/Ejector Bases Refill, 12 pcs
ATK25	Articulation Keepers Refill, 25 pcs
ATK300	Articulation Keepers Refill, 300 pcs
ATS20	Articulation Spacers, 20 pcs
ATT12	Acrylic frame Refill, 12 pcs
ATT576	Acrylic frame Refill, 576 pcs
ATN50	Transport System 1 × 50 pcs Disposable Carrying Cases 1 × 50 pcs Plastic Retainers

PINDEX

Uncompromising model and die precision which assures parallel channels for quick, easy and reliable pin placement. The laser drilling machine offers more accurate drill positioning for better die registration and the visible beam reduces eye fatigue. The focus-free laser outlasts conventional bulbs, and manual drilling to pre-set depth is possible.





P113L

The PINDEX System is considered to be the «benchmark» used in model making. It provides a simple way to fabricate precise parallel drill channels. Laser Drilling Machine- Assures parallel channels for quick, easy, and reliable pin placement.


- More accurate drill positioning for better die registration
- Visible beam in all lighting for reduced eye fatigue
- Focus-free laser outlasts conventional bulbs
- Manual drilling to pre-set depth System Accessories
- Durable, high quality carbide drill

REF	Item
P113LCE	Laser Drilling Machines 230V Europe 1 x Hand drill 1 x Set of flexible socketformer 2 x Tools for drill changing 1 x Socket wrench 1 x Pin introduction set
P113LUK	Laser Drilling Machines 230V Uk 1 x Hand drill 1 x Set of flexible socketformer 2 x Tools for drill changing 1 x Socket wrench 1 x Pin introduction set
PX112S	Medium Pins & Split Sleeves, 1000 pcs
PX112T	Medium Pins/Sleeve Brn, 1 pcs
PX113L	Laser System 115V 1 x Hand drill 1 x Set of flexible socketformer 2 x Tools for drill changing 1 x Socket wrench 1 x Pin introduction set
PX114S	Long Pins, 1000 pcs
PX115	Short Pins & Plastic Sleeves, Bulk, 1000 pcs
PX116	Long Pin Kit, 100 pcs
PX117	Short Pins & Plastic Sleeves, 100 pcs
PX119	Carbide Drills, 3 pcs
PX120	Reamer, 1 pcs
PX125	Laser Retro-Fit Kit, 1 pcs
PX144	Self-Articulating Pins & Plastic Sleeves, 1000 pcs
PX145	Self-Articulating Pins & Plastic Sleeves, Bulk, 100 pcs
PX145S	Plastic Sleeves For Articulating Pins, Bulk, 1000 pcs

HYGENIC Waxes

Small details, big impact

The HYGENIC wax portfolio offers different wax sizes and shapes for various indications. Users can chose between bite wafers, u-shaped occlusal rim wax, square wax ropes, round wax strips, wax sticks which can be used for undercuts of individual trays, bite registration, bite plane appliances etc.



HYGENIC yellow bite wax



H00807

Economical alternative to beeswax with many of the same characteristics.

REF	Item
H00807	HYGENIC Yellow Bite Wax, 1 × 450 g
H00813	HYGENIC Yellow Bite Wax, 1 pcs

HYGENIC bite-wafers



H00825

This bite-mask contains a thin aluminium foil within two wax layers. U-shaped and designed for a fast precise retain of the situation.

REF	Item
H00825	HYGENIC Bite-Wafers, Double, 24 pcs
H00839	HYGENIC Bite-Wafers, Yellow, Soft, 50 pcs
H03365	HYGENIC Bite-Wafers, Green, Hard, 50 pcs

HYGENIC boxing wax



H00815

Designed for boxing impressions, withstands extreme heat given off by dental stone during setting.

REF	Item
H00815	HYGENIC Boxing Wax, Regular, 1 x 450 g
H00816	HYGENIC Boxing Wax, Extra Fine, 1 x 450 g
H00837	HYGENIC Boxing Wax, Regular, 1 pcs
H00838	HYGENIC Boxing Wax, Extra Fine, 1 pcs

HYGENIC medium-soft no. 3 pink/red wax



H00805

Easily softened and cooled to minimize shrinkage. Soft modeling.

REF	Item
H00806	HYGENIC Medium Soft,# 3 Pink Wax, 1 x 450 g
H00805	HYGENIC Red Wax, 1 x 450 g
H00811	HYGENIC Red Wax, 1 x 2.25 kg
H00812	HYGENIC Medium Soft,# 3 Pink Wax, 1 x 2.25 kg

HYGENIC u-shaped occlusal rim wax



H04030

Recessed base aids adaption to base plate. Can be reshaped at room temperature to fit individual arch. If it's too rigid, just warm the wax slightly in warm water.

REF	Item
H00822	HYGENIC U-Shaped Occlusal Rim Wax, Standard Package, 100 pcs
H04030	HYGENIC U-Shaped Occlusal Rim Wax, Lab Package, 21 pcs

HYGENIC orthodontic tray wax



H00827

Sticky, soft and narrow wax ropes.

REF	Item
H00827	HYGENIC Orthodontic Tray Wax, White, 48 pcs

HYGENIC wax ropes



H00819

Sticky, soft and narrow wax ropes.

REF	Item
H00817	HYGENIC Round Wax Strips, Red, 80 pcs
H00818	HYGENIC Round Wax Strips, White, 80 pcs
H00819	HYGENIC Square Wax Ropes, Red, 1 pcs
H00820	HYGENIC Square Wax Ropes, White, 1 pcs

HYGENIC wax sticks



H00826

Readily formed without heating, sticks are softer and tackier than „Round Wax Strips“ to minimize tissue displacement.

REF	Item
H00826	HYGENIC Wax Sticks, Scented, 48 pcs

HYGENIC extra tough wax



H00803

Durable, all-purpose wax that is a carefully prepared blend of exclusive waxes. Tough, pliable and sufficiently tacky to hold teeth well. Takes polish and boils out.

REF	Item
H00803	HYGENIC Extra Tough Wax, 1 × 450 g
H00809	HYGENIC Extra Tough Wax, 1 pcs

Laboratory Accessories

COLTENEs laboratory accessories are made for laboratory needs. The HYGENIC Flexibole® mixing bowls, as well as the nylon spatula are designed so that the mass to be mixed does not stick to them.



Laboratory Accessories



H00563


Plaster, stone, investments, and alginates will not stick to these self-cleaning Flexibole® mixing bowls or the nylon spatula. More flexible than rubber bowls for more efficient mixing. The high-quality trimmer for silicones is used to remove any excess gingival mask material during the restoration process.

REF	Item
H00563	Mixing Bowl Small, 1 pcs
H00564	Mixing Bowl Medium, 1 pcs
H00565	Mixing Bowl Large, 1 pcs
H00566	Mixing Bowl Extra Large, 1 pcs
H00567	Mixing Bowl Jumbo, 1 pcs
H00676	3/16 x 3/32 x 50 Ft Amber Tube, 1 pcs
H00682	1/4 x 1/16 x 50 Ft Amber Tube, 1 pcs
H00684	1/4 x 1/16 x 50 Ft Black Tube, 1 pcs
H00685	1/4 x 3/32 Amber, 1 pcs
H00688	1/4 x 1/8 x 50 Ft Amber Tube, 1 pcs
H00697	5/16 X 3/32 X 50 Ft Amber Tube, 1 pcs
H00757	Saliva Eject Basketscreen A238, 1 pcs
H06012	Spatula, 3 pcs
H06541	100' Clean-Flex 1/2Hve Tubing, 1 pcs
8075	Trimmer for Silicones, 1 pcs

NOTES









LOCATIONS

CA	<p>SciCan Ltd. 1440 Don Mills Road Toronto, Ontario, M3B 3P9 / Canada T +1 416 445 1600 F +1 416 445 272</p>	CH	<p>Coltene/Whaledent AG Feldwiesenstrasse 20 9450 Altstätten / Switzerland T +41 71 757 53 00 F +41 71 757 53 01 info.ch@coltene.com</p>
CN	<p>Coltene/Whaledent China Trading CO. Ltd. Rm. 1103/1104, Bldg. B2, Wu Dong Plaza #9 Chegongzhuang Street, West District Beijing 100044 / P.R. China T +86 10 6212 3251 F +86 10 6212 3252 info.cn@coltene.com</p>	DE	<p>Coltene/Whaledent GmbH & Co. KG Raiffeisenstraße 30 89129 Langenau / Germany T +49 7345 805 0 F +49 7345 805 201 info.de@coltene.com</p>
ES PT	<p>Coltene Iberia S.L.U. c/Doctor Fleming, 3 - Bajo A 28220 Majadahonda, Madrid / Spain T +34 91 63 86 492 F +34 91 63 86 279 info.ib@coltene.com</p>	FR	<p>Coltene/Whaledent S.a.r.l. Synergie Park 3 avenue Louis Néel 59260 Lezennes / France T +33 320 64 64 00 F +33 320 64 64 03 info.fr@coltene.com</p>
FR	<p>Micro-Mega SA 12 rue du Tunnel 25000 Besancon / France T +33 381 54 42 42 F +33 381 54 42 39</p>	IN	<p>Coltene/Whaledent Private Ltd. 106/107, Hallmark, Vasant Oscar L.B.S Marg, Mulund West Mumbai 400080, Maharashtra / India T +91 22 25923626 F +91 22 25623165 info.in@coltene.com</p>
IT	<p>Coltene Italia S.r.l c/o Polo Tecnologico della Brianza Via Alessandro Volta, 94 - Cond. Bianchi 20832 Desio (MB) / Italy T +39 0362 62 67 32 F +39 0362 63 25 80 info.it@coltene.com</p>	JP	<p>Coltene Japan LLC 2-2-28 2nd Floor, Nishikicho 190-0022 Tachikawa Tokyo / Japan T +81 42 595 6945 F +81 42 595 6946 info.jp@coltene.com</p>
UK	<p>Coltene/Whaledent Ltd. The President Suite Kendal House, Burgess Hill Victoria Way West Sussex, RH15 9NF / UK T +44 1444 235486 F +44 1444 870640 info.uk@coltene.com</p>	US	<p>Coltene/Whaledent Inc. 235 Ascot Parkway Cuyahoga Falls, Ohio 44223 / USA T +1 330 916 8800 F +1 330 916 7077 info.us@coltene.com</p>



© COLTENE - www.coltene.com

