



Radiography Holders	216
Digital Sensor Holders	217
Film & Phosphor Plate Holders	222
Radiography Aprons	225
Film Mounts	226
Small Equipment	228
Kodak Film & Chemicals	230







X-Ray Accessories

Radiography Holders

Which radiography holder is right for your equipment, and for you?

for Digital Sensors	Universal	Autoclavaki	Disposable	Paralleling	Bisecting de	Anteriors	Posteriors	Bitewings	Endodontic	Page #
XCP-DS [®] Sensor Holders		ullet		ullet		ullet	ullet			219
 XCP-DS FIT™ Sensor Holders									٠	217
 XCP-0RA™ One Ring & Arm System										217
Uni-Grip® 360 Sensor Holder						•				218
										218
				•						221
for Film and Phosphor Plates										
XCP® Film & PSP Holder										222
XCP-ORA [™] One Ring & Arm System										222
EndoRay®II Film & PSP Holder									٠	224
 Uni-Bite™ Film & PSP Holder	٠									224
Uni-Grip® Film & PSP Holder										222
Snap-A-Ray® Xtra Film & PSP Holder										222
		\bullet								224
RINN® Bitewing Loops										223

Digital Sensor Holders



XCP-ORA[™] - One Ring & Arm Positioning System

The XCP-ORA - One Ring & Arm Positioning System reduces the number of components needed for accurately positioning radiographic media with RINN[®] XCP[®] or XCP-DS[®] bite-blocks for anterior, posterior and bitewings to just one arm, one ring.

XCP-ORA One Ring & Arm, 1-pack XCP-ORA One Ring & Arm, 2-pack 550771 550772

Product Code







XCP-DS FIT™

Reduce loading time during your digital radiograph procedures with the new XCP-DS FIT Universal Autoclavable Sensor Holder.

- Universal design provides a custom fit to your sensor
- One holder works with both size 1 and size 2 sensors for less inventory
- Autoclavable for lower cost per use
- Rigid sensor support keeps sensor positioned accurately
- Silicone-cushioned edges help provide patient comfort and prevent sensor cover perforations
- Use with the XCP-ORA 3-in-1 Positioning System for easy paralleling

XCP-DS FIT Hygiene Kit

 XCP-ORA Arm and Ring, 8 XCP-DS FIT bite-blocks
(2 each: anterior, posterior, horizontal bitewing, vertical bitewing)

XCP-DS FIT Endo Kit

• 1 x Endo Ring, 1 Endo Arm,

• 2 x XCP-DS FIT Endo Bite-blocks

XCP-DS FIT Complete Kit

559909

559900

559908


• XCP-ORA Arm and Ring, Endo Ring and Arm, 10 XCP-DS FIT bite-blocks (2 each: anterior, posterior, horizontal bitewing, vertical bitewing, endo)

XCP-DS FIT Bite-block Refill Packs, 2 bite-blocks per pack

Anterior	559901
Posterior	559902
Horizontal Bitewing	559903
Vertical Bitewing	559904
Endo	559905



Digital Sensor Holders



Product Code Uni-Grip[®] 360 Universal Digital Disposable Sensor Holder

Adjustable bite-block for anterior, posterior, horizontal bitewings, or vertical bitewings in the one holder, allowing the user to take a very fast full mouth series without the need to use multiple bite-blocks or holders. Universal fit for all digital sensors which locks in any size sensor cover for maximum stability - no need for tight fitting bags that make sensor removal hard.

Uni-Grip 360 Disposable Holder for Digital Sensors Includes:

- 50 x Uni-Grips
- 50 x Digital Sensor Barriers

Uni-Grip 360 Combo Kit for Digital Sensors Includes:

550252

550052

- 2 x 550052 kits (50 Uni-Grips and SD Digital Sensor Barrier)
- 1 x Uni-GripAR

Uni-GripAR® 3-in-1 Arm & Ring Positioning System

Works with any Uni-Grip holder. Colour-coded for easy assembly & fully autoclavable.

Uni-GripAR Positioning Arm

550100

Snap-A-Ray[®] DS

Allows ability to quickly load and place digital sensor holders for accurate anterior and posterior periapical radiographs. An optional arm and ring is available for easy alignment.

Snap-A-Ray DS 1 pack (does not include the optional arm and ring)550321Snap-A-Ray DS 3 pack (does not include the optional arm and ring)550323Arm & Ring for Snap-a-Ray DS & DS Endo550322

Snap-A-Ray[®] DS Endo

Similar to the Snap-A-Ray DS with a cut-out bite surface to accommodate rubber dam clamps.

Snap-A-Ray DS Endo 1 pack,	550421
(does not include the optional arm and ring)	
Snap-A-Ray DS Endo 3 pack,	550423
(does not include the optional arm and ring)	
Arm & Ring for Snap-A-Ray DS & DS Endo	550322









Digital Sensor Holders

Product Code



Suni*

Acteon*

XDR*

Sigma Biomedics*

Cygnus Black*

551604

551020

551800

551010



X-Ray Accessories

Digital Sensor Holders

Product Code

RINN® XCP Kits For Digital Sensors (continued)

Each kit is suitable for a number of different digital sensors. Please refer to the table below.

XCP-DS Kit Name	Other systems the XCP-DS® Kit is suitable for
Schick CDR*	Schick CDR PlusWire, CDR Wireless, Belmont, BelSensor
Gendex e2v*	Gendex VisualiX/GX-S/HDI, Accent, Acteon Sopro, Cyber Medical XDR, Cygnus, Dentrix, Dürr, e2v, HDX, Mediadent, Progeny MPS/MPSe, and Visiodent systems
Trophy Kodak RVG 5000/6000*	RVG Access, RVG Ultimate
Trophy Kodak RVG 5*	RVG ui, RVG Reference, RVG Gold
Suni*	Suni metal-style only; not for Suni-Ray or Dr. Suni Plus Also suitable for CamSight, DMD, Escape, Lightyear, and Multimedia

For full details and details of XCP-DS replacement baskets and bitepieces for specific digital sensors, please contact your Local DENTSPLY Territory Manager.

RINN® Digital Sensor Loops for Bitewing Radiographs

RINN Digital Sensor Loops are designed to slip over digital sensors, providing a handy tab to correctly position sensors for bitewing radiographs.

Sensor Size	Size 0	Size 1	Size 2
Schick, Dr. Suni Plus, Sirona XIOS*	559011	559012	559021
Gendex HDI, Instrumentarium, * Schick Wireless, Kodak RVG 5000/6000,* Suni Ray*		559012	559022
DEXIS*			559022
DEXIS Platinum*		559021	
Kodak 5100/6100*		559021	559023
Dent-X, Gendex eHD, RVG 5*		559021	559023



*Schick, Gendex, Kodak, Suni, Sirona XIOS, Instrumentarium, DEXIS, DEXIS Platinum and Dent-Xare not registered trademarks of the DENTSPLY International Group of Companies.

Digital Sensor Holders

Product Code

RINN[®] Digital Sensor Loops for Bitewing Radiographs (continued)

X-Small Sensor Loop, (Pkt 500)	559011
Small Sensor Loop, (Pkt 500)	559012
Large Sensor Loop, (Pkt 500)	559021
X-Large Sensor Loop, (Pkt 500)	559022
XX-Large Sensor Loop, (Pkt 500)	559023
Assorted Pack (50 of each size)	559099



XCP-DS[®] Disposable Cord Holder

Secure virtually any digital sensor cord to the RINN® indicator arm quickly and easily.

XCP-DS Disposable Cord Holder, (Pkt 50)

550613



Universal Digital Sensor Covers

Protect your sensor with RINN's premium quality sensor covers.

- Fits all #1 and #2 size sensors
- Rounded edges for comfort

Universal Digital Sensor Covers, (Pkt 500)



X-Ray Accessories

Film and Phosphor Plate Holders







Uni-Grip[®] Disposable Radiography Holder for Film/PSP

Accommodates ANY #1 or #2 radiograph media; attaches in seconds. Reposition bitepieces in just a moment. Colour-coded channels ensure correct positioning. Compact size and rounded edges are easy on the patient.

550055 Uni-Grip Disposable Holder for Film & Phosphor Plates

- 50 x Uni-Grips
- 50 x Digital Sensor Barriers

550255 Uni-Grip Combo Kit for Film and Phosphor Plates Includes:

- 2 x 550055 kits (50 Uni-Grips and SD Digital Sensor Barrier)
- 1 x Uni-GripAR

Uni-GripAR[®] 3-in-1 Arm & Ring Positioning System

Works with any Uni-Grip holder, colour-coded for easy assembly & fully autoclavable.

Uni-GripAR Positioning Arm

550100

Product Code



Snap-A-Ray[®] Xtra - Film & Phosphor Plate Holder

The Snap-A-Ray Xtra Film & Phosphor Plate Holder brings an extra measure of efficiency and comfort to anterior and posterior bisecting angle technique radiographs.


Snap-A-Ray Xtra Film & Phosphor Plate Holder, 1-pack	540501
Snap-A-Ray Xtra Film & Phosphor Plate Holder, 3-pack	540503

XCP[®] Kits for Film

Although the colour-coded XCP Arms and Rings have been superseded by the XCP-ORA™ (One Ring&Arm), they are still available to purchase along with the colour-coded XCP bite-blocks for film.

XCP Kit 2000 with Bitewing

- Includes:
- 4 x Arms: 1 Anterior, 2 Posterior, 1 x Bitewing
- 4 x Aiming Rings: 2 Anterior/Bitewing, 2 x Posterior
- 17 x Bite Blocks: 4 Anterior, 8 Posterior, 5 x Bitewing (2 x #2H, 1 x #0/1H, 1 x #2V, 1 x #1V)
- 1 x Instruction Poster



Film and Phosphor Plate Holders





•••••		Product Code
XCP [®] Kits for Film	(continued)	
Arms sold individually	Anterior (Blue)	540857
	Posterior (Yellow)	540858
	Bitewing (Red)	540927
	Endodontic (Green)	550597
Rings sold individually	Anterior (Blue)	540865
	Posterior (Yellow)	540860
	Bitewing (Red)	540934
	Endodontic (Green)	550598
Bite-blocks	Anterior (x25) (Blue)	540861
	Posterior (x25) (Yellow)	540862
	Ant/Post (x10/10) (Blue/Yellow)	540863
	Endodontic (x3) (Green)	540935
Bitewing	#2 Horizontal (x6) (Red)	540929
	#1 Horizontal (x6) (Red)	540933
	#2 Vertical (x6) (Red)	542042
	#1 Vertical (x6) (Red)	540922

Bitewing Bite-blocks for Kodak SURESOFT Film Packets

Dual-slot design accommodates both SURESOFT and standard #2 film packets.

XCP Horizontal Bitewing Bite-block for Kodak SURESOFT Film Packets (Pkt 2)	540629
XCP Vertical Bitewing Bite-block for Kodak SURESOFT Film Packets (Pkt 6)	540630

Disposable Bite-Blocks, Bitewing Holders & Mouth Props



Disposable XCP [®] Bite-blocks, 100 pieces #0, #1 or #2 film	540870
Bitewing Loops, 500 pieces #1 film	141026
Bitewing Loops, 500 pieces #2 film	141027
Bitewing Tabs, 500 pieces	140500
EZ-Prop® Mouth Props, Styrofoam x 100 pieces accommodates dental floss	540900



X-Ray Accessories

Film and Phosphor Plate Holders



EndoRay® II Film Holder

Prevent cone cutting with precise angulation and alignment of film to the x-ray cone. Fits over files, clamps and rubber dams.

EndoRay II, 2 pack	540303
Replacement Parts (1 per pack)	
EndoRay II Arm	540310
EndoRay II Basket	540311
EndoRay II Ring	540312



Uni-Bite[™] Film Holder

Take perfect radiographs with this universal film holder. Great for anterior bitewings and digital phosphor plates.

Uni-Bite, (Pkt 1)

540401

Product Code



RINN[®] FLIP-RAY[™] -Universal Film & Phosphor Plate Holder

The RINN FLIP-RAY's rotating media clip lets the clinician take bitewing radiographs horizontally or vertically with a simple flip. The aiming ring allows the clinician to capture accurate radiographs the first time, every time.

RINN FLIP-RAY Universal Bitewing Holder, (Pkt 5)

Radiography Aprons

Product Code

Protective Aprons

Soothe-Guard Air[™] lead-free aprons have all the features of our traditional Soothe-Guard[®] aprons, plus an innovative lead-free, alloy lining which provides protection equal to lead with less weight, in fact, up to 37% lighter than lead-lined aprons.

Soothe-Guard Air Lead-Free Aprons, 0.3mm Lead Equivalency

Adult Apron		8611049
Pano-Cape Api	on	8611053
Technician Apr	on	8611056



Soothe-Guard®

Lead lined aprons have a unique silky-smooth material with an exquisite feel. Sturdy sewn-in hanging loops provide easy storage and a new collar design fits better, for less gapping and greater patient comfort. Soothe-Guard material front and back, makes these aprons easy to clean (just wipe with damp cloth or disinfectant wipe).

Soothe-Guard Lead Lined Aprons 0.3mm Lead Equivalency

	ipi ene elemin Leau Equinaterio)	
	Adult Apron with Collar	6611050
	Child Apron with Collar	6611047
	Protective Collar	6611051
	Pano-Cape Apron	6611053
0.5mm Lead Equivalency*		
	Adult Apron	5611049
	Pano-Cape Apron	5611053

*0.5mm lead equivalency aprons are made-to-order (not returnable). 4 to 6 weeks lead time.



Apron Hanger Heavy-duty, white epoxy-coated steel. Mounting hardware is included. Dimensions: 32 x 378 x 35mm

Apron Hanger



X-Ray Accessories

Film Mounts







Product Code Top quality film mounts from DENTSPLY RINN®, the leading supplier of film mounts worldwide since 1923.

Masked vinyl pocket

- EZ-View[®] Masked Pocket 21, 22, 23 Series
- Available in Black (21, 23 Series) or Green (22 Series)
- Masking provides the most comfortable viewing, and enhances viewing of fine details
- 23 Series Open-Window Masked Mounts let you have it both ways the easyloading of a pocket mount, the clarity of an open window
- Can be duplicated with exceptional clarity
- Patient name/date even appear on the duplicate
- Can be written on with pen and pencil

Open window paperboard Eezeemount® Paperboard 10 Series

- Films slide in easily, stay securely in die-cut windows
- Sturdy, reusable paperboard construction
- Accepts any pen or pencil

Open Window Vinyl EZ-Tab® Plastic 30 Series

- X-rays snap quickly into tabs for easy loading
- Chemical and water-resistant vinyl permits processing in the mount
- Recessed windows keep radiographs safe

Notched/Die-Cut TransView® 25 Series

- Light-diffusing vinyl offers full visibility of radiographs
- Films held securely by easy-opening die-cuts
- Thin profile requires minimum filing space

See <u>next page</u> for the full range of available RINN Film Mounts.

EZ-View[®] FasTear[®] Tear-Apart Mount - Pull out the desired number of pockets and tear into desired configuration.

EZ-View Fas Iear Universal Mounts	200500R
EZ-View [®] Universal Cut-Apart Mount - Pull out the desired number and cut into desired configuration. The clear vinyl pocket style offers protection, visibility, and convenience.	
EZ-View Universal Cut-Apart Mounts, (Pkt 1000) 2	2001021000
Protective Mount Envelopes - Heavy-duty, end-opening grey paper e	nvelones

Protective Mount Envelopes - Heavy-duty, end-opening grey paper envelopes for filing x-ray mounts. Large enough for full mouth surveys.

RINN Mount Envelopes, Size #20, (Pkt 100)

500849

.

Marking Pencils - White - A unique white pencil that writes easily on film mounts and x-rays.

Marking Pencil, (Pkt 3)

Film Mounts

			Eezeemount [®] Paperboard 10: Grey	EZ-Tab [®] Open Window 30: Grey	TransView [®] Notched 25: Translucent	EZ-View [®] Masked Pocket 21: Black	EZ-View [®] Masked Pocket 22: Green	EZ-View [®] Masked Pocket 23: Open Window
IH #2			100101 76 x 111mm			216011 67 x 95mm	226011 67 x 95mm	236011 67 x 95mm
V #2			100113 76 x 111mm					
2H		#0	100121 76 x 111mm					
		#1	100105 76 x 111mm					
		#2	100119 76 x 111mm	306021 70 x 101mm		216021 67 x 95mm	226021 67 x 95mm	236021 67 x 95mm
		#3					227021 67 x 121mm	
3H #2			100134 76 x 165mm					
3V #2			100135 76 x 143mm					
4H #2			100142 76 x 213mm					236040 67 x 181mm
4H #2							226045 92 x 101mm	
4V #2								236042 70 x 140mm
6H #2			100144 111 x 160mm					
6V #2							228062 120 x 120mm	
Full-Mouth	Series Mounts (FM)	()						
14 Windows 8H #2 6V #2					250184 14WT 108 x 289mm			
16 Windows 6H #2 6V #2 4BW #2							228162 136 x 206mm	
18 Windows 8H #2 6V #2 4BW #2						216180 [2 folds] 121 x 292mm		
20 Windows 3H #2 3V #1 (5/3) 4VBW #2						216203 [2 folds] 121 x 327mm		
	istory Mounts							
10BW #2							222345 229 x 92mm	232345 229 x 92mm
	ilm Mounts		100107					
1H OCC			100107 111 x 101mm					





Small Equipment





Universal Collimator

Protect your patients by limiting the size of the primary x-ray beam to #2 film.

Quickly applies into open end of round PID

X-Ray Accessories

- Vertical or horizontal use
- Use with any RINN[®] Arm & Ring positioning device for accurate alignment

Universal Collimator

540853

530105

530108

Product Code

RINN® Chairside Darkrooms

The Chairside Darkroom allows film to be conveniently processed in the surgery. The Chairside Darkroom features a light filtering top, allowing clear vision, an ergonomic design and a chemical-resistant surface.

Chairside Darkroom

Includes:

- 1 x Ruby Light Filter
- 4 x Single-Clip film Hangers
- 4 x 240ml Beakers with attached snap-on lids

Mini-Chairside Darkroom

Includes:

- 1 x Amber Light Filter
- 2 x Single-Clip film Hangers
- 3 x 240ml Beakers with attached snap-on lids

Replacement Beakers for both units (4 x 240ml) 530255



Film Processing Hangers

Single Clip (1)	length 57mm	510258
Single Clip (12)	57mm	51025812
2 Clips	124mm	510972
4 Clips	176mm	510973
6 Clips	228mm	510260
10 Clips	214mm	510261
12 Clips	276mm	510262
14 Clips (Short)	283mm	510263



Small Equipment



Mini-Ray Duplicator

Duplicate your intraoral films with this handy battery-powered film duplicator. For use with the RINN $^{\circ}$ Chairside Darkrooms or daylight loaders of automatic processors.

- Duplicates in 8-10 seconds.
- Accommodates #0, #1, #2 and #3 films
- Includes a low battery indicator light.

Mini-Ray Duplicator

721220

670520

670542

Product Code



EZ-View® Auto Light Flatscreen Viewer

Save time and energy with the auto-on feature. Automatically turns the viewer on and off as film mounts or x-rays are slid under and out from the holder bar.

- Cold Cathode Tube light allows hours of life
- 30 adjustable brightness settings
- Mount on the wall or countertop
- Viewing area: 6 x 13.5 (152 x 343mm)
- Dimensions: 8 x 16.5 x 1.25 (216 x 419 x 32mm)
- Enhance analysis of radiographs with the 2x magnifier

EZ-View Auto Light Flatscreen Viewer EZ-View Auto Light Flatscreen Viewer with Magnifier 230V

Universal Viewer with Magnifier

- Standard fluorescent light source
- 2 x magnifier
- Easy to clean surfaces that will not rust when your diagnostic film isn't fully dry

Universal View with Magnifier, 230-240 V	670445
Magnifier Assembly	003704





X-Ray Accessories

Kodak Film and Chemicals



		Product Code
Kodak* Film		
Occlusal		
	DF 49	1666122
	DF 50	1666163
	IO 41P	1169143
Periapical		
	DF 54	1228840
	DF 56	1273747
	DF 57	1753664
	DF 58	1658194
	DF 58C	1717131
	DF 54C	1445360
	DF 57C	8987877
	IP 01P	8675332
	IP 21C	1539931
	IP 21P	1163401
	IP 22C Double Film	1560390
Bitewing		
	DF 40	8393043
	DF 42	1296771
	DF 48	8348658
	IB 31P	8298929
T-MAT		
	T-MAT G Film (50) 15x30cm	1869080
Pano		4000746
	Panoramic Film 15x30 DF76 50SH	1990712

Developer / Fixer

Ready Matic Developer / Fixer	1028869
Rapid Developer / Fixer Pack	1838374
Developer GBX 500ml	4037180
Developer GBX 5ltr	4037206
Fixer GBX 500ml	4037214
Fixer GBX 5ltr	4037230

*Kodak is not a registered trademark of the DENTSPLY International Group of Companies.

Notes

Client Services: Australia: 1300 55 29 29 | New Zealand: 0800 DENTSPLY (33 68 77)