

Bonding materials

VOCO GmbH, Department of Knowledge Communication

Anton-Flettner-Str. 1-3

Postfach 767

D-27472 Cuxhaven

Tel.: +49 (0)4721-719-0

Fax: +49 (0)4721-719-109

info@voco.de

www.voco.de



VOCO provides the dentist with a number of different bonds. The bonds range from a 4th generation, 2-bottle total etch bond to a 7th generation, 1-bottle self-etch. A summary of the different bonds is presented below.

Product	Etching strategy	Layers	Polymerisation	Compatibility
Solobond Plus	Total-Etch	2 (primer, then adhesive)	LC	All composites (A second layer of adhesive is recommended with DC/SC composites)
Solobond M	Total-Etch	1	LC	All light-curing composites
Admira Bond	Total-Etch	1	LC	All light-curing composites
Futurabond NR	Self-Etch	1 (primer and adhesive are mixed in advance)	LC	All light-curing composites
Futurabond DC	Self-Etch	1 (primer and adhesive are mixed in advance)	DC	All composites
Futurabond M	Self-Etch	1	LC	All light-curing composites

Conclusion: The range of 4 generations of bonds offers every dentist a bond that is his taste.