

# *elite cord*



**Clinical**  
Retraction System

**Zhermack**   
Dental

# elite cord

simply beyond precision



Zhermack, leader in the production of avant-garde dental material for impression taking, presents the new range of gingival displacement cord, **elite cord**, in combination with Zhermack silicones, able to optimise impression precision simply and quickly.



## EASY TO APPLY

Thanks to the knitted structure, the cord remains in position, enabling ideal non-traumatic displacement to be obtained rapidly

## AVAILABLE IN FIVE DIFFERENT THICKNESSES (000, 00, 0, 1, 2)

Enables gingival opening to be obtained in a range of conditions and with different gum types

## HIGH ABSORBENCY

The 100% cotton fibre keeps the marginal gum tissue temporarily dry during restoration, margin finishing and adhesive cementation

## VERSATILE (can be impregnated as required)


Can be used non impregnated, or impregnated with the preferred haemostatic agent for optimal control of bleeding

## CONTRASTING COLOURS


Optimise visibility of the operating field and prevent oversights

## PRECISION


The first time, every time



**STEP 1** Choose the technique and sizes to be used  
**STEP 2** Take the cord  
**STEP 3** Position the cord in the gingival sulcus



**STEP 4** Leave the cord in the sulcus for the time necessary  
**STEP 5** Wash before removing the cord  
**STEP 6** Remove the cord from the sulcus



## Packaging

<b>C610000</b>	Elite Cord # 000 (1 x 275 cm)	<b>C610001</b>	Elite Cord Kit # 000 (5 x 275 cm)	<b>C610051</b>	Elite Cord Kit # 000-00-0-1-2 (5 x 275 cm)
<b>C610010</b>	Elite Cord # 00 (1 x 275 cm)	<b>C610011</b>	Elite Cord Kit # 00 (5 x 275 cm)		
<b>C610020</b>	Elite Cord # 0 (1 x 275 cm)	<b>C610021</b>	Elite Cord Kit # 0 (5 x 275 cm)		
<b>C610030</b>	Elite Cord # 1 (1 x 275 cm)	<b>C610031</b>	Elite Cord Kit # 1 (5 x 275 cm)		
<b>C610040</b>	Elite Cord # 2 (1 x 275 cm)	<b>C610041</b>	Elite Cord Kit # 2 (5 x 275 cm)		

## BRANCHES

### GERMANY

Zhermack GmbH Deutschland  
 Öhlmühle 10 - D-49448 Marl

T +49 (0) 54 43 / 20 33 - 0  
 F +49 (0) 54 43 / 20 33 - 11  
 info@zhermack.de  
 www.zhermack.de

### POLAND

Zhermapol Sp. z o.o.  
 ul. Augustówka 14  
 02 - 981 Warszawa - Polska

T +48 22 858 82 72  
 +48 22 292 71 17  
 F +48 22 642 07 14  
 biuro@zhermack.com  
 www.zhermapol.pl

### USA

Zhermack Inc.  
 P.O. Box 4195, River Edge  
 New Jersey 07661 - 4195

T +1 (877) 819 6206  
 +1 (732) 389 8540  
 F +1 (732) 389 8543  
 info@zhermackusa.com  
 www.zhermackusa.com

Zhermack S.p.A. | Via Bovazecchino, 100 | 45021 Badia Polesine (RO) ITALY  
 T +39 0425 597611 | F +39 0425 53596 | info@zhermack.com | www.zhermack.com

Zhermack