## **HIRO ROCK®**

Type 4 modified synthetic dental die stone with high strength.

Type 4 modified synthetic dental die stone with high strength, consisting of a blend of refined selected alpha hemihydrates and special resin, for implant work, sectioned, master and static models, all kinds of prosthetics with high precision etc.

- Top price-performance ratio;
- Thixotropic Exceptional hardness;
- Extreme easy processing, smooth and creamy mixing;
- Dimensional stability;
- It can be scanned with several scanning systems;
- Glossy surface Low expansion.



• Water/Powder Ratio: 20ml/100g

• Setting Time: 10 - 12 min.

• Expansion after 2 h:  $\leq 0.1\%$ 

Compressive Strength after 1 hr:

• Compressive Strength [MPa; N/mm<sup>2</sup>] - 50MPa (= 7251psi) Compressive Strength after 24 hrs:

• Compressive Strength [MPa; N/mm<sup>2</sup>] - 89MPa (=12908 psi) ISO 6873, Data referred to 20ml/100g



**Packaging:** 2.5kg | 5kg | 10kg | 20kg



- Sunglow Yellow > Golden Brown
- Pistachio Green

## **HIRO HARD ROCK®**

Type 4 modified synthetic dental die stone with high strength.

Type 4 modified synthetic dental die stone with high strength, consisting of a blend of refined selected alpha hemihydrates and special resin, perfect for implant work, sectioned and master models, all kinds of prosthetics with high precision, inlay onlay etc.

• Thixotropic

• High innitial compression strength, resists chipping and splintering • Low expansion • Can be scanned with several scanning systems • Dimensional stability and glossy surface

- Smooth and creamy mix Glossy surface
- Ideal also for plastic basing systems
- Ideal for CAD/CAM or other accuracy required purposes

Compressive Strength after 1 hr: • Compressive Strength [MPa; N/mm<sup>2</sup>] - 55MPa (=7977psi) Compressive Strength after 24 hrs:

• Compressive Strength [MPa; N/mm<sup>2</sup>] - 99MPa (=14358psi) ISO 6873, Data referred to 20ml/100g



## **Color:**

- ▶ Golden Brown ▶ Vanilla Yellow
- Pistachio Green



Water/Powder Ratio: 20ml/100g

- Setting Time: 10 12 min.
- Expansion after 2 h:  $\leq$  0.08 %

**VDXTD** 

www.mutsumidental.eu