


**tempolink®**  
proper formulation for specific requirements



**3 new types of temporary cements for cementing temporary crowns & bridges:**

**tempolink® clear** is an eugenol-free, resin based for trial cementing restorations (crowns, bridges, inlays and onlays). The colorless, transparent material especially suits for cementation in the anterior region. By building a film it seals the restoration reliably, yet ensures easy removal when desired. The dual curing systems guarantees flexible processing. Excesses can easily be removed from the restoration and preparation.

**tempolink® single & multi** are self-curing, eugenol-free, polycarboxylate-based cements, intended for both short- and long term crowns & bridges temporaries, inlays and onlays. The hydrophilic material with excellent fluidity enables the prothesis to be easily and completed placed, yet affords easy removal when desired.

**tempolink® single** is designed especially for the cementation of single crowns and smaller applications, **tempolink® multi** should be used for polynominal provisionals resp. larger restorations.

The 3 types are antibacterial, odour- and tasteless and available in functional mini-mix cartridges. Each standard package contains one cartridges of 5 ml and 8 mixing canulas.

