


mis[®] | C1

A Conical Connection Implant

Enhanced
Geometry



Esthetics






Simplicity



Enhanced geometry


Conical connection




- Friction fit
- Ideal seal
- Reduced micro-movements

Enhanced geometry

Thread design



Thick threads:
Bone condensing
properties



Sharp threads:
Bone cutting
properties




Enhanced geometry

Thread design


Implant apical shape



Anterior Socket




Posterior Socket



Enhanced geometry

Platform switch


Platform switch design
Designed for soft tissue
preservation and growth



Esthetics

Prosthetics


- Broad range of all restorative solutions
- Concave emergence profile abutments
- Uncompromised accuracy
- Gold shade
- Color-coded



Esthetics

Prosthetics

A consistent emergence profile throughout all prosthetic components, allows an efficient restorative procedure; leading to favorable esthetic results.



Simplicity

Range

Narrow Platform

Ø 3.30



Ø 2.50

Standard Platform

Ø 3.75



Ø 3.10

Ø 4.20



Ø 3.60

Wide Platform


Ø 5



Ø 4.50

Simplicity

Final drill



Simplicity

Combo package



Simplicity

Surgical kit

Minimum tools for
maximum efficiency.



Simplicity

Step drills

These drills are designed with a "step" matching the diameter of the previous drill used in the sequence, ensuring a highly accurate drilling procedure.



mis[®]
THANK YOU

