

glassbite – the new clear class

from 0 to 80 in 60 sec.



At the IDS 2007 DETAX presents the first clear-transparent bite registration materials, with its highest final shore hardness of 80 Shore A – and such in just 60 seconds.

From the Automix2 cartridge **glassbite** is applied directly on the toothrow, being maximum comfort for patient and dentist. The optimal consistency prevents resistance when the patient closes his mouth – thus excluding displacements during checkbite. Only 60 seconds in the mouth makes it most convenient for the patient; due to the extrem high shore hardness it is easy to trim (not flexible, not brittle) to eliminate excess material. The thixotropic material is designed not to flow into interdental spaces, for positioning control onto occlusal.

glassbite fast-set accurately reproduces the fundamental occlusal details. It can be also used as matrix material for the light curing of composites, for positioning of auxiliaries or for orthodontic applications.

glassbite is packaged in a standard box (2 Automix2 cartridges of 50 ml and, 4 mixing cannulas), also available as eco-pack (4 Automix2 cartridges, 16 mixing cannulas). Free in each packages are contouring tips for a precise and material-saving application.