



RESULT PRODUCT TEST

# ROEKO GuttaFlow® bioseal



## Barometer | Testphase Readers Product Test – The Original



### ACTIVE REGENERATION IN THE ROOT CANAL

GuttaFlow bioseal is an innovative filling material. Upon contact with fluids, it provides bioactive material such as calcium and silicate and activates biochemical processes. By formation of hydroxylapatite crystals the bioactive material supports adhesion and regeneration of bone and dentin tissue.


### 3-IN-1-OBTURATION MATERIAL

- Combines free-flow gutta-percha with a sealer and bio ceramic at room temperature.
- Uncomplicated, easy flowing filling material.
- Tight seal.
- Working time: 5 minutes; Curing: 12 – 16 minutes.
- Excellent radiopacity.

### PRODUCT RATING


#### Most Convincing Features for the Readers

- Simple Handling
- Flowability
- Support for Regeneration
- Radiopacity
- Overall Result
- Short Setting Time




### Consistency

- Ideal
- Acceptable
- Highly Viscous (Hard)
- Lowly Viscous (Soft)




### Flow Properties

- Excellent
- Good
- Sufficient
- No Details Available




### Treatment Results

- Very Satisfied
- Satisfied
- Less Satisfied
- No Details Available



### Further Usage

- Yes
- Perhaps
- No



### CONCLUSION

The smart product features of GuttaFlow bioseal have convinced the testers. Many users were particularly impressed by the simple handling and the excellent consistency and flow properties. The good radiopacity and bioactive features also received positive ratings from the dentists. 85% of the GuttaFlow bioseal testers were either very satisfied or satisfied with the treatment results. 97% could envisage the further usage of GuttaFlow bioseal in their practice.