



Class B, simply profitable

NEW



01

Automatic digital
cycle recording
on USB stick



02

Type B cycle
in less than
30 minutes

ECO-B

OK

03

Automatic
water filling
function

lina

Technical data

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

Lina sterilization cycles

		B-Universal 134			B-Prion 134			B-Universal 121
Temperature		134 °C			134 °C			121 °C
Holding time		3.5'			18'			15'
Total cycle duration including drying		ECO-B	Full Load	Typical Load	ECO-B	Full Load	Typical Load	Full Load
	Lina 17	29'	54'	49'	44'	69'	64'	75'
	Lina 22	33'	64'	55'	48'	79'	70'	79'
Max solid: ECO/full/typical loads/ Max porous load		Lina 17: 0.5kg / 4kg / 2kg / 1.5kg • Lina 22: 0.5kg / 5kg / 2kg / 1.8kg						
Test cycles		Helix/Bowie Dick/Vacuum test						

Your advantages at a glance

- > The ECO B cycle reduces the cycle time by half
- > Low-noise operation with the patented water separation system
- > User-friendly keypad for simple navigation to all sterilization options
- > Built-in dust filter to protect internal components
- > Programmable cycle start
- > Excellent value for money
- > Integrated automatic water filling valve (optional)
- > Automatic digital cycle recording on USB stick (optional)



	Lina 17	Lina 22
Type:	PR013-003-17	PR013-003-22
Chamber size:	17 l	22 l
Power supply:	200 – 240 Vac 50 / 60 Hz; 10 A	
Power consumption	2 kW	
Dimensions overall (W x H x D):	450 x 435 x 599 mm	
Weight:	39 kg	40 kg
Main/Used water tanks:	4 / 3.5 litres	
Working range:	5 to 10 cycles	
Usable space in chamber (W x H x D):	195 x 195 x 297 mm	195 x 195 x 390 mm
Connection type:	1 serial port, 1 USB + 1 serial (optional), integrated automatic water filling valve (optional)	

Lina sterilizers were designed, certified and validated with the most stringent Directives and Standards:

93/42/EEC	Medical Device Directive	EN 61326	Electromagnetic compatibility
97/23/EEC	Pressure equipment Directive	EN 61010-1	Safety requirements
ISO 17665-1	Validation and routine control	EN 61010-2-040	Specific requirements for steam sterilizer
EN 13060	Small steam sterilizer		

W&H Sterilization Srl.
Via Bolgara, 2
24060 Brusaporto (BG), Italy
t + 39 035-66 63 000
f + 39 035-50 96 988
office.sterilization@wh.com
wh.com