Sintraumal+Surgi-Fil by roeko

Atraumatic needle-suture sets and sterile single sutures



Sintraumal + Surgi-Fil rоеко

Surgi-Fil + Sintraumal are gentle on the tissue:

- The sutures have a hydrophobic coating and are non-capillary. As neither secretion nor water is absorbed by the suture, there is far less risk of infection.
- Sintraumal sutures do not cut into the tissue due to their special coating, thus reducing trauma.


Surgi-Fil + Sintraumal are safe and easy to use:


- Sterile in a double peel-pack
- Easy and secure knotting
- The polyester suture is green and the silk black for easy identification in the tissue
- The tangle free suture can be quickly removed from the spool.

Sintraumal

Atraumatic needle-suture sets

Sintraumal is a combination of precision ground, eyeless needles and non-resorbable silk or polyester sutures. The stainless steel needles are smoothly bound to the sutures, so pain and tissue damage often caused by needles with eyes, cannot occur.







Surgi-Fil

Surgical sutures without needles

<u>Surgi-Fil silk sutures</u> are made from specially braided silk, a natural fibre. The tissue compatible impregnation makes the suture non-capillary.

The specially braided thread is particularly suitable for sub-cutaneous stitches, ligatures and purse-string ligatures as well as for closing wounds. The thread is coloured black with haematein for improved visibility.





<u>Surgi-Fil polyester sutures</u>, Flexafil, are braided. Their hydrophobic coating makes them highly tissue compatible and prevents wicking.

The excellent properties of Flexafil sutures make it suitable for ligatures, purse-string ligatures, vessel occlusion and for closing wounds. They are coloured green, to differentiate them from the silk sutures.

| Indication | Needle 1:1 True sizes | Needle type | Gauge metric | USP | Suture 45 cm | REF |
|---|--------------------------|-------------|------------------------|-----|-----------------|---------|
| Operations in firm gum tissue | | HS 20 | 1.5 | 4/0 | Silk, black | 490 301 |
| Extractions | | HS 25 | 2.0 | 3/0 | Silk, black | 490 302 |
| Apical extractions and flap operations | | DS 20 | 1.5 | 4/0 | Silk, black | 490 303 |
| Extractions | | HS 25 | 2.5 | 2/0 | Silk, black | 490 305 |
| Extractions, e.g. Wisdom teeth, Serial extractions, Stiches in firm tissues | | HRT 18 | 2.0 | 3/0 | Flexafil, green | 490 306 |
| Flap operations, Periodontal surgery, Gingivectomy | | GR 20 | 1.5 | 4/0 | Silk, black | 490 307 |
| Extractions, e. g. Wisdom teeth, Flap operations | | AS 30 | 2.0 | 3/0 | Silk, black | 490 308 |

Contents: 12 pc.

Sintraumal is available in 7 combinations, the right combination for each case: Available needles are round bodied needle, cutting needle, round bodied cutting needle (trocar finish). Available as half circle, 3/8 circle, straight and fish hook.





roeко - optimizing surgical procedures



resorbable gelatine sponge with colloidal silver



Roeko-Swab sterile gauze swabs

Order form

Sintraumal Surgi-Fil **Needle-suture sets Surgical sutures** metric USP Cont. REF Qty. metric USP Cont. REF Qty. HS 20 1.5 4/0 Silk, black 12 St. 490 301 4.0 Silk, black 12 St. 490 409 1 HS 25 2.0 3/0 Silk, black 12 St. 490 302 3.5 0 Silk, black 12 St. 490 410 Silk, black DS 20 1.5 4/0 Silk, black 12 St. 490 303 2.5 12 St. 490 411 2/0 490 412 HS 25 2.5 2/0 Silk, black 12 St. 490 305 2.0 3/0 Silk, black 12 St. HRT18 2.0 3/0 Flexafil, green 12 St. 490 306 1.5 4/0 Silk, black 12 St. 490 413 GR 20 1.5 4/0 Silk, black 12 St. 490 307 2.0 3/0 Flexafil, green 12 St. 490 421 Silk, black 4/0 490 422 AS 30 2.0 3/0 12 St. 490 308 1.5 Flexafil, green 12 St.

length 45 cm

length 70 cm

www.coltenewhaledent.com © 2004 Coltène/Whaledent

> **Coltène/Whaledent GmbH + Co. KG** Raiffeisenstraße 30

Rainesensuade 30 89129 Langenau/Germany Tel. +49 (0)7345 805 0 Fax +49 (0)7345 805 201 productinfo@coltenewhaledent.de Coltène/Whaledent Inc.

235 Ascot Parkway Cuyahoga Falls, Ohio 44223/USA Tel. +1 330 916 8800 Fax +1 330 916 7077 productinfo@coltenewhaledent.com

Coltène AG

Feldwiesenstraße 20 9450 Altstätten/Switzerland Tel. +41 (0)71 757 53 00 Fax +41 (0)71 757 53 01 office@coltenewhaledent.ch

